

<b>REF</b>	501-2 6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期: 2024-09 1 / 6
------------	----------------	------------	------------	------------	---------	------------------------

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(Li)锂	(CRBM)卡马西平	(DGXN)地高辛
(PHYT)苯妥英	(THEO)茶碱	(VALP)丙戊酸
(TT3)总三碘甲状腺氨酸	(FT3)游离T3	(TT4)总T4
(FT4 )游离T4	(TSH)促甲状腺素	(Cortisol)皮质醇
(FSH)促卵泡刺激素	(LH)促黄体生成素	(PRL)催/泌乳素
(TES)睾酮	(E2)雌二醇	(C-P)C-肽
(FA)叶酸	(INS)胰岛素	(VB12)维生素B12
(IgG)免疫球蛋白 IgG	(IgM)免疫球蛋白 IgM	(IgA )免疫球蛋白 IgA
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(Ferritin)铁蛋白
(F-PSA)游离前列腺特异性抗原	(FE3/uE3)游离雌三醇	(17-OHP)17-羟孕酮
(ASD)雄烯二酮	(ALD)醛固酮	(DHEA-S)硫酸脱氢表雄酮
(HGH)人生长激素	(25-OH-Vit D)25羟基维生素D	(ACTH)促肾上腺皮质激素
(AK)阿米卡星	(Gentamicin)庆大霉素	(Phenobarbital)苯巴比妥
(Tobramycin)妥布霉素	(Vancomycin)万古霉素	(CT)降钙素
(Acetaminophen)对乙酰氨基酚	(cPSA)复合前列腺特异性抗原	(Salicylate)水杨酸盐

<b>REF</b>	501-2	6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期:	2024-09
								2 / 6

批号: 2FCCEA6			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
<b>(17-OHP)17-羟孕酮</b>			
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.06	1.64-2.48
<b>(25-OH-Vit D)25羟基维生素D</b>			
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	ng/mL	26.3	21-31.6
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA )	nmol/L	65.8	52.6-79
<b>(Acetaminophen)对乙酰氨基酚</b>			
Abbott Alinity C Series Abbott \ Colorimetry	mg/L	34.2	27.4-41
<b>(ACTH)促肾上腺皮质激素</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pg/mL	197	158-236
<b>(AFP)甲胎蛋白</b>			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	87.3	69.8-104.8
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	97	77.6-116.4
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	97.8	78.2-117.4
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	92.1	73.7-110.5
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	94	75.2-112.8
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	89.8	71.8-107.8
<b>(AK)阿米卡星</b>			
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	mg/L	13.5	10.8-16.2
<b>(ALD)醛固酮</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pg/mL	262	210-314
<b>(ASD)雄烯二酮</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.61	1.29-1.93
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.05	0.83-1.27
<b>(CEA)癌胚抗原</b>			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	16.7	13.4-20
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	16.7	13.4-20
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	19.4	15.5-23.3
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	18.6	14.9-22.3
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	14.7	11.8-17.6
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	17.6	14.1-21.1
<b>(Cortisol)皮质醇</b>			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	162	130-194
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	150	120-180
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	159	127-191
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	188	150-226
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	182	146-218
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	171	137-205
<b>(C-P)C-肽</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.83	5.46-8.2
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.3	3.44-5.16
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.93	4.76-7.1
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	6.14	4.92-7.36
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.47	4.37-6.57

<b>REF</b>	501-2	6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期:	2024-09
								3 / 6

<b>(cPSA)复合前列腺特异性抗原</b>			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/L	780	624-936
<b>(CRBM)卡马西平</b>			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	µg/mL	7.93	6.34-9.52
<b>(CT)降钙素</b>			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/L	195	156-234
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/L	179	143-215
<b>(DGXN)地高辛</b>			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.2	1.76-2.64
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.2	1.76-2.64
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.37	1.89-2.85
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	ng/mL	2.19	1.75-2.63
<b>(DHEA-S)硫酸脱氢表雄酮</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1410	1128-1692
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1450	1160-1740
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1740	1392-2088
<b>(E2)雌二醇</b>			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pmol/L	569	455-683
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	587	470-704
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	848	678-1018
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	895	716-1074
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	602	482-722
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	815	652-978
<b>(FA)叶酸</b>			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	10.1	8.1-12.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	9.84	7.87-11.81
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	11.1	8.9-13.3
<b>(FE3/uE3)游离雌三醇</b>			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.94	3.16-4.72
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	4.36	3.48-5.24
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.94	3.14-4.74
<b>(Ferritin)铁蛋白</b>			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	150	120-180
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	140	112-168
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	185	148-222
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	169	135-203
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	271	217-325
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	224	179-269
<b>(F-PSA)游离前列腺特异性抗原</b>			
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.65	1.32-1.98
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.9	1.52-2.28
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.54	1.24-1.84
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.48	1.18-1.78
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.2	1.76-2.64

<b>REF</b>	501-2	6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期:	2024-09
								4 / 6

(FSH)促卵泡刺激素			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	25.1	20.1-30.1
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	21	16.8-25.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	25.1	20.1-30.1
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/L	27.1	21.7-32.5
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	IU/L	35.7	28.6-42.8
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	23.6	18.9-28.3
(FT3)游离T3			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pmol/L	15.4	12.3-18.5
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	12.3	9.8-14.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	11.8	9.4-14.2
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	9.73	7.78-11.68
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	11.7	9.4-14
(FT4 )游离T4			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pmol/L	33.6	26.9-40.3
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	pmol/L	39.8	31.8-47.8
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pmol/L	29.3	23.4-35.2
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	pmol/L	26.4	21-31.8
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	32.2	25.8-38.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pmol/L	30.2	24-36.4
(Gentamicin)庆大霉素			
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L	8674	6922-10426
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	10408	8326-12490
(hCG)人绒毛膜促性腺激素			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	IU/L	34.1	27.3-40.9
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	28.8	23-34.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	31.8	25.4-38.2
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	31.5	25.2-37.8
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	35.4	28.3-42.5
(HGH)人生长激素			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	µg/L	7.98	6.38-9.58
(IgA)免疫球蛋白 IgA			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	mg/dL	176	141-211
(IgE)免疫球蛋白 IgE			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/mL	99.7	79.8-119.6
(IgG)免疫球蛋白 IgG			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	mg/dL	719	575-863
(IgM)免疫球蛋白 IgM			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	mg/dL	72.7	58.2-87.2
(INS)胰岛素			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	mIU/L	68.6	54.9-82.3
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	64.1	51.3-76.9
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	47.3	37.8-56.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	mIU/L	69.5	55.6-83.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	51.4	41.1-61.7

<b>REF</b>	501-2	6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期:	2024-09
								5 / 6

(LH)促黄体生成素			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	IU/L	23.2	18.6-27.8
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	IU/L	17.6	14.1-21.1
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	IU/L	15.1	12.1-18.1
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	IU/L	26.2	21-31.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	18	14.4-21.6
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	IU/L	18	14.4-21.6
(Li)锂			
Abbott Alinity C Series Abbott \ ISE (direct)	mmol/L	1.16	0.92-1.4
(Phenobarbital)苯巴比妥			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	nmol/L	121983	97586-146380
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	862	690-1034
(PHYT)苯妥英			
Abbott Alinity C Series Abbott \ Enzyme immunoassay (EIA)	µg/mL	12.3	9.8-14.8
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	µg/mL	15.2	12.2-18.2
(PRL)催/泌乳素			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	24	19.2-28.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	16.2	13-19.4
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	17.1	13.7-20.5
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	14.6	11.7-17.5
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	22.3	17.8-26.8
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	23	18.4-27.6
(PSA)前列腺特异性抗原			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.8	2.24-3.36
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.32	1.86-2.78
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.72	2.18-3.26
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.51	2.01-3.01
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.28	1.82-2.74
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.54	2.04-3.04
(Salicylate)水杨酸盐			
Abbott Alinity C Series Abbott \ Enzymatic	mg/dL	15.4	12.3-18.5
(TES)睾酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	0.003	0.002-0.004
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.85	3.08-4.62
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.69	2.95-4.43
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	3.91	3.13-4.69
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	5.39	4.31-6.47
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.05	1.63-2.47
(THEO)茶碱			
Abbott Alinity C Series Abbott \ Enzyme immunoassay (EIA)	µg/mL	14.2	11.4-17
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	µg/mL	13.3	10.6-16
(Tobramycin)妥布霉素			
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L	6338	5058-7618
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	5011	4029-5993

<b>REF</b>	501-2	6x5mL	<b>EXP</b>	2027-06-11	<b>LOT</b>	2FCCEA6	更新日期:	2024-09
								6 / 6

(TSH)促甲状腺素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	4.52	3.62-5.42
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	mIU/L	4.84	3.88-5.8
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	mIU/L	5.66	4.52-6.8
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	mIU/L	5.24	4.2-6.28
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	mIU/L	4.15	3.31-4.99
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	mIU/L	5.25	4.19-6.31
(TT3)总三碘甲状原氨酸			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.8	1.44-2.16
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.01	1.61-2.41
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	1.85	1.48-2.22
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.19	1.75-2.63
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	2.05	1.63-2.47
(TT4)总T4			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	ng/mL	113	90-136
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	116	93-139
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	ng/mL	114	91-137
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA )	ng/mL	73.4	58.9-87.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	ng/mL	114	91-137
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	ng/mL	136	109-163
(VALP)丙戊酸			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	µg/mL	81.3	65-97.6
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	µg/mL	79	63.2-94.8
(Vancomycin)万古霉素			
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L	9867	7894-11840
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	12903	10322-15484
(VB12)维生素B12			
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA )	pg/mL	598	478-718
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pg/mL	456	365-547
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA )	pg/mL	658	526-790
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )	pg/mL	447	358-536
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA )	pg/mL	588	470-706

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。