

REF	501-1 6x5mL	EXP	2027-06-11	LOT	2FCA796	更新日期: 2024-09 1 / 6
------------	----------------	------------	------------	------------	---------	------------------------

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(Li)锂	(CRBM)卡马西平	(DGXN)地高辛
(PHYT)苯妥英	(THEO)茶碱	(VALP)丙戊酸
(TT3)总三碘甲状腺氨酸	(FT3)游离T3	(TT4)总T4
(FT4)游离T4	(TSH)促甲状腺素	(Cortisol)皮质醇
(FSH)促卵泡刺激素	(LH)促黄体生成素	(PRL)催/泌乳素
(TES)睾酮	(E2)雌二醇	(C-P)C-肽
(FA)叶酸	(INS)胰岛素	(VB12)维生素B12
(IgG)免疫球蛋白 IgG	(IgM)免疫球蛋白 IgM	(IgA)免疫球蛋白 IgA
(IgE)免疫球蛋白 IgE	(AFP)甲胎蛋白	(CEA)癌胚抗原
(hCG)人绒毛膜促性腺激素	(PSA)前列腺特异性抗原	(Ferritin)铁蛋白
(F-PSA)游离前列腺特异性抗原	(FE3/uE3)游离雌三醇	(17-OHP)17-羟孕酮
(ASD)雄烯二酮	(ALD)醛固酮	(DHEA-S)硫酸脱氢表雄酮
(HGH)人生长激素	(25-OH-Vit D)25羟基维生素D	(ACTH)促肾上腺皮质激素
(AK)阿米卡星	(Gentamicin)庆大霉素	(Phenobarbital)苯巴比妥
(Tobramycin)妥布霉素	(Vancomycin)万古霉素	(CT)降钙素
(Acetaminophen)对乙酰氨基酚	(cPSA)复合前列腺特异性抗原	

REF 501-1	6x5mL	EXP 2027-06-11	LOT 2FCA796	更新日期: 2024-09
				2 / 6

批号: 2FCA796				
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD	
(17-OHP)17-羟孕酮				
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.54	0.432-0.648	
(25-OH-Vit D)25羟基维生素D				
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20.2	16.2-24.2	
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	nmol/L	50.5	40.4-60.6	
(Acetaminophen)对乙酰氨基酚				
Abbott Alinity C Series Abbott \ Colorimetry	mg/L	10.6	8.5-12.7	
(ACTH)促肾上腺皮质激素				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pg/mL	24.6	19.7-29.5	
(AFP)甲胎蛋白				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.3	14.6-22	
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.5	14-21	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.2	13.8-20.6	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.1	15.3-22.9	
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18	14.4-21.6	
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.2	13.8-20.6	
(AK)阿米卡星				
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	mg/L	5	4-6	
(ALD)醛固酮				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pg/mL	118	94-142	
(ASD)雄烯二酮				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.92	0.74-1.1	
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.54	0.432-0.648	
(CEA)癌胚抗原				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.26	1-1.52	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.39	1.11-1.67	
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.06	0.84-1.28	
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.1	1.68-2.52	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.14	0.92-1.36	
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.11	0.89-1.33	
(Cortisol)皮质醇				
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	32.1	25.7-38.5	
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.3	19.5-29.1	
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	32.9	26.3-39.5	
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	27.2	21.8-32.6	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	30.1	24.1-36.1	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.8	21.4-32.2	
(C-P)C-肽				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.1	0.88-1.32	
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.97	0.77-1.17	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.91	0.73-1.09	
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.98	0.78-1.18	
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.73	0.59-0.87	

REF	501-1	6x5mL	EXP	2027-06-11	LOT	2FCA796	更新日期:	2024-09
								3 / 6

(cPSA)复合前列腺特异性抗原			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/L	270	210-330
(CRBM)卡马西平			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	µg/mL	2.72	2.18-3.26
(CT)降钙素			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/L	49.9	39.9-59.9
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/L	45.9	36.7-55.1
(DGXN)地高辛			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.74	0.592-0.888
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	ng/mL	0.68	0.54-0.82
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.74	0.592-0.888
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.85	0.68-1.02
(DHEA-S)硫酸脱氢表雄酮			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	720	576-864
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	715	572-858
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	685	548-822
(E2)雌二醇			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	352	282-422
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	281	225-337
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	272	218-326
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pmol/L	259	207-311
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	275	220-330
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	374	299-449
(FA)叶酸			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.71	2.97-4.45
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.63	3.71-5.55
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.86	3.08-4.64
(FE3/uE3)游离雌三醇			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.48	1.18-1.78
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.35	1.09-1.61
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.35	1.09-1.61
(Ferritin)铁蛋白			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	57.5	46-69
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	51.9	41.5-62.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	66.9	53.5-80.3
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	85.5	68.4-102.6
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	71.8	57.4-86.2
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	55.4	44.3-66.5
(F-PSA)游离前列腺特异性抗原			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.44	0.36-0.52
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.34	0.272-0.408
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.41	0.33-0.49
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.49	0.39-0.59
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.38	0.3-0.46

REF	501-1	6x5mL	EXP	2027-06-11	LOT	2FCA796	更新日期:	2024-09
								4 / 6

(FSH)促卵泡刺激素				
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L		6.72	5.38-8.06
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	IU/L		10.1	8.1-12.1
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L		7.42	5.94-8.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L		8.2	6.56-9.84
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L		8.37	6.7-10.04
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L		8.2	6.56-9.84
(FT3)游离T3				
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L		4.57	3.65-5.49
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L		4.05	3.25-4.85
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L		4.05	3.25-4.85
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pmol/L		4.14	3.31-4.97
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L		4.17	3.34-5
(FT4)游离T4				
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L		10.6	8.5-12.7
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L		12.4	9.8-15
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L		10.9	8.7-13.1
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L		10.9	8.7-13.1
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pmol/L		12.4	9.8-15
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L		10.2	8.1-12.3
(Gentamicin)庆大霉素				
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L		4870	3906-5834
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L		3825	3075-4575
(hCG)人绒毛膜促性腺激素				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	IU/L		6.05	4.83-7.27
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L		7.22	5.78-8.66
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L		6.24	5-7.48
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L		5.62	4.5-6.74
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L		5.46	4.36-6.56
(HGH)人生长激素				
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	µg/L		2.87	2.29-3.45
(IgA)免疫球蛋白 IgA				
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	mg/dL		198	158-238
(IgE)免疫球蛋白 IgE				
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL		348	278-418
(IgG)免疫球蛋白 IgG				
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	mg/dL		777	622-932
(IgM)免疫球蛋白 IgM				
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	mg/dL		75.2	60.2-90.2
(INS)胰岛素				
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L		18.8	15-22.6
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L		19.1	15.3-22.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L		14.7	11.8-17.6
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	mIU/L		19.2	15.4-23
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L		12.3	9.8-14.8

REF	501-1	6x5mL	EXP	2027-06-11	LOT	2FCA796	更新日期:	2024-09
								5 / 6

(LH)促黄体生成素			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.82	1.46-2.18
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.32	1.06-1.58
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.32	1.06-1.58
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.94	1.56-2.32
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.14	0.92-1.36
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	1.47	1.17-1.77
(Li)锂			
Abbott Alinity C Series Abbott \ ISE (direct)	mmol/L	0.39	0.31-0.47
(Phenobarbital)苯巴比妥			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	51293	41034-61552
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	30905	24724-37086
(PHYT)苯妥英			
Abbott Alinity C Series Abbott \ Enzyme immunoassay (EIA)	µg/mL	6	4.8-7.2
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	µg/mL	7.04	5.65-8.43
(PRL)催/泌乳素			
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.66	6.11-9.21
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.69	4.55-6.83
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.33	4.27-6.39
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.76	7.01-10.51
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.05	4.83-7.27
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.31	5.85-8.77
(PSA)前列腺特异性抗原			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.78	0.619-0.941
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.88	0.7-1.06
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.83	0.671-0.989
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.64	0.52-0.76
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.78	0.619-0.941
(TES)睾酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.001	0.001-0.001
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.52	1.22-1.82
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.76	1.41-2.11
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.39	1.11-1.67
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.57	0.456-0.684
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.22	0.98-1.46
(THEO)茶碱			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	µg/mL	4.25	3.39-5.11
Abbott Alinity C Series Abbott \ Enzyme immunoassay (EIA)	µg/mL	4.1	3.28-4.92
(Tobramycin)妥布霉素			
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L	1906	1521-2291
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	1435	1148-1722
(TSH)促甲状腺素			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.58	0.464-0.696
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.45	0.36-0.54
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.39	0.31-0.47

REF	501-1	6x5mL	EXP	2027-06-11	LOT	2FCA796	更新日期:	2024-09
								6 / 6

Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.51	0.41-0.61
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.63	0.504-0.756
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.46	0.368-0.552
(TT3)总三碘甲状腺氨酸			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.77	0.61-0.93
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.91	0.73-1.09
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.711	0.569-0.853
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.83	0.664-0.996
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.87	0.691-1.049
(TT4)总T4			
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	41.4	33.2-49.6
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	65.1	52.1-78.1
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	54.7	43.7-65.7
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	58.5	46.7-70.3
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	55.5	44.3-66.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	55.5	44.3-66.7
(VALP)丙戊酸			
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	µg/mL	33.6	26.9-40.3
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	µg/mL	36.9	29.5-44.3
(Vancomycin)万古霉素			
Abbott Alinity C Series Abbott \ Immunoturbidimetry (ITA)	nmol/L	5058	4046-6070
Siemens ADVIA Centaur CP Siemens \ Dipstick	nmol/L	3919	3135-4703
(VB12)维生素B12			
Siemens ADVIA Centaur CP Siemens \ Chemiluminescence Immunoassay (CLIA)	pg/mL	372	298-446
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	278	222-334
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	262	210-314
Abbott Alinity i Series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	345	276-414
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pg/mL	449	359-539

注: 此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来, 仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。