

<b>REF</b> QLA-0062 6x4.5mL	<b>EXP</b> 2024-12-05	<b>LOT</b> 2407A0062	更新日期： 2024-08 1 / 3
-----------------------------	-----------------------	----------------------	------------------------

本报告含以下项目的示范实验室数据

分析项目(Analyte)	分析项目(Analyte)	分析项目(Analyte)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：产品用作非定值质控品时，可用于检测分析物列表所列项目的各类仪器设备。

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

<b>REF</b> QLA-0062	6x4.5mL	<b>EXP</b> 2024-12-05	<b>LOT</b> 2407A0062	更新日期: 2024-08
				2 / 3

批号: 2407A0062				
项目\仪器\试剂方法	单位	均值	+ / - 2 SD	
<b>(BASO %)嗜碱性粒细胞 %</b>				
Mindray BS series		5	0-10	
Mindray \ Flow cytometry	%			
Sysmex XE series		50	30-70	
sysmex \ Calculate	%			
Sysmex XS series		10	0-20	
sysmex \ Calculate	%			
Sysmex XT series		55	25-85	
sysmex \ Calculate	%			
Sysmex XN series		5	0-10	
sysmex \ Calculate	%			
<b>(EOS %)嗜酸性粒细胞 %</b>				
Mindray BC-6000 series		10	0-20	
Mindray \ Flow cytometry	%			
Sysmex XE series		15	0-30	
sysmex \ Calculate	%			
Sysmex XS series		10	0-20	
sysmex \ Calculate	%			
Sysmex XT series		15	0-30	
sysmex \ Calculate	%			
Sysmex XN series		10	0-20	
sysmex \ Calculate	%			
<b>(HCT)血细胞比容/压积</b>				
Mindray BC-6000 series		35.7	27.9-43.5	
Mindray \ Light scattering	%			
Sysmex XE series		35.5	27.9-43.1	
sysmex \ Light scattering	%			
Sysmex XS series		35.7	27.9-43.5	
sysmex \ Light scattering	%			
Sysmex XT series		35.1	27.5-42.7	
sysmex \ Light scattering	%			
Sysmex XN series		35.5	27.9-43.1	
sysmex \ Light scattering	%			
<b>(HGB/ Hb)血红蛋白</b>				
Mindray BC-6000 series		133	107-159	
Mindray \ Colorimetry	g/L			
Sysmex XE series		129	103-155	
sysmex \ Light scattering	g/L			
Sysmex XS series		132	108-156	
sysmex \ Light scattering	g/L			
Sysmex XT series		134	108-160	
sysmex \ Light scattering	g/L			
Sysmex XN series		132	106-158	
sysmex \ Light scattering	g/L			
<b>(LYM %)淋巴细胞 %</b>				
Mindray BC-6000 series		35	20-50	
Mindray \ Flow cytometry	%			
Sysmex XE series		30	15-45	
sysmex \ Calculate	%			
Sysmex XS series		35	20-50	
sysmex \ Calculate	%			
Sysmex XT series		35	20-50	
sysmex \ Calculate	%			
Sysmex XN series		40	20-60	
sysmex \ Calculate	%			
<b>(MCH)平均红细胞血红蛋白含量</b>				
Mindray BC-6000 series		30.2	24.4-36	
Mindray \ Light scattering	pg			
Sysmex XE series		31	25.4-36.6	
sysmex \ Light scattering	pg			
Sysmex XS series		30	24.2-35.8	
sysmex \ Light scattering	pg			
Sysmex XT series		30	24.2-35.8	
sysmex \ Light scattering	pg			
Sysmex XN series		31.3	25.5-37.1	
sysmex \ Light scattering	pg			
<b>(MCHC)平均红细胞血红蛋白浓度</b>				
Mindray BC-6000 series		348	284-412	
Mindray \ Calculate	g/L			
Sysmex XE series		347	285-409	
sysmex \ Light scattering	g/L			
Sysmex XS series		343	279-407	
sysmex \ Light scattering	g/L			
Sysmex XT series		347	283-411	
sysmex \ Light scattering	g/L			
Sysmex XN series		344	282-406	
sysmex \ Light scattering	g/L			
<b>(MCV)平均红细胞体积</b>				
Mindray BC-6000 series		87.9	71.9-104	
Mindray \ Light scattering	fL			

注: 产品用作非定值质控品时, 可用于检测分析物列表所列项目的各类仪器设备。

注: 此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。

<b>REF</b> QLA-0062	6x4.5mL	<b>EXP</b> 2024-12-05	<b>LOT</b> 2407A0062	更新日期: 2024-08
				3 / 3

Sysmex XE series sysmex \ Light scattering	fL	86.9	70.9-103
Sysmex XS series sysmex \ Light scattering	fL	89	73-105
Sysmex XT series sysmex \ Light scattering	fL	86	70-102
Sysmex XN series sysmex \ Light scattering	fL	87.9	71.9-104
<b>(MONO %)单核细胞 %</b>			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	5	0-10
Sysmex XS series sysmex \ Calculate	%	5	0-10
Sysmex XT series sysmex \ Calculate	%	10	0-20
Sysmex XN series sysmex \ Calculate	%	5	0-10
<b>(NEUT %)嗜中性粒细胞 %</b>			
Mindray BC-6000 series Mindray \ Flow cytometry	%	50	30-70
Sysmex XE series sysmex \ Calculate	%	60	40-80
Sysmex XS series sysmex \ Calculate	%	55	35-75
Sysmex XT series sysmex \ Calculate	%	50	30-70
Sysmex XN series sysmex \ Calculate	%	55	35-75
<b>(PLT)血小板</b>			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E+9/L)	240	190-290
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E+9/L)	249	199-299
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E+9/L)	247	197-297
Sysmex XT series sysmex \ Light scattering	Giga/L (1E+9/L)	242	192-292
Sysmex XN series sysmex \ Light scattering	Giga/L (1E+9/L)	245	195-295
<b>(RBC)红细胞</b>			
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E+12/L)	4.35	3.43-5.27
Sysmex XE series sysmex \ Light scattering	T/L (1E+12/L)	4.41	3.49-5.33
Sysmex XS series sysmex \ Light scattering	T/L (1E+12/L)	4.41	3.53-5.29
Sysmex XT series sysmex \ Light scattering	T/L (1E+12/L)	4.36	3.44-5.28
Sysmex XN series sysmex \ Light scattering	T/L (1E+12/L)	4.4	3.48-5.32
<b>(WBC)白细胞</b>			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E+9/L)	7.3	5.84-8.76
Sysmex XE series sysmex \ Light scattering	Giga/L (1E+9/L)	7.33	5.89-8.77
Sysmex XS series sysmex \ Light scattering	Giga/L (1E+9/L)	7.24	5.82-8.66
Sysmex XT series sysmex \ Light scattering	Giga/L (1E+9/L)	7.32	5.88-8.76
Sysmex XN series sysmex \ Light scattering	Giga/L (1E+9/L)	7.26	5.8-8.72

注: 产品用作非定值质控品时, 可用于检测分析物列表所列项目的各类仪器设备。

注: 此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同, 或因制造商检测方法的改变, 均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求, 实验室须遵循相关技术规范确立自己的均值和可接受范围。