

REF	539-2	6×3mL	EXP	2023-01-30	LOT	336C9FD	更新日期：	2024-04
1 / 5								

批号：336C9FD					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.77	10.4%	<5	<100
(AFP)甲胎蛋白					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	121	6.5%	19	401
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	119	5.3%	11	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	139	4.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	132	9.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	9.4%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	10.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	126	9.5%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	41.2	10%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	59.3	9.9%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	30.2	9.3%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.3	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.8	10.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	28.7	2.3%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.7	8.8%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.5	0.4%	<5	<100
(ASD)雄烯二酮					
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.01	11.4%	<5	<100
(CA 125)糖类抗原 125					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	66.8	9.9%	18	374
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	86.3	6.2%	11	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	56.5	11.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	60.4	10.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	78.4	10.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	62.8	4.2%	<5	<100
(CA 15-3)糖类抗原 15-3					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	29.7	7.3%	9	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	49.6	12.1%	17	347
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	41.4	9.9%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	40.2	9.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.5	10.4%	<5	<100
(CA 19-9)糖类抗原 19-9					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	34.4	10.8%	11	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	48.4	9.9%	<5	<100

REF	539-2	6×3mL	EXP	2023-01-30	LOT	336C9FD	更新日期：	2024-04
2 / 5								

		批号：336C9FD			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	33.9	12.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	38.1	10.2%	<5	<100
(CEA)癌胚抗原					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.5	7.2%	17	384
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	17.5	10.9%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.87	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.7	11.2%	<5	<100
(Cortisol)皮质醇					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	239	10.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	225	9.7%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	242	10.3%	<5	<100
(C-P)C-肽					
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.62	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.33	11.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.92	9.9%	<5	<100
(E2)雌二醇					
Abbott ARCHITECT i series OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	994	2.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	935	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1098	9.8%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	976	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	896	2.1%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1015	10%	<5	<100
(FA)叶酸					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.22	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.46	10.7%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.67	9.9%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5	4.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	10.4%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.87	10.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.98	9.8%	<5	<100
(FSH)促卵泡刺激素					
Abbott ARCHITECT i series OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	33.1	1.6%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	32.5	7.1%	11	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	24.9	9.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.7	10.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	25.2	2%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	33.1	10.6%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2023-01-30	LOT	336C9FD	更新日期：	2024-04
3 / 5								

		批号：336C9FD			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.2	10.2%	<5	<100
(FT3)游离T3					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	11.4	4.8%	12	265
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	12.1	5.2%	11	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.6	11.3%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	9.8%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	13.8	9.4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.4	10.5%	<5	<100
(FT4)游离T4					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	27.5	5.2%	18	393
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	31	5.7%	11	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	33.6	10.7%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	35.9	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.6	9.8%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	31.3	9.9%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.5	11.1%	<5	<100
(HCY)同型半胱氨酸					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	21.9	11%	<5	<100
(INS)胰岛素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.7	10.4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	33.1	10.3%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.6	2.4%	<5	<100
(LH)促黄体生成素					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	20.6	5.6%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.7	10.1%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	19.8	9.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.7	10.7%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.1	10.5%	<5	<100
(PRL)催/ 泌乳素					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	17.7	7.9%	7	<100
Abbott ARCHITECT i series OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	16.5	4.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	9.5%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.1	3.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.5	9.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.2	9.9%	<5	<100
(PROG)孕酮					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.8	10.5%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.09	10%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2023-01-30	LOT	336C9FD	更新日期：	2024-04
4 / 5								

批号：336C9FD					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.94	2.6%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.2	11.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.7	10.5%	<5	<100
(PSA)前列腺特异性抗原					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.5	11.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.8	12%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.1	9.9%	<5	<100
(PTH)甲状旁腺素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.27	9.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.45	10.1%	<5	<100
(TES)睾酮					
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.534	10.3%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.43	9.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.516	10.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.536	10.2%	<5	<100
(Tg)甲状腺球蛋白					
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	10.4%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.44	10.5%	<5	<100
(THEO)茶碱					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	24.5	9.6%	<5	<100
(TSH)促甲状腺素					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	10.5	3%	<5	<100
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	11.1	10.8%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.74	10.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.86	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.88	10.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.85	2%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	9.63	4%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.3	10.7%	<5	<100
(TT3)总三碘甲状腺氨酸					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.34	10.6%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.18	9.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.61	8.9%	<5	<100
(TT4)总T4					
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	137	5.5%	12	249
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	151	6.1%	11	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	95.6	8.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.009	10.4%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	88.5	2.6%	<5	<100

REF	539-2	6×3mL	EXP	2023-01-30	LOT	336C9FD	更新日期：	2024-04
5 / 5								

		批号：336C9FD			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	103	9.8%	<5	<100
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	9%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	97.5	2.2%	<5	<100
(VB12)维生素B12					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	627	10.4%	<5	<100
(β2M)β-2-微球蛋白					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1430	11.2%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	129	10.9%	<5	<100
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.134	10.4%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。