

REF	671-2	6×1mL	EXP	2023-12-25	LOT	40033BA	更新日期：	2024-04
1 / 2								

批号：40033BA					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	Seconds	58.7	8.3%	88	1397
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	83.8	3.1%	8	214
BE Thrombolyzer X/XRC/XRM Dedicated \ Light scattering	Seconds	53.9	5%	5	<100
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	73.9	8.4%	29	401
Sysmex CS 5100 Siemens Dade Actin FSL(APTT) \ Light scattering	Seconds	78.7	5.4%	17	467
Sysmex CS 5100 sysmex \ Light scattering	Seconds	67.7	5.5%	16	296
Stago STA series Stago \ Magnetic Beads	Seconds	76.2	12%	<5	<100
(AT-III)抗凝血酶III					
BE Thrombolyzer X/XRC/XRM Dedicated \ Magnetic Beads	%	65.5	8.8%	5	<100
Sysmex CA 1500/6000/7000 Dedicated \ Nephelometry	%	42	13.4%	29	394
Sysmex CS 5100 sysmex \ Nephelometry	%	32.2	12.4%	16	283
Sysmex CS 5100 OTHER \ Nephelometry	%	40.3	13.3%	13	203
Stago STA series Stago \ Magnetic Beads	%	42.5	12%	<5	<100
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	%	34	20.4%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	3.53	12%	<5	<100
(FDP)纤维蛋白(原)降解产物					
BE Thrombolyzer X/XRC/XRM Dedicated \ Magnetic Beads	mg/L	5.27	24.1%	5	<100
Sysmex CA 1500/6000/7000 Dedicated \ Nephelometry	mg/L	10.5	6.1%	28	342
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	mg/L	5.06	9.3%	<5	<100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	8.38	12%	<5	<100
(FIB)纤维蛋白原					
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	g/L	1.91	7.6%	87	1493
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	g/L	2	5.8%	8	217
BE Thrombolyzer X/XRC/XRM Dedicated \ Light scattering	g/L	2.23	8%	5	<100
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	g/L	1.88	6.7%	25	368
Sysmex CS 5100 sysmex \ Light scattering	g/L	1.96	5.6%	17	431
Sysmex CS 5100 Dedicated \ Magnetic Beads	g/L	1.91	6.8%	15	441
Sysmex CS 2000i/2100i Dedicated \ Magnetic Beads	g/L	1.95	5.5%	10	226
Stago STA series Stago \ Magnetic Beads	g/L	2.11	12%	<5	<100
(PT INR)PT国际标准化比值					
Sysmex CA 1500/6000/7000 sysmex \ Calculate	INR	2.44	6%	8	216
BE Thrombolyzer X/XRC/XRM Dedicated \ Calculate	INR	2.54	6.6%	5	<100
BE Thrombolyzer X/XRC/XRM OTHER \ Calculate	INR	2.57	4.1%	45	814
Sysmex CS 2000i/2100i Dedicated \ Calculate	INR	2.39	6.1%	<5	<100
Sysmex CS 5100 Dedicated \ Calculate	INR	2.53	3.6%	<5	222
BE Thrombolyzer X/XRC/XRM Behnk Elektronik (BE) \ Calculate	INR	2.61	4.9%	34	488
Sysmex CS 5100 sysmex \ Calculate	INR	2.27	16.2%	33	734
Sysmex CA 1500/6000/7000 Dedicated \ Calculate	INR	2.73	15.6%	11	143
Stago STA series Stago \ Calculate	INR	2.53	12%	<5	<100
(PT)血浆凝血酶原时间					
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	Seconds	29.7	4.7%	86	1216
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	28.7	5.2%	8	214
BE Thrombolyzer X/XRC/XRM Dedicated \ Magnetic Beads	Seconds	28.6	6.5%	5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	25.1	16.6%	33	735
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	29.9	10%	29	402

REF	671-2	6×1mL	EXP	2023-12-25	LOT	40033BA	更新日期：	2024-04
2 / 2								

		批号：40033BA			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Stago STA series Stago \ Magnetic Beads	Seconds	27.4	12%	<5	<100
(TT)凝血酶时间					
BE Thrombolyzer X/XRC/XRM OTHER \ OTHER	Seconds	10.5	4.2%	86	1308
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	15.6	2.8%	8	214
BE Thrombolyzer X/XRC/XRM Dedicated \ Light scattering	Seconds	10.3	2.8%	5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	19.2	17.5%	33	762
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	19.1	15.1%	29	399
Sysmex CS 5100 Dedicated \ Magnetic Beads	Seconds	18.3	3.7%	14	442
Sysmex CS 2000i/2100i Dedicated \ Magnetic Beads	Seconds	17.5	5.6%	10	227
Stago STA series Stago \ Magnetic Beads	Seconds	14.8	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。