

REF	530-3	6×3mL	EXP	2023-05-28	LOT	T329343	更新日期：	2024-04
1 / 1								

批号：T329343					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CK-MB mass)肌酸激酶-MB质量					
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	ng/mL	55.7	11.7%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	32	7.5%	5	105
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	39.1	10%	<5	<100
(CRP)C反应蛋白					
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	mg/dL	1.1	10.7%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Roche ACCU-CHEK series Roche \ Nephelometry	mg/L	1.55	10%	<5	<100
(hs-CRP)超敏C反应蛋白					
Beckman AU series Beckman Coulter \ immunoturbidimetric	mg/L	7.89	10%	<5	<100
(Myo / Mb)肌红蛋白					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	738	2.5%	8	113
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	381	10%	<5	<100
Perkinelmer SuperFlex PerkinElmer \ Chemiluminescence Immunoassay (CLIA)	μg/L	125	9.8%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	pg/mL	2356	9%	8	133
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	pmol/L	277	9%	8	133
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1393	10%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	164	10%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.14	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。