

|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 510-1 | 6×3mL | EXP | 2024-06-18 | LOT | 30A6143 | 更新日期： | 2024-04 |
| 1 / 3 |       |       |     |            |     |         |       |         |

| 批号：30A6143                                                                               |       |         |       |          |      |
|------------------------------------------------------------------------------------------|-------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                              | 单位    | 均值 \ CV |       | Lab \ 点数 |      |
| (ADA)腺苷脱氨酶                                                                               |       |         |       |          |      |
| AB Sciex API Series<br>AllergScreen \ Peroxidase                                         | U/L   | 2       | 0%    | <5       | <100 |
| (AFP)甲胎蛋白                                                                                |       |         |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)  | ng/mL | 4.21    | 2.9%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | ng/mL | 4.05    | 3%    | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 3.31    | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                    | ng/mL | 4.04    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )           | ng/mL | 4.14    | 2.8%  | <5       | <100 |
| (Alb)白蛋白                                                                                 |       |         |       |          |      |
| Beckman AU series<br>Beckman Coulter \ Bromocresol green                                 | g/L   | 1.05    | 6.7%  | <5       | <100 |
| (ALP)碱性磷酸酶                                                                               |       |         |       |          |      |
| AB Sciex API Series<br>ABX \ dry chemistry                                               | U/L   | 11.2    | 0%    | <5       | <100 |
| (Anti -HAV IgM)抗-甲肝病毒 IgM                                                                |       |         |       |          |      |
| maglumi system<br>Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)                   | IU/mL | 5       | 0%    | <5       | <100 |
| (AST/ SGOT)天门冬氨酸氨基转移酶                                                                    |       |         |       |          |      |
| AB Sciex API Series<br>ABX \ dry chemistry                                               | U/L   | 11.2    | 0%    | <5       | <100 |
| (CA 125)糖类抗原 125                                                                         |       |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)             | IU/mL | 9.86    | 94.1% | 35       | 752  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)  | IU/mL | 23      | 1.6%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | IU/mL | 31.7    | 4.1%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 18.8    | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                    | IU/mL | 26.8    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )           | IU/mL | 34.3    | 3.3%  | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                          | IU/mL | 57.2    | 9.1%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)            | IU/mL | 27.6    | 0.1%  | <5       | <100 |
| Abbott Alinity i Series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )                | IU/mL | 70.9    | 0.2%  | <5       | <100 |
| (CA 15-3)糖类抗原 15-3                                                                       |       |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)             | IU/mL | 10.5    | 43.7% | 36       | 764  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)  | IU/mL | 15.8    | 4.5%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | IU/mL | 14.8    | 4.5%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                    | IU/mL | 15.8    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )           | IU/mL | 17.5    | 3.5%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 11.8    | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)            | IU/mL | 22.1    | 0.1%  | <5       | <100 |
| (CA 19-9)糖类抗原 19-9                                                                       |       |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)             | IU/mL | 9.54    | 42.8% | 36       | 761  |
| Mindray CL-2000i/1000i/6000i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA ) | IU/mL | 31.2    | 14.7% | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)  | IU/mL | 12.3    | 4.1%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | IU/mL | 28.4    | 5.6%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)            | IU/mL | 29.4    | 0.2%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )           | IU/mL | 26.9    | 6.6%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 11.6    | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                    | IU/mL | 14.1    | 12%   | <5       | <100 |
| (CA 50)糖类抗原 50                                                                           |       |         |       |          |      |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                          | IU/mL | 11.7    | 6.4%  | <5       | <100 |

质控品名： 肿瘤标志物液体非定值质控品

|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 510-1 | 6×3mL | EXP | 2024-06-18 | LOT | 30A6143 | 更新日期： | 2024-04 |
| 2 / 3 |       |       |     |            |     |         |       |         |

| 批号：30A6143                                                                              |        |         |       |          |      |
|-----------------------------------------------------------------------------------------|--------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位     | 均值 \ CV |       | Lab \ 点数 |      |
| (CA 72-4)糖类抗原 72-4                                                                      |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL  | 1.87    | 4.9%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL  | 2.03    | 4.2%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | IU/mL  | 2.51    | 12%   | <5       | <100 |
| (CEA)癌胚抗原                                                                               |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)            | ng/mL  | 5.43    | 4.5%  | 8        | 112  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 2.41    | 5.6%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 2.76    | 4.4%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 3.72    | 1%    | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )  | ng/mL  | 1.18    | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | ng/mL  | 2.46    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | ng/mL  | 2.79    | 4.1%  | <5       | <100 |
| (Cyfra 21-1)细胞角蛋白19片段                                                                   |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 2.02    | 5.8%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | ng/mL  | 2.38    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | ng/mL  | 1.29    | 5.7%  | <5       | <100 |
| (DGXN)地高辛                                                                               |        |         |       |          |      |
| ABX Pentra Series<br>ImmunoCAP \ Chemiluminescence enzyme immunoassay (CLEIA)           | ng/mL  | 1.56    | 0%    | <5       | <100 |
| (Ferritin)铁蛋白                                                                           |        |         |       |          |      |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 62.8    | 3.6%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 58.5    | 0%    | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | ng/mL  | 66.6    | 2.6%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Roche \ Immunonephelometry (INA)                          | ng/mL  | 67.8    | 12%   | <5       | <100 |
| (F-PSA)游离前列腺特异性抗原                                                                       |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)            | ng/mL  | 0.994   | 6.7%  | 8        | 112  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 0.898   | 2.1%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | ng/mL  | 0.96    | 12.1% | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )  | ng/mL  | 0.85    | 12%   | <5       | <100 |
| (GLU)葡萄糖                                                                                |        |         |       |          |      |
| AB Sciex API Series<br>ABX \ dry chemistry                                              | mmol/L | 62.2    | 0%    | <5       | <100 |
| (PSA)前列腺特异性抗原                                                                           |        |         |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 1.86    | 2.2%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 1.75    | 3.5%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )  | ng/mL  | 1.71    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | ng/mL  | 1.57    | 2.5%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | ng/mL  | 2.23    | 12%   | <5       | <100 |
| (SCC)鳞状细胞癌抗原                                                                            |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 0.991   | 27.6% | 9        | 167  |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL  | 1.39    | 5.5%  | <5       | <100 |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 0.95    | 3.2%  | <5       | <100 |
| Mindray CL-2000i/1000i/6000i series<br>Mindray \ Chemiluminescence Immunoassay (CLIA )  | ng/mL  | 1.12    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | ng/mL  | 1.04    | 7.6%  | <5       | <100 |
| (β-hCG)人绒毛膜促性腺激素β亚基                                                                     |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Chemiluminescence enzyme immunoassay (CLEIA)            | IU/mL  | 0.01    | 6.6%  | 8        | 112  |

|                |       |       |                |            |                |         |       |         |
|----------------|-------|-------|----------------|------------|----------------|---------|-------|---------|
| <div>REF</div> | 510-1 | 6×3mL | <div>EXP</div> | 2024-06-18 | <div>LOT</div> | 30A6143 | 更新日期： | 2024-04 |
| 3 / 3          |       |       |                |            |                |         |       |         |

|                                                                                         |       |            |       |          |      |
|-----------------------------------------------------------------------------------------|-------|------------|-------|----------|------|
|                                                                                         |       | 批号：30A6143 |       |          |      |
| 项目\仪器\试剂\方法                                                                             | 单位    | 均值 \ CV    |       | Lab \ 点数 |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL | 0.013      | 3.8%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Chemiluminescence Immunoassay (CLIA )                   | IU/mL | 0.013      | 11.9% | <5       | <100 |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。