

REF	671-1	6×1mL	EXP	2023-01-13	LOT	4000C47	更新日期：	2024-04
1 / 1								

		批号：4000C47			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CS-1300 sysmex \ Light scattering	Seconds	33.6	3.4%	8	213
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	30.6	4.4%	49	1307
Sekisui Coapresta 3000 SekiSui \ OTHER	Seconds	28.8	22.3%	<5	<100
Werfen ACL series Dedicated \ Light scattering	Seconds	35.5	6.1%	14	331
Werfen ACL series Werfen \ Light scattering	Seconds	29.6	9.8%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	36.5	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	32.6	9.8%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	31.8	10.1%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CS 5100 Dedicated \ Nephelometry	%	93.1	4.5%	37	1048
Sekisui Coapresta 3000 SekiSui \ Nephelometry	%	94.2	19.6%	<5	<100
Werfen ACL series Dedicated \ Nephelometry	%	99.9	3.8%	15	398
(D-D(DDU)) D-二聚体(DDU)					
Stago STA series SekiSui \ Nephelometry	mg/L	1.06	12%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CS 5100 Dedicated \ Nephelometry	mg/L	1.55	7.7%	<5	<100
Sysmex CA 1500/6000/7000 SekiSui \ Nephelometry	mg/L	1.58	8.9%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Stago STA series SekiSui \ Nephelometry	mg/L	4.07	12%	<5	<100
Werfen ACL series SekiSui \ Nephelometry	mg/L	3.59	8.9%	<5	<100
Sysmex CA 1500/6000/7000 SekiSui \ Nephelometry	mg/L	4.23	9%	<5	<100
(FIB)纤维蛋白原					
Werfen ACL series Dedicated \ Light scattering	g/L	3.6	3.7%	<5	<100
Sysmex CS 5100 Dedicated \ Light scattering	g/L	3.61	4.1%	10	298
Werfen ACL series Werfen \ Light scattering	g/L	3.87	8.3%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	3.88	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	3.35	8.4%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	g/L	3.16	9.2%	<5	<100
(PT INR)PT国际标准化比值					
Stago STA series Stago \ Calculate	INR	1.02	12%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	12.4	4.8%	32	871
Werfen ACL series Dedicated \ Light scattering	Seconds	12	2.2%	14	370
Werfen ACL series Werfen \ Light scattering	Seconds	11.6	9.7%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	12.8	9.4%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	13.2	12%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	13.1	9.2%	<5	<100
(TT)凝血酶时间					
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	17.4	9.8%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	16.5	12%	<5	<100
Werfen ACL series SUNBIO \ Light scattering	Seconds	14.6	8.2%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。