

REF	671-2	6×1mL	EXP	2023-01-13	LOT	4003357	更新日期：	2024-04
1 / 2								

		批号：4003357			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CS-1300	Seconds	84.8	11.2%	8	211
sysmex \ Light scattering					
Sekisui Coapresta 3000	Seconds	75.1	10.9%	8	118
SekiSui \ OTHER					
Sysmex CS 5100	Seconds	70.8	8.4%	49	1310
Dedicated \ Light scattering					
Werfen ACL series	Seconds	68.4	9.6%	<5	<100
Werfen \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	96.3	9.6%	<5	<100
Siemens \ Light scattering					
Stago STA series	Seconds	75	12%	<5	<100
Stago \ Magnetic Beads					
Sysmex CS 5100	Seconds	98.6	9.3%	<5	<100
Siemens \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	72.5	0%	<5	<100
SUNBIO \ Light scattering					
(AT-III)抗凝血酶III					
Sekisui Coapresta 3000	%	44	6.6%	5	<100
SekiSui \ Nephelometry					
Sysmex CS 5100	%	46.2	7.3%	37	1049
Dedicated \ Nephelometry					
Werfen ACL series	%	46.9	9.2%	<5	<100
Dedicated \ Nephelometry					
(D-D(DDU)) D-二聚体(DDU)					
Werfen ACL series	mg/L	0.689	10.4%	<5	<100
Dedicated \ Nephelometry					
Stago STA series	mg/L	2.68	12%	<5	<100
SekiSui \ Nephelometry					
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CS 5100	mg/L	2.88	9.7%	<5	<100
Dedicated \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	2.96	9.5%	<5	<100
SekiSui \ Nephelometry					
(FDP)纤维蛋白(原)降解产物					
Sekisui Coapresta 3000	mg/L	9.31	7.2%	7	<100
SekiSui \ OTHER					
Sysmex CS 5100	mg/L	10.8	9.5%	43	946
Dedicated \ Nephelometry					
Werfen ACL series	mg/L	4970	9.9%	<5	<100
SekiSui \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	5.89	10%	<5	<100
SekiSui \ Nephelometry					
Stago STA series	mg/L	7.24	12%	<5	<100
SekiSui \ Nephelometry					
(FIB)纤维蛋白原					
Sysmex CA 1500/6000/7000	g/L	1.96	9.2%	<5	<100
Siemens \ Light scattering					
Stago STA series	g/L	2.54	12%	<5	<100
Stago \ Magnetic Beads					
Sysmex CS 5100	g/L	1.87	9.1%	<5	<100
Siemens \ Light scattering					
Werfen ACL series	g/L	2.21	8.6%	<5	<100
Werfen \ Light scattering					
(PT INR)PT国际标准化比值					
Sysmex CS 5100	INR	2.92	5.4%	<5	<100
Dedicated \ Calculate					
Stago STA series	INR	2.9	12%	<5	<100
Stago \ Calculate					
(PT)血浆凝血酶原时间					
Sysmex CS 5100	Seconds	30.6	5.1%	45	1218
Dedicated \ Light scattering					
Sysmex CS-1300	Seconds	26.9	5.4%	<5	<100
sysmex \ Light scattering					
Werfen ACL series	Seconds	30.6	9.5%	<5	<100
Werfen \ Light scattering					
Sysmex CS 5100	Seconds	32.3	9%	<5	<100
Siemens \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	32.1	9.7%	<5	<100
Siemens \ Light scattering					
Stago STA series	Seconds	29.6	12%	<5	<100
Stago \ Magnetic Beads					
Sysmex CA 1500/6000/7000	Seconds	35.2	0%	<5	<100
SUNBIO \ Light scattering					
(TT)凝血酶时间					
Sysmex CS 5100	Seconds	17.8	3.4%	28	699
Dedicated \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	17.6	9.7%	<5	<100
Dedicated \ Light scattering					
Stago STA series	Seconds	14.6	12%	<5	<100
Stago \ Magnetic Beads					
Werfen ACL series	Seconds	17.2	9.9%	<5	<100
SUNBIO \ Light scattering					

报告类别： 汇总实验室数据报告



质控品名： 凝血七项非定值质控品

REF	671-2	6×1mL	EXP	2023-01-13	LOT	4003357	更新日期：	2024-04
2 / 2								

		批号：4003357	
项目\仪器\试剂\方法	单位	均值 \ CV	Lab \ 点数

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。