

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
1 / 12								

批号：33828C6					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(17-OHP)17-羟孕酮					
Agilent agilent \ MS	ng/mL	0.759	11.9%	<5	<100
(25-OH-VD2)25羟基维生素D2					
OTHER Dedicated \ OTHER	ng/mL	0.03	3.7%	<5	<100
(25-OH-VD3)25羟基维生素D3					
OTHER Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.1	4.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20	11.9%	<5	<100
(25-OH-Vit D)25羟基维生素D					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.7	17.1%	20	452
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	61.7	17.1%	20	452
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.8	11.9%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	54.4	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20.7	9.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	51.6	9.8%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.9	10.3%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	64.8	10.3%	<5	<100
(AFP)甲胎蛋白					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.84	5.1%	9	213
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.32	2.8%	7	177
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.44	8%	68	1297
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.96	4.4%	5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.84	4.6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.17	5.6%	36	928
OTHER Abbott \ OTHER	ng/mL	7.25	8%	30	584
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.41	12%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.6	5%	<5	<100
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.69	5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.34	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.26	3%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.58	6.5%	18	468
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.51	12%	<5	<100
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	134	9.9%	<5	<100
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	pg/mL	70.3	11.9%	<5	100
(Anti-Tg)抗甲状腺球蛋白抗体					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	11.2	5.3%	8	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	19.1	37.1%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	21	9.1%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	10.3	21.2%	29	591
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	1.68	14%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	3.55	11.6%	18	480
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	17.4	1%	17	459

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
2 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	9.91	9.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	1.45	8.3%	<5	<100
OTHER Abbott \ OTHER	IU/mL	6.49	0%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	12.8	8.2%	9	113
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	7.95	7.8%	8	243
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15.7	16.3%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14.6	10.2%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	15.4	17.6%	<5	128
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	11.2	18.9%	29	590
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	4.78	7.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	6.7	9.1%	18	470
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	9.2	2.7%	17	459
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	8.31	9.9%	<5	<100
(BAP (Ost))骨碱性磷酸酶 (骨泌素)					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.1	6.7%	5	128
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.9	2.2%	<5	<100
(CA 125)糖类抗原 125					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.3	3.2%	5	<100
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	37.2	3.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	36.9	5.3%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	32	9.1%	36	941
OTHER Abbott \ OTHER	IU/mL	39.2	7.5%	30	593
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	34.8	2.8%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	34.2	7.7%	24	469
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.3	17%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	48.5	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	29.8	1.8%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	34.9	5%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	43.1	12%	<5	<100
(CA 15-3)糖类抗原 15-3					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.1	4.8%	<5	<100
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	28.8	6.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	14	9.5%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.7	4%	36	942
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	30.6	7.8%	24	449
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.4	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	18.7	17.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	2.9%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	26.7	4.9%	<5	<100
(CA 19-9)糖类抗原 19-9					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	26.7	8.5%	7	103

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
3 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	35.4	4%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	8%	36	932
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	24.2	3.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	8.5%	24	469
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	20.6	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	21.5	21.5%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	28.2	2.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	40.8	5.7%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	33.6	5.1%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.3	12%	<5	<100
(CEA)癌胚抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.32	8.9%	66	1326
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.22	4.1%	5	<100
OTHER Abbott \ OTHER	ng/mL	2.94	14.7%	30	593
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.63	11.1%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.86	8.9%	29	697
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.9	11.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	15.1%	15	385
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.05	1.3%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.99	12%	<5	<100
(CK)肌酸激酶					
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	130	10%	<5	<100
(CK-MB mass)肌酸激酶-MB质量					
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.91	9.3%	10	144
(CK-MB)肌酸激酶-MB (活性)					
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.93	6.6%	9	252
Beckman AU series Dedicated \ Immune suppression	U/L	9.5	9.1%	7	123
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.1	14.9%	44	871
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.79	26.5%	<5	<100
TOSHIBA TBA MAKER \ Immune suppression	U/L	8.51	7.2%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.49	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.25	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	7.09	8.2%	18	476
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.87	9.7%	18	352
(Cortisol)皮质醇					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	102	5.5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	92.4	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	115	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	134	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	94.3	8.3%	18	475
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	99	0.7%	17	458

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
4 / 12								

批号：33828C6					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(C-P)C-肽					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.9	12.2%	7	104
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.903	3.8%	6	105
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.814	4.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.06	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.648	11.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.409	14%	19	482
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.997	2.8%	17	456
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.956	3.4%	10	266
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.93	9.7%	<5	<100
(DGXN)地高辛					
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.932	15.6%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.88	10%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.87	11.9%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3364	3.7%	<5	115
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3390	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3416	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3906	10.4%	<5	<100
(E2)雌二醇					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	342	8.5%	9	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	292	4.7%	6	355
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	268	12.4%	6	102
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	288	5.1%	<5	123
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	346	10.7%	<5	103
OTHER Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	343	5.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	335	3.3%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	342	4.2%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	285	9.2%	29	594
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	331	2.4%	24	491
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	263	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	246	11.9%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	312	4.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	298	6.9%	18	468
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	275	11.4%	15	327
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	6.77	17.3%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	6.14	6.1%	18	470
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	6.76	9.9%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
5 / 12								

		批号：33828C6			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(FA)叶酸					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.52	23.9%	5	107
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.79	8.6%	<5	<100
maglumi system Snibe \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.81	1.8%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.47	4.9%	<5	123
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.66	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.87	6.7%	19	531
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.04	8.2%	18	479
(Ferritin)铁蛋白					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	76.9	8.8%	6	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	66.5	7.3%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	73.9	13.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	80.5	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	119	9.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	106	6.7%	19	468
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	72.8	7.4%	18	474
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	105	1.1%	17	458
DiaSorin Liaison /XL Diasorin \ Chemiluminescence Immunoassay (CLIA)	ng/mL	95.9	12%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	69.2	9.8%	<5	<100
(F-PSA)游离前列腺特异性抗原					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.97	4.1%	6	102
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.75	5.5%	41	924
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.16	4.8%	29	599
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.06	11.9%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.48	6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.71	4.9%	19	474
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.93	6.5%	18	474
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.81	12%	<5	<100
OTHER Abbott \ OTHER	ng/mL	2.05	3.1%	<5	<100
(FSH)促卵泡刺激素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.37	14.3%	8	109
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.89	3%	6	184
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.37	10.6%	6	100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.47	10.4%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.66	1.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.66	6.8%	<5	123
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.85	4.6%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.21	5.1%	29	600
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.92	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.8	11.9%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.4	3.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	6.34	7.3%	18	471

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
6 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.3	3.1%	18	462
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.77	4.8%	12	284
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.51	12%	<5	<100
(FT3)游离T3					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.81	6.1%	8	243
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.92	7.1%	6	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.1	5.9%	5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.13	3.4%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.66	2.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.07	8.3%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.36	7.1%	23	462
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.43	5.7%	21	1001
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.5	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.18	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	4.95	3.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	4.78	7%	18	949
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	6.72	3%	17	458
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.99	12%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	6.51	2%	<5	<100
(FT4)游离T4					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.8	6.6%	8	243
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.8	9%	6	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.7	4.1%	5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.4	5.1%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.1	4.3%	<5	103
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.3	9.1%	<5	100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.7	5.4%	23	468
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.4	5.8%	21	986
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	14.1	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.8	2.5%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	11.6	7.1%	18	944
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	18.5	1%	17	458
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	12%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	16	3.7%	<5	<100
(hCG)人绒毛膜促性腺激素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	0.031	3.8%	<5	<100
OTHER Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.62	3.8%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.42	9.6%	23	457
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.22	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	0.029	7.9%	18	481
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.74	45.4%	14	339

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
7 / 12								

批号：33828C6					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(HCY)同型半胱氨酸					
TOSHIBA TBA Bio-assay \ Enzymatic Cycling Assays Method	μmol/L	11.7	3.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	μmol/L	10.8	6.6%	<5	<100
Beckman AU series Dedicated \ Enzymatic Cycling Assays Method	μmol/L	12.1	15.8%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.9	9.9%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	12.1	9.9%	<5	<100
(HGH)人生长激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	1.51	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	1.99	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	1.41	5.5%	18	471
(HS-CTnl) (高) 超敏肌钙蛋白I					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.061	11.8%	<5	<100
(HS-CTnT)超敏肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.019	9.3%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	131	10%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	160	11.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	161	8%	19	440
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	218	9.5%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	17.7	11.9%	<5	<100
(INS)胰岛素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	7.34	12.5%	8	107
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.39	3.6%	6	104
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	6.72	5.7%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.33	5.1%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.98	2.8%	<5	117
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.01	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	6.95	11.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.36	4.7%	19	484
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	8.15	0.8%	17	457
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	12	4.7%	10	265
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.83	10%	<5	<100
(IntactPTH)全段甲状旁腺激素					
Abbott ARCHITECT i series Abbott \ OTHER	pg/mL	2.98	11.6%	6	<100
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	3.62	10.3%	<5	<100
(LH)促黄体生成素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.51	6.9%	7	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.4	4.4%	6	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.38	3.6%	6	188
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.62	10.1%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.81	5.4%	<5	124
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.41	4.9%	<5	<100

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
8 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.54	8.6%	29	600
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.47	1.6%	24	495
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.62	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.79	11.9%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.92	6.1%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.43	2.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	3.76	8.4%	18	466
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	5.17	5.4%	12	288
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.93	12%	<5	<100
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.66	11.6%	<5	<100
(Myo / Mb)肌红蛋白					
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	40.7	3.5%	9	256
TOSHIBA TBA OTHER \ OTHER	μg/L	42.8	15.1%	5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	28.9	6.3%	44	877
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	26.9	8.9%	<5	103
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	29.1	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	34.2	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	28.5	7%	18	473
J&J Vitros ECL Ortho \ Chemiluminescence Immunoassay (CLIA)	μg/L	39.3	6.6%	18	319
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	41.8	5.8%	10	146
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	182	5.1%	42	843
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	21.4	5.1%	42	843
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	149	4.5%	<5	<100
Ortho Vitros 3600/5600 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.5	4.5%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	202	9.9%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	23.8	9.9%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	184	5.1%	11	249
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	21.6	5.1%	11	249
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pg/mL	234	2.8%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	pmol/L	27.5	2.8%	<5	<100
J&J Vitros ECL Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	133	0%	<5	<100
J&J Vitros ECL Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.6	0%	<5	<100
(PCT)降钙素原					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.463	6.7%	10	227
(PRL)催/ 泌乳素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.15	6%	6	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.72	2.3%	6	187
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.34	11.4%	6	101
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.25	2.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.89	5.4%	<5	123

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
9 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	5%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.42	5%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.69	4.9%	29	587
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.17	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.01	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.83	9.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.3	5.9%	18	472
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.24	0.8%	17	457
Siemens Immulite1000/2000 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.64	5.2%	12	284
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.12	12%	<5	<100
(PROG)孕酮					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.983	32%	8	134
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.04	4.1%	7	140
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.94	6.2%	6	302
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.848	7.1%	5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.023	6.8%	5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.11	14.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.916	31.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.27	10.6%	<5	126
OTHER Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.15	8.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.45	9.3%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.29	10.2%	<5	<100
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.869	16.7%	28	602
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.22	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.01	11.9%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.21	18.9%	18	464
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.78	2.1%	17	457
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.17	8.3%	12	283
(PSA)前列腺特异性抗原					
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.17	11.2%	7	<100
Autolumo A2000/A2000 PLUS Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.99	3.1%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.55	11.6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.62	2.7%	20	485
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.08	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.16	1.4%	<5	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.52	2.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.22	6.1%	19	519
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.05	6.9%	18	479
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.71	10.4%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.33	12%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
10 / 12								

批号：33828C6					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
OTHER Abbott \ OTHER	ng/mL	4.16	6.8%	<5	<100
(PTH)甲状旁腺素					
Siemens ADVIA Centaur XP/XPT	pmol/L	1.83	13.3%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	pmol/L	1.91	6.7%	<5	<100
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Tosoh AIA 900/1800	pmol/L	1.94	11.9%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	pmol/L	2.05	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	pmol/L	2	8.4%	18	472
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Roche cobas 4000/6000/8000	pmol/L	1.95	1.1%	17	459
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	pmol/L	1.97	9.4%	12	270
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	pmol/L	2.26	9.8%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	1.96	9.7%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
OTHER	pmol/L	1.88	7.7%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
(SHBG)性激素结合球蛋白					
Beckman UniCel Dxl 600/800	nmol/L	33.2	12%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
(TES)睾酮					
Abbott ARCHITECT i series	ng/mL	6.39	9.5%	8	117
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	ng/mL	6.71	2.4%	6	230
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e411	ng/mL	6.02	3.9%	5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	ng/mL	6.27	4.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	ng/mL	5.53	3.7%	<5	<100
Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman AU series	ng/mL	5.29	4.6%	<5	<100
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
OTHER	ng/mL	6.28	7.4%	29	596
ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	ng/mL	5.13	5%	22	594
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	5.28	11.9%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	6.17	10.2%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Tosoh AIA 900/1800	ng/mL	6.61	11.9%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	6.13	0.9%	17	457
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens Immulite1000/2000	ng/mL	6.03	7%	12	285
Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)					
(Tg)甲状腺球蛋白					
Beckman UniCel Dxl 600/800	ng/mL	9.97	3.8%	8	242
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	10.7	9.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	ng/mL	8.26	11.3%	18	474
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Abbott ARCHITECT i series	ng/mL	9.88	9.9%	<5	<100
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	8.04	10.2%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
(THEO)茶碱					
Agilent agilent \ MS	µg/mL	8.7	11.9%	<5	<100
(TnI)肌钙蛋白I					
Ortho Vitros 3600/5600	ng/mL	0.079	8.5%	9	136
Dedicated \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	ng/mL	0.066	35.8%	<5	143
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
J&J Vitros ECL	ng/mL	0.075	8.1%	18	337
Ortho \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	ng/mL	0.07	9.7%	10	276
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
11 / 12								

		批号：33828C6			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(TnT)肌钙蛋白T					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02	3.7%	8	165
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02	4.4%	42	854
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.02	10.2%	<5	<100
(TSH)促甲状腺素					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.76	9.3%	9	138
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.12	4%	8	242
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.34	5.2%	6	191
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.5	3.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3	3.4%	<5	<100
OTHER Abbott \ OTHER	mIU/L	2.47	6%	30	596
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.48	4.8%	30	562
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.96	5%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.39	4.9%	21	983
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.9	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.87	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.32	1.3%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.83	5.2%	18	949
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.34	2.1%	17	458
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.63	7.6%	13	272
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3	12%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.86	4.4%	<5	<100
(TT3)总三碘甲状腺氨酸					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.13	6.7%	8	243
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.17	7.9%	5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.45	7.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.27	9.5%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.58	4%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.2	5.4%	21	561
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.14	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.23	3.6%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.07	9.5%	18	476
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.18	17%	10	123
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.956	12%	<5	<100
OTHER Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.32	5.5%	<5	<100
(TT4)总T4					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	64.2	5%	<5	<100
OTHER Abbott \ OTHER	ng/mL	73.1	6.7%	30	588
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	93.5	15.2%	30	582
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	75.4	7.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	69.2	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	74.8	2.7%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	76.7	8.9%	18	476

REF	540-1	6×3mL	EXP	2025-02-16	LOT	33828C6	更新日期：	2024-04
12 / 12								

		批号：33828C6			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	59.7	5.5%	17	436
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	58.1	10.2%	10	139
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	55.4	12%	<5	<100
(VB12)维生素B12					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	347	9.9%	5	112
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	485	5.9%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pg/mL	610	11.4%	<5	<100
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	389	5.1%	<5	123
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pg/mL	312	21.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	324	11.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	509	11.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	528	4.8%	19	514
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	383	6.5%	18	474
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	1020	9.9%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1170	11.9%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.006	7.7%	9	132
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.006	4%	7	136
OTHER ABX \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	8.4%	6	132
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.003	10.3%	5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.004	15.1%	<5	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.004	4.4%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.003	3.4%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.003	4.6%	23	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.003	11.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	2.6%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.005	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。