

REF	600-1	30×1.7mL	EXP	2024-03-10	LOT	393B5E3	更新日期：	2024-04
1 / 1								

		批号: 393B5E3			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Ca)钙					
Radiometer ABL Series Dedicated \ ISE (direct)	mmol/L	0.82	2.7%	6	131
OTHER OTHER \ OTHER	mmol/L	1.31	0.8%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ ISE (direct)	mmol/L	1.47	8.8%	<5	<100
Radiometer ABL Series Rayto \ ISE (indirect)	mmol/L	1.48	8.8%	<5	<100
Roche cobas b 121/123/221 Dedicated \ ISE (indirect)	mmol/L	1.44	9%	<5	<100
(Cl)氯					
Radiometer ABL Series Rayto \ ISE (direct)	mmol/L	75.8	9.8%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ ISE (direct)	mmol/L	82.3	9%	<5	<100
IL Gem Premier 3000 series Dedicated \ ISE (direct)	mmol/L	81.2	9.5%	<5	<100
(K)钾					
OTHER OTHER \ OTHER	mmol/L	2.15	2.4%	<5	<100
Siemens RAPIDPoint 400/500 Siemens \ ISE (direct)	mmol/L	2.35	8.9%	<5	<100
IL Gem Premier 3000 series Werfen \ ISE (direct)	mmol/L	2.29	9.6%	<5	<100
Radiometer ABL Series Rayto \ ISE (direct)	mmol/L	2.36	8.9%	<5	<100
(Na)钠					
OTHER OTHER \ OTHER	mmol/L	115	0.6%	<5	<100
IL Gem Premier 3000 series Dedicated \ ICP-MS	mmol/L	109	8.3%	<5	<100
Radiometer AQTSeries Dedicated \ ISE (direct)	mmol/L	113	8%	<5	<100
Radiometer ABL Series Rayto \ ISE (direct)	mmol/L	112	8%	<5	<100
(PCO2)二氧化碳分压					
Radiometer ABL Series Radiometer \ Electrode	mmHg	62	5.5%	9	109
OTHER OTHER \ Electrode	mmHg	64.7	1.4%	<5	<100
Radiometer AQTSeries Rayto \ Electrode	mmHg	62.1	9.5%	<5	<100
Roche cobas b 121/123/221 Roche \ Electrode	mmHg	62.8	8.1%	<5	<100
Nova Stat Profile pHox series Dedicated \ Electrode	mmHg	59.6	9.8%	<5	<100
Radiometer ABL Series Rayto \ Electrode	mmHg	63.8	9.2%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	61.1	9.6%	<5	<100
(pH)酸碱度					
Radiometer ABL Series Dedicated \ Electrode	No Units	7.63	0.1%	6	144
OTHER OTHER \ Electrode	No Units	7.14	0.2%	<5	<100
Radiometer ABL Series Radiometer \ Electrode	No Units	7.19	0.2%	10	133
Roche cobas b 121/123/221 Roche \ Electrode	No Units	7.31	9.4%	<5	<100
Radiometer ABL Series Rayto \ Electrode	No Units	7.22	9.6%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	No Units	7.32	9.7%	<5	<100
Radiometer AQTSeries Rayto \ Electrode	No Units	7.39	9.3%	<5	<100
(PO2)氧分压					
Radiometer ABL Series Radiometer \ Electrode	mmHg	66.5	6.3%	9	112
OTHER OTHER \ Electrode	mmHg	79.5	1.2%	<5	<100
Nova Stat Profile pHox series Dedicated \ Electrode	mmHg	86.5	7.8%	<5	<100
IL Gem Premier 3000 series Dedicated \ Electrode	mmHg	73.6	9.8%	<5	<100
Radiometer ABL Series Rayto \ Electrode	mmHg	71.5	9.9%	<5	<100
Roche cobas b 121/123/221 Roche \ Electrode	mmHg	72.8	9.3%	<5	<100
Siemens RAPIDPoint 400/500 Dedicated \ Electrode	mmHg	84.2	8%	<5	<100
Radiometer AQTSeries Dedicated \ Electrode	mmHg	69.7	9.8%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。