

质控品名： 多专业免疫分析非定值质控品

REF	539-1	6×3mL	EXP	2023-11-05	LOT	336A1CB	更新日期：	2024-05
1 / 3								

批号：336A1CB					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(AFP)甲胎蛋白					
Roche cobas E601/602	ng/mL	7.89	6.7%	21	602
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	ng/mL	6.76	4.4%	13	326
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS	pg/mL	161	11.2%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(Anti-Tg)抗甲状腺球蛋白抗体					
Roche cobas E601/602	IU/mL	17	14.7%	21	573
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	IU/mL	6.79	6.9%	13	378
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	0.75	12%	<5	<100
Mindray \ Chemiluminescence Immunoassay (CLIA)					
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Roche cobas E601/602	IU/mL	11.4	18.4%	21	552
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	IU/mL	8.56	5.3%	13	379
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	9.05	12%	<5	<100
Mindray \ Chemiluminescence Immunoassay (CLIA)					
(CA 125)糖类抗原 125					
Roche cobas E601/602	IU/mL	9.43	6.3%	21	595
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	IU/mL	10.4	4.3%	13	315
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
(CA 15-3)糖类抗原 15-3					
Roche cobas E601/602	IU/mL	12.1	6.8%	21	589
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	IU/mL	11.9	12.7%	13	289
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
(CA 19-9)糖类抗原 19-9					
Roche cobas E601/602	IU/mL	5.2	6.1%	21	621
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(CEA)癌胚抗原					
Roche cobas E601/602	ng/mL	3.18	5.2%	21	616
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	ng/mL	3.13	5.3%	13	343
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
(Cortisol)皮质醇					
Roche cobas E601/602	ng/mL	109	16.6%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(C-P)C-肽					
Roche cobas E601/602	ng/mL	1.09	6.5%	21	204
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	ng/mL	0.778	4.5%	12	351
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(DHEA-S)硫酸脱氢表雄酮					
Roche cobas E601/602	ng/mL	3130	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(E2)雌二醇					
Roche cobas E601/602	pmol/L	238	7.8%	21	224
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	pmol/L	231	6.4%	13	374
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	239	4.9%	11	<100
OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)					
(F-PSA)游离前列腺特异性抗原					
Roche cobas E601/602	ng/mL	0.116	7.1%	21	223
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(FSH)促卵泡刺激素					
Abbott ARCHITECT i series	IU/L	5.02	2.3%	8	<100
OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)					
Roche cobas E601/602	IU/L	4.91	6.4%	21	206
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	IU/L	4.94	3.1%	13	378
Abbott \ Chemiluminescence Immunoassay (CLIA)					
(FT3)游离T3					
Roche cobas E601/602	pmol/L	5.93	5%	21	591
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Abbott ARCHITECT i series	pmol/L	4.85	4.6%	13	380
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	pmol/L	4.72	4.7%	13	325
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					

质控品名： 多专业免疫分析非定值质控品

REF	539-1	6×3mL	EXP	2023-11-05	LOT	336A1CB	更新日期：	2024-05
2 / 3								

		批号：336A1CB			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.57	12%	<5	<100
(FT4)游离T4					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.8	4.7%	21	589
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.7	3.1%	13	375
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	12.7	3.2%	13	345
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.9	12%	<5	<100
(hCG)人绒毛膜促性腺激素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	3.5	6.4%	20	439
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.6	6.7%	13	383
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	4.14	7.3%	13	297
Ortho Vitros 3600/5600 Ortho \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.64	12.1%	<5	<100
(HCY)同型半胱氨酸					
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.6	10%	<5	<100
(HGH)人生长激素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	2.06	10%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Immage 800 Roche \ Immunoturbidimetry (ITA)	IU/mL	90.9	10%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	215	10%	<5	<100
(INS)胰岛素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.15	6.8%	21	221
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	6.61	3.2%	13	381
(LH)促黄体生成素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.54	6.3%	21	214
(PRL)催/ 泌乳素					
Abbott ARCHITECT i series OTHER \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	8.64	2.4%	7	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.27	5.9%	21	203
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.22	2.6%	13	372
(PROG)孕酮					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.785	6.7%	21	477
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.636	5.5%	13	376
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.738	11.9%	13	286
(PSA)前列腺特异性抗原					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.606	7.1%	21	209
(PTH)甲状旁腺素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.02	13.4%	21	199
(TES)睾酮					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.46	7.5%	21	226
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.75	4.6%	13	375
(Tg)甲状腺球蛋白					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.2	5.6%	19	231
(TSH)促甲状腺素					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	1.94	4.2%	21	623
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	1.58	3.3%	13	377
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	1.51	4%	13	338

REF	539-1	6×3mL	EXP	2023-11-05	LOT	336A1CB	更新日期：	2024-05
3 / 3								

批号：336A1CB					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.25	12%	<5	<100
(TT3)总三碘甲状腺氨酸					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.25	5.1%	21	235
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.02	3.5%	13	377
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.972	4.9%	13	342
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.97	12%	<5	<100
(TT4)总T4					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	64.1	5.8%	21	245
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	58.3	3.2%	13	377
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	58.1	4.5%	13	340
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	63.5	12%	<5	<100
(VB12)维生素B12					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	391	9.8%	21	210
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	855	10%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.004	10%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。