

REF	521-2	6×3mL	EXP	2026-05-22	LOT	31B52C1	更新日期：	2024-04
1 / 3								

批号：31B52C1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Olympus AU400 Olympus \ Bromocresol green	g/L	0.038	10.1%	<5	<100
(ASO)抗链球菌溶血素 O 抗体					
Abbott ARCHITECT c series BSBE \ OTHER	IU/mL	219	1.7%	5	102
Mindray BS series MedicalSystem \ OTHER	IU/mL	288	1.5%	<5	<100
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	213	2.1%	<5	<100
Beckman AU series BSBE \ OTHER	IU/mL	238	22%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	272	10%	<5	<100
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	242	1.3%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	166	4.8%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	IU/mL	335	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	IU/mL	265	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	251	0.4%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	IU/mL	264	10%	<5	<100
(AT-III)抗凝血酶III					
Beckman Image 800 Beckman Coulter \ Nephelometry	g/L	0.236	10%	<5	<100
(C3)补体 C3					
Mindray BS series MedicalSystem \ OTHER	mg/dL	129	1.8%	<5	<100
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	152	4.3%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	144	10%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	144	10%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	142	10%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	136	3.5%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	1.47	4.2%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	134	10%	<5	<100
(C4)补体 C4					
Mindray BS series MedicalSystem \ OTHER	mg/dL	32.5	1.3%	<5	<100
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	39.3	4.7%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	39.8	4.3%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	39	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	0.342	3.9%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	33	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	39.3	10%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	34	10%	<5	<100
(Cer)铜蓝蛋白					
Olympus AU400 Olympus \ Immunonephelometry (INA)	g/L	0.483	10%	<5	<100
(CRP)C反应蛋白					
Mindray BS series MedicalSystem \ OTHER	mg/dL	2.83	1.4%	<5	<100
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	2.8	3%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	3.01	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	2.98	9.7%	<5	<100
TOSHIBA TBA TOSHIBA \ Immunoturbidimetry (ITA)	mg/dL	2.2	6.1%	<5	<100
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	2.76	0.4%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	2.88	1.8%	<5	<100

REF	521-2	6×3mL	EXP	2026-05-22	LOT	31B52C1	更新日期：	2024-04
2 / 3								

批号：31B52C1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Ferritin)铁蛋白					
Siemens BN II Siemens \ Immunoturbidimetry (ITA)	ng/mL	116	10%	<5	<100
Olympus AU400 Olympus \ Immunoturbidimetry (ITA)	ng/mL	135	10%	<5	<100
(Hap (HPT))触珠蛋白 (结合珠蛋白)					
Abbott ARCHITECT c series chongqiong zybio \ OTHER	mg/L	1321	5%	9	223
Beckman AU series chongqiong zybio \ OTHER	mg/L	1428	10.9%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/L	1350	10%	<5	<100
(IgA)免疫球蛋白 IgA					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	326	2.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	262	2.5%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	293	10%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	292	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	271	10%	<5	<100
Abbott ARCHITECT c series DiaSys \ Immunoturbidimetry (ITA)	mg/dL	267	10%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	326	10%	<5	<100
(IgE)免疫球蛋白 IgE					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	48.5	2.6%	<5	<100
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	48.8	4%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	47.1	10%	<5	<100
(IgG)免疫球蛋白 IgG					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	1564	4.5%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	1607	2%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1471	10%	<5	<100
Abbott ARCHITECT c series DiaSys \ Immunoturbidimetry (ITA)	mg/dL	1415	10%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	1240	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	1380	10%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	1520	10%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1411	6.3%	<5	<100
(IgM)免疫球蛋白 IgM					
Siemens BN II Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/dL	210	2.9%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	182	4.2%	<5	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	mg/dL	172	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Roche \ Chemiluminescence enzyme immunoassay (CLEIA)	mg/dL	161	1.6%	<5	<100
Abbott ARCHITECT c series DiaSys \ Immunoturbidimetry (ITA)	mg/dL	163	10%	<5	<100
Olympus AU400 Olympus \ Immunonephelometry (INA)	mg/dL	165	10%	<5	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	179	10%	<5	<100
KHB Polaris i2400 KHB \ Chemiluminescence Immunoassay (CLIA)	mg/dL	262	10%	<5	<100
(Kappa Light Chain)κ轻链抗体					
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	10.6	10%	<5	<100
(Lambda Light Chain)λ 轻链抗体					
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	5.68	10%	<5	<100
(PAB)前白蛋白					
Mindray BS series OTHER \ OTHER	g/L	0.325	1.5%	<5	<100
Beckman AU series BSBE \ OTHER	g/L	3.46	18.4%	<5	<100
Abbott ARCHITECT c series BSBE \ OTHER	g/L	3.09	4.9%	10	218

质控品名： 特定蛋白液体非定值质控品

REF	521-2	6×3mL	EXP	2026-05-22	LOT	31B52C1	更新日期：	2024-04
3 / 3								

批号：31B52C1					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
TOSHIBA TBA	g/L	0.357	1.1%	<5	<100
TOSHIBA \ Immunoturbidimetry (ITA)					
KHB Polaris i2400	g/L	0.321	10%	<5	<100
KHB \ Chemiluminescence Immunoassay (CLIA)					
Beckman AU series	g/L	0.32	10%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
Olympus AU400	g/L	0.315	10%	<5	<100
Olympus \ Immunonephelometry (INA)					
Siemens ADVIA Centaur XP/XPT	g/L	287	0.4%	<5	<100
Roche \ Chemiluminescence enzyme immunoassay (CLEIA)					
(RBP)视黄醇结合蛋白					
Siemens BN II	mg/dL	4.73	9.9%	<5	<100
Siemens \ Immunonephelometry (INA)					
(RF)类风湿因子					
Abbott ARCHITECT c series	IU/mL	77.2	2.7%	5	105
BSBE \ OTHER					
Beckman AU series	IU/mL	71.6	10.6%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
Mindray BS series	IU/mL	80.1	1.6%	<5	<100
MedicalSystem \ OTHER					
Siemens BN II	IU/mL	54.1	6.5%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman Image 800	IU/mL	67.1	10%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
Siemens BN ProSpec	IU/mL	55.2	3.3%	<5	<100
Siemens \ Immunonephelometry (INA)					
Mindray BS series	IU/mL	60.1	10%	<5	<100
Mindray \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	IU/mL	74.2	0.7%	<5	<100
TOSHIBA \ Immunoturbidimetry (ITA)					
KHB Polaris i2400	IU/mL	76.9	10%	<5	<100
KHB \ Chemiluminescence Immunoassay (CLIA)					
(TRF)转铁蛋白					
Siemens BN II	mg/dL	237	2.8%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Olympus AU400	mg/dL	235	10%	<5	<100
Olympus \ Immunoturbidimetry (ITA)					
(α-1-AG)α-1-酸性糖蛋白					
Siemens BN II	g/L	1.1	7.9%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman AU series	g/L	0.698	1.7%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
Olympus AU400	g/L	0.772	10%	<5	<100
Olympus \ Immunoturbidimetry (ITA)					
(α-1-AT)α-1-抗胰蛋白酶					
Beckman AU series	mg/dL	152	5.4%	<5	<100
Beckman Coulter \ Immunoturbidimetry (ITA)					
Olympus AU400	mg/dL	149	10%	<5	<100
Olympus \ Immunoturbidimetry (ITA)					
(α2M)α-2-巨球蛋白					
Beckman Image 800	mg/L	2180	10%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					
(β2M)β-2-微球蛋白					
Mindray BS series	ng/mL	4181	1.1%	<5	<100
MedicalSystem \ OTHER					
Beckman AU series	ng/mL	3596	14.2%	<5	<100
BSBE \ OTHER					
Abbott ARCHITECT c series	ng/mL	3621	5.6%	10	217
BSBE \ OTHER					
Olympus AU400	ng/mL	4280	9.3%	<5	<100
Olympus \ Immunonephelometry (INA)					
KHB Polaris i2400	ng/mL	3340	10%	<5	<100
KHB \ Chemiluminescence Immunoassay (CLIA)					
Mindray BS series	ng/mL	3970	10%	<5	<100
Mindray \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	ng/mL	2828	1%	<5	<100
TOSHIBA \ Immunoturbidimetry (ITA)					
Beckman AU series	ng/mL	3980	2.8%	<5	<100
Beckman Coulter \ Immunonephelometry (INA)					

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。