

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
1 / 10								

		批号：33876EF			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(25-OH-Vit D)25羟基维生素D					
Roche cobas E601/602	ng/mL	70.2	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	nmol/L	175	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(AFP)甲胎蛋白					
Beckman UniCel Dxl 600/800	ng/mL	196	8.6%	6	<100
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
Abbott ARCHITECT i series	ng/mL	226	6.4%	46	905
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	253	4.6%	29	799
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	ng/mL	235	6.4%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	254	0.7%	17	457
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens Atellica	ng/mL	243	7%	15	449
Siemens \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	ng/mL	252	3.4%	14	309
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	ng/mL	319	8.1%	<5	<100
Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens Immulite1000/2000	ng/mL	230	12%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	ng/mL	233	12%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Caris200	ng/mL	191	12%	<5	<100
WanTai \ Chemiluminescence Immunoassay (CLIA)					
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS	pg/mL	5360	12%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(Anti-Tg)抗甲状腺球蛋白抗体					
Beckman UniCel Dxl 600/800	IU/mL	2.33	7.5%	6	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
YHLO iFlash Analyzer	IU/mL	73.6	8.1%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	IU/mL	19.7	17.5%	25	541
Abbott \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	IU/mL	52.7	4.6%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	IU/mL	32.5	0.7%	17	459
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	IU/mL	35.2	8.6%	14	501
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	IU/mL	1.83	14.8%	10	224
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Simens IM1600	IU/mL	27.2	8.1%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	28.7	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Beckman UniCel Dxl 600/800	IU/mL	14	3.9%	6	<100
Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)					
YHLO iFlash Analyzer	IU/mL	22.2	3.3%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	IU/mL	25.5	4.8%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	IU/mL	24.5	7.4%	25	526
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	IU/mL	14.5	1%	17	458
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas 4000/6000/8000	IU/mL	21.2	12.5%	14	610
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman UniCel Dxl 600/800	IU/mL	14.9	8.6%	10	253
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	13.8	12%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(CA 125)糖类抗原 125					
Abbott ARCHITECT i series	IU/mL	179	37.6%	33	644
Abbott \ Chemiluminescence Immunoassay (CLIA)					
Beckman UniCel Dxl 600/800	IU/mL	113	3.4%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	94.5	0.3%	17	458
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Siemens Atellica	IU/mL	107	4%	14	425
Siemens \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	IU/mL	88.7	3.1%	14	318
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
2 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	130	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	96.3	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	166	12%	<5	<100
(CA 15-3)糖类抗原 15-3					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	96.7	21.5%	34	709
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	38.9	3.1%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	74.8	1.2%	17	459
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.8	15.1%	15	429
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	81.1	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	76.8	12%	<5	<100
(CA 19-9)糖类抗原 19-9					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	666	20.7%	34	763
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	207	3.6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	114	1.6%	17	459
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	265	9.5%	15	448
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	98	4.8%	14	304
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	78.7	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	194	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	106	12%	<5	<100
(CEA)癌胚抗原					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	30.8	8.1%	47	1006
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.4	3.8%	29	834
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	25.6	3.7%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	28.7	1.2%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23.3	11.1%	15	446
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.9	2.3%	14	308
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.4	3.2%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.3	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	43	12%	<5	<100
(CK)肌酸激酶					
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	117	12%	<5	<100
(CK-MB)肌酸激酶-MB (活性)					
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	30.5	13.1%	8	168
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	32.8	6.4%	7	<100
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	24	4.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	28.2	8.1%	24	613
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	27.5	6.5%	20	532
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	29.5	12%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	27.6	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	46.8	12%	<5	<100
(Cortisol)皮质醇					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	283	6.4%	8	216

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
3 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	299	10.2%	7	207
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	299	8.3%	6	<100
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	370	7.7%	<5	<100
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	147	8.4%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	275	3.3%	24	644
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	321	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	364	12%	<5	<100
(C-P)C-肽					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.76	3.2%	7	205
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.74	9.2%	7	128
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.25	9.1%	21	411
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.5	1.4%	17	457
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.88	6.1%	15	452
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	4.57	5.2%	14	270
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.29	3.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.98	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.68	12%	<5	<100
(Cyfra 21-1)细胞角蛋白19片段					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.3	5.2%	10	191
(DGXN)地高辛					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.98	12%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8023	7.9%	5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8821	4.2%	14	473
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9875	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9285	12%	<5	<100
(E2)雌二醇					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1793	2.9%	7	232
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1630	6.3%	7	202
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	498	5.7%	7	165
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1637	7.8%	6	<100
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1855	6%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1794	2.9%	28	1316
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1715	6.1%	<5	<100
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1762	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1454	10.4%	19	366
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1798	1%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1836	7%	15	446
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	1394	7.6%	15	300
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1461	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1711	1.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1510	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1850	12%	<5	<100

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
4 / 10								

		批号：33876EF			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(EPO)促红细胞生成素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	128	5%	9	209
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/mL	130	5.6%	6	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	149	12%	<5	<100
(FA)叶酸					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.1	5%	9	209
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.4	6.6%	8	215
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	14.9	5.7%	6	<100
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	16.2	6.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.2	4.5%	18	505
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.6	8.3%	13	179
(Ferritin)铁蛋白					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	388	6.3%	9	209
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	386	7.2%	6	<100
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	619	4.4%	<5	105
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	507	24.4%	33	623
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	545	0.5%	17	459
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	574	3%	14	<100
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	448	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	420	12%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	10.5	6.5%	6	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.14	34.2%	33	661
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.6	4.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.36	3.1%	18	487
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.83	12%	<5	<100
(FSH)促卵泡刺激素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	50.7	7.9%	8	217
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	48.3	7.2%	7	206
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	54.1	4%	7	204
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	73.6	5%	7	169
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	47.6	4.9%	6	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	56.2	10.5%	33	664
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	85.8	5.8%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	55.8	3%	28	807
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	52.3	1.1%	17	456
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	69.1	12.5%	15	444
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	63.8	6%	15	307
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	55.4	4.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	45.5	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	39.2	12%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	65.6	12%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
5 / 10								

批号：33876EF					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	56.7	12%	<5	<100
(FT3)游离T3					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17.4	5.7%	7	205
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.9	3.7%	7	185
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	16.1	9%	6	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	28	9.1%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	23.2	7.7%	35	723
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	16	5%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17	6%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	22	6.9%	17	878
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	28.4	0.6%	17	458
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	26.5	2.7%	14	416
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	24.9	6.1%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.9	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	19.9	12%	<5	<100
(FT4)游离T4					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	46.1	5.2%	7	<100
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	47.6	4.8%	7	207
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	44.7	4%	7	185
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	55.7	9.8%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.8	7.7%	35	746
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	31.6	3.8%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	48.9	2.7%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.5	4.8%	17	864
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	68.8	1.7%	17	458
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	53.6	6.2%	14	419
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	48.8	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	47.3	3%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	pmol/L	40.8	12%	<5	<100
(hCG)人绒毛膜促性腺激素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	376	4.9%	6	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	406	2.8%	6	145
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	422	2.8%	<5	<100
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	411	4.3%	<5	<100
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	310	16.8%	15	435
Abbott ARCHITECT i series Abbott \ OTHER	IU/L	402	8.6%	15	369
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	648	6.9%	14	297
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	615	12%	<5	<100
(HCY)同型半胱氨酸					
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	39.6	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	40.3	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	32.5	2.5%	<5	<100

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
6 / 10								

		批号：33876EF			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(HE4)人附睾蛋白4					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	199	23.2%	21	413
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	185	2.6%	14	309
(HGH)人生长激素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	7.5	5.6%	6	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	9.66	12%	<5	<100
(IgE)免疫球蛋白 IgE					
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	321	3%	7	208
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	251	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	273	12%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	202	12%	<5	<100
(INS)胰岛素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	129	3.9%	7	217
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	112	4.2%	7	172
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	132	3.1%	21	614
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	98.4	6.7%	20	407
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	136	1.7%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	114	4.1%	15	444
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	148	5.9%	14	284
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	118	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	126	12%	<5	<100
(IntactPTH)全段甲状旁腺激素					
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33	15.2%	9	<100
Abbott ARCHITECT i series Abbott \ OTHER	pg/mL	31.9	13.1%	11	218
Abbott ARCHITECT i series Abbott \ OTHER	pmol/L	33.3	4.1%	<5	<100
(LH)促黄体生成素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	42.2	5.2%	8	217
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	43.5	5.4%	7	205
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	55.8	2.6%	7	205
Siemens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	53.3	3.1%	7	173
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	39.4	5.4%	5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	43.4	10.1%	35	758
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	53.2	8.6%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	55.9	3.9%	28	796
Tosoh AIA 900/1800 OTHER \ Chemiluminescence Immunoassay (CLIA)	IU/L	53.8	3.6%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	55.9	1.1%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	52.9	5%	15	442
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/L	75.6	6.4%	15	312
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	56.6	3.5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	56.4	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/L	78.1	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	47.5	12%	<5	<100

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
7 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	72.5	12%	<5	<100
(Myo / Mb)肌红蛋白					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	238	5.3%	9	209
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	255	3.5%	8	170
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	230	3.1%	7	<100
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	μg/L	317	3.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	280	4.2%	24	607
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	247	4.4%	20	538
Ortho Vitros 3600/5600 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	338	7.6%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μg/L	239	5.6%	<5	<100
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	μg/L	345	6.5%	15	412
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	289	12%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2351	4.9%	24	613
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	277	4.9%	24	613
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2444	6.4%	20	530
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	287	6.4%	20	530
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2593	4.2%	11	208
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	305	4.2%	11	208
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	2210	12%	<5	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	260	12%	<5	<100
(PRL)催/ 泌乳素					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	31	4.8%	7	209
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	37.2	2.9%	7	195
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.5	5%	7	171
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	27.4	4.6%	6	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	35	2.9%	28	797
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	38	13.3%	18	389
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	34.8	1.8%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23.5	9.4%	15	428
Siemens Immulite1000/2000 Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	31.3	5%	15	310
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	32.7	1.8%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	30.1	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	38.2	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28.8	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26.8	3.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	36.9	12%	<5	<100
(PROG)孕酮					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.8	7.2%	7	203
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.78	3.5%	7	134
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	20.8	10%	6	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20.3	9.6%	5	124

质控品名： 复合免疫分析非定值质控品

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
8 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.5	13%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	21.7	2.5%	27	1159
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	22	2.8%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.6	7.3%	19	380
Roche cobas E601/602 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	21.5	0.8%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	15.8	5.2%	15	445
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	17.1	6.8%	15	281
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19	3.5%	<5	<100
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.1	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	24.2	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.5	12%	<5	<100
(proGRP)胃泌素释放肽前体					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/L	8680	7.8%	16	334
(PSA)前列腺特异性抗原					
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	23.8	8%	6	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.1	19.7%	33	678
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	24.6	3.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	20	5.8%	18	504
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.7	6.3%	15	450
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28.5	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.9	12%	<5	<100
(PTH)甲状旁腺素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	23.3	4.9%	6	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	19.4	0.9%	17	458
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	21.3	11.6%	11	257
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	25.8	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	20.2	12%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.2	26.8%	<5	<100
(SCC)鳞状细胞癌抗原					
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.7	8.6%	16	332
(TES)睾酮					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.65	4.6%	7	210
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.41	4.8%	7	205
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.479	4.4%	7	171
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.53	4.7%	6	<100
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.88	11.4%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.52	5.2%	28	807
Tosoh AIA 900/1800 OTHER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.49	3.3%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.42	11.5%	19	411
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.6	3.1%	17	458
Siemens Atellica Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.71	8.2%	15	436
Siemens Immulite1000/2000 Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.95	8.4%	15	307
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.49	12%	<5	<100

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
9 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.29	3.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.54	12%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.99	12%	<5	<100
(Tg)甲状腺球蛋白					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	6.74	4.8%	6	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.89	5.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.3	12%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.71	4.1%	24	604
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.81	6.5%	20	529
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.63	3.4%	11	208
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.98	12%	<5	<100
(TSH)促甲状腺素					
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	27.4	5.2%	7	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.5	5%	7	206
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31.8	4.6%	7	181
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	36.4	9.3%	<5	<100
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	37.3	4.1%	<5	<100
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	34.2	3.5%	28	778
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	25.5	4.7%	28	623
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29.8	3.9%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	24.2	4.8%	17	866
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.6	1.3%	17	458
Abbott ARCHITECT i series Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	26.5	6.3%	14	280
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.9	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	29.4	12%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.6	2.7%	<5	<100
(TT3)总三碘甲状腺原氨酸					
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.01	5.1%	8	164
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.31	7.8%	7	207
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	3.1	7%	6	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.46	8.1%	37	817
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.6	2.7%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.61	4%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.83	7.7%	17	445
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.67	1.4%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.67	12%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.54	12%	<5	<100
(TT4)总T4					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	204	10.5%	7	211
Simens IM1600 Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	185	5.6%	7	185
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	164	9.5%	7	<100

REF	540-3	6×3mL	EXP	2025-11-16	LOT	33876EF	更新日期：	2024-05
10 / 10								

		批号：33876EF			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	166	8.9%	37	752
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	ng/mL	162	3.2%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	176	6.4%	<5	<100
Abbott Alinity ci-series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	162	5.6%	17	436
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	162	3.8%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	ng/mL	133	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	143	12%	<5	<100
(VB12)维生素B12					
maglumi system Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	386	6.2%	9	272
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	712	6.3%	9	208
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	670	5.1%	6	<100
Abbott ARCHITECT i series Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	797	4.2%	<5	104
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1067	5.3%	18	497
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pg/mL	852	10.7%	15	222
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	840	7.9%	10	234
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	976	12%	<5	<100
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	2770	12%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2810	12%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Beckman AU series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.41	4.1%	7	205
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.419	2%	6	139
Siemens ADVIA Chemistry XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.342	7.5%	<5	<100
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.395	2.3%	<5	<100
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.395	4.8%	<5	<100
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.413	1%	17	458
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.396	2.4%	14	514
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.471	1.9%	<5	<100
Caris200 WanTai \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.413	12%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.522	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。