

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 1 / 12 |       |       |     |            |     |         |       |         |

| 批号：33828CF                                                                              |        |         |       |          |      |
|-----------------------------------------------------------------------------------------|--------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位     | 均值 \ CV |       | Lab \ 点数 |      |
| (17-OHP)17-羟孕酮                                                                          |        |         |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 0.798   | 5.8%  | <5       | <100 |
| Agilent<br>agilent \ MS                                                                 | ng/mL  | 0.759   | 12%   | <5       | <100 |
| (25-OH-VD3)25羟基维生素D3                                                                    |        |         |       |          |      |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 22      | 8.3%  | 12       | 178  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 23.3    | 2.7%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 20      | 12%   | <5       | <100 |
| (25-OH-Vit D)25羟基维生素D                                                                   |        |         |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 21      | 6.8%  | 8        | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | nmol/L | 52.5    | 6.8%  | 8        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 20.7    | 8.4%  | <5       | 137  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | nmol/L | 51.6    | 8.4%  | <5       | 137  |
| Beckman Access<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )               | ng/mL  | 21.8    | 12%   | <5       | <100 |
| Beckman Access<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )               | nmol/L | 54.4    | 12%   | <5       | <100 |
| (AFP)甲胎蛋白                                                                               |        |         |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 5.79    | 5.5%  | 7        | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 6.96    | 6%    | 7        | 130  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 7.42    | 4.8%  | 30       | 884  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 6.26    | 5.3%  | 29       | 516  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 7.29    | 3.9%  | 14       | 309  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 8.03    | 9.8%  | 11       | <100 |
| maglumi system<br>Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)                 | ng/mL  | 10      | 3.3%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 7.34    | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL  | 6.51    | 12%   | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 5.91    | 11.8% | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 6.65    | 12%   | <5       | <100 |
| (ALD)醛固酮                                                                                |        |         |       |          |      |
| DiaSorin Liaison /XL<br>Diasorin \ Chemiluminescence Immunoassay (CLIA )                | pg/mL  | 70.3    | 12%   | <5       | <100 |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | pg/mL  | 134     | 10%   | <5       | <100 |
| (AMH)抗缪勒氏管激素                                                                            |        |         |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 4.73    | 4.5%  | 7        | <100 |
| (Anti-Tg)抗甲状腺球蛋白抗体                                                                      |        |         |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL  | 1.71    | 9%    | 6        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL  | 16.9    | 27.6% | <5       | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL  | 21.9    | 4.6%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL  | 11.8    | 21.4% | 24       | 538  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL  | 18.6    | 14.1% | 22       | 532  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL  | 1.48    | 20.7% | 11       | 260  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL  | 9.97    | 13.5% | 10       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | IU/mL  | 20.9    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)      | IU/mL  | 9.91    | 9.3%  | <5       | <100 |

质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 2 / 12 |       |       |     |            |     |         |       |         |

| 批号：33828CF                                                                              |       |         |       |          |      |
|-----------------------------------------------------------------------------------------|-------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位    | 均值 \ CV |       | Lab \ 点数 |      |
| (Anti-TPO)抗甲状腺过氧化物酶抗体                                                                   |       |         |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | IU/mL | 6.89    | 5.6%  | 6        | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL | 14.4    | 9.9%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 15      | 15.6% | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 12.4    | 6.3%  | 25       | 517  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL | 12.9    | 19.8% | 21       | 641  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL | 24.4    | 8.2%  | 10       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 7.39    | 10.2% | 10       | 263  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | IU/mL | 15      | 8%    | <5       | <100 |
| (ASD)雄烯二酮                                                                               |       |         |       |          |      |
| Simens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                        | ng/mL | 2.55    | 10.1% | 10       | <100 |
| (CA 125)糖类抗原 125                                                                        |       |         |       |          |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL | 29.2    | 3.4%  | 7        | 129  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 33.3    | 5.1%  | <5       | 102  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 39.7    | 3.5%  | 19       | 396  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL | 28.3    | 3.2%  | 14       | 318  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL | 34.9    | 2.6%  | 11       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | IU/mL | 43.1    | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 33.9    | 10.6% | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/mL | 48.5    | 12%   | <5       | <100 |
| (CA 15-3)糖类抗原 15-3                                                                      |       |         |       |          |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL | 29.6    | 4.4%  | 7        | 125  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 28.5    | 7.1%  | <5       | 111  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 31.6    | 6.7%  | 19       | 372  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL | 29.8    | 10.4% | 11       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 16.2    | 11.7% | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/mL | 29.4    | 12%   | <5       | <100 |
| (CA 19-9)糖类抗原 19-9                                                                      |       |         |       |          |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL | 27.8    | 10.3% | 8        | 137  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 29.4    | 6.5%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 170     | 8.4%  | 19       | 396  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL | 24.7    | 5.1%  | 14       | 304  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL | 23.3    | 5.9%  | 10       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | IU/mL | 45.3    | 12%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/mL | 20.6    | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 26.2    | 10.7% | <5       | <100 |
| (CEA)癌胚抗原                                                                               |       |         |       |          |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL | 3.53    | 9.4%  | 7        | 129  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 3.43    | 5.6%  | 30       | 866  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 3.39    | 10.2% | 30       | 581  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL | 3.43    | 4.8%  | 14       | 308  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | ng/mL | 3.67    | 4.9%  | 11       | <100 |
| maglumi system<br>Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)                 | ng/mL | 3.43    | 18.4% | <5       | <100 |

质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 3 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |       |            |       |          |      |
|-----------------------------------------------------------------------------------------|-------|------------|-------|----------|------|
|                                                                                         |       | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位    | 均值 \ CV    |       | Lab \ 点数 |      |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL | 2.99       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 2.55       | 12%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 4.9        | 12%   | <5       | <100 |
| (CK)肌酸激酶                                                                                |       |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Roche \ Urease-UV-rate                                    | U/L   | 130        | 10.1% | <5       | <100 |
| (CK-MB)肌酸激酶-MB (活性)                                                                     |       |            |       |          |      |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 4.13       | 13.3% | 8        | 317  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL | 3.83       | 8%    | 7        | <100 |
| Roche Hitachi/Hitachi Systems<br>MAKER \ Immune suppression                             | U/L   | 7.56       | 14.9% | 6        | <100 |
| TOSHIBA TBA<br>MAKER \ Immune suppression                                               | U/L   | 7.52       | 10.4% | <5       | <100 |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)      | ng/mL | 2.96       | 3.7%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 3.78       | 8.6%  | 25       | 653  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL | 3.79       | 5.7%  | 19       | 507  |
| Roche cobas 4000/6000/8000<br>Roche \ Immunoturbidimetry (ITA)                          | U/L   | 9.25       | 10%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 5.08       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 4.5        | 12%   | <5       | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL | 4.49       | 10%   | <5       | <100 |
| (Cortisol)皮质醇                                                                           |       |            |       |          |      |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | ng/mL | 97.3       | 7.3%  | 9        | 115  |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL | 98.3       | 10%   | 7        | 203  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)     | ng/mL | 117        | 4.9%  | 5        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 137        | 6.7%  | <5       | 109  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL | 109        | 5.4%  | 14       | 194  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL | 108        | 5.8%  | 13       | 114  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 93.9       | 7.9%  | 11       | 286  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 115        | 12%   | <5       | <100 |
| (C-P)C-肽                                                                                |       |            |       |          |      |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 0.879      | 3%    | 7        | 113  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 0.827      | 10.4% | 20       | 382  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 1.12       | 22.9% | <5       | 104  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL | 1.03       | 3.9%  | 15       | 297  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)     | ng/mL | 0.955      | 3.5%  | 14       | 269  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL | 1.09       | 12.9% | 13       | 176  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | ng/mL | 0.696      | 15.1% | 11       | 207  |
| maglumi system<br>Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)                 | ng/mL | 0.401      | 2.4%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 0.648      | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL | 0.858      | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 0.93       | 9.7%  | <5       | <100 |
| (Cyfra 21-1)细胞角蛋白19片段                                                                   |       |            |       |          |      |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 0.93       | 8.4%  | 17       | 342  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | ng/mL | 0.973      | 10.5% | 11       | <100 |

质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 4 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        | 批号：33828CF |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| (DGXN)地高辛                                                                               |        |            |       |          |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 0.88       | 10%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 0.87       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL  | 0.859      | 12%   | <5       | <100 |
| (DHEA-S)硫酸脱氢表雄酮                                                                         |        |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 3192       | 4.6%  | 14       | 473  |
| Simens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                        | ng/mL  | 2640       | 10.3% | 10       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 3390       | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 3416       | 10%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 3906       | 10.4% | <5       | <100 |
| (E2)雌二醇                                                                                 |        |            |       |          |      |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | pmol/L | 309        | 6.1%  | 7        | 205  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | pmol/L | 314        | 12.4% | <5       | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | pmol/L | 332        | 2.4%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | pmol/L | 330        | 4.5%  | 28       | 1316 |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pmol/L | 318        | 7.5%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | pmol/L | 319        | 9.1%  | 17       | 332  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | pmol/L | 245        | 11.6% | 15       | 306  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | pmol/L | 322        | 3.5%  | 13       | 111  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | pmol/L | 342        | 4.6%  | 12       | 177  |
| Simens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                        | pmol/L | 87.6       | 8.7%  | 11       | 284  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | pmol/L | 246        | 12%   | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | pmol/L | 331        | 1.7%  | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | pmol/L | 275        | 12%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pmol/L | 351        | 10%   | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | pmol/L | 327        | 1.8%  | <5       | <100 |
| (EPO)促红细胞生成素                                                                            |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | mIU/mL | 5.95       | 4.8%  | 13       | 113  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | mIU/mL | 7.28       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | mIU/mL | 6.76       | 9.9%  | <5       | <100 |
| (FA)叶酸                                                                                  |        |            |       |          |      |
| maglumi system<br>Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)                 | ng/mL  | 12.1       | 4%    | 9        | 270  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 3.65       | 13.4% | 9        | 231  |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 8.45       | 10%   | 7        | 104  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 3.51       | 20.1% | 5        | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL  | 3.32       | 13%   | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 4.45       | 22%   | 13       | 297  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 3.77       | 9.6%  | 13       | 113  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 3.76       | 12.8% | 10       | 140  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 3.66       | 12%   | <5       | <100 |

质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 5 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        | 批号：33828CF |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| (Ferritin)铁蛋白                                                                           |        |            |       |          |      |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 78.3       | 7.2%  | 7        | 102  |
| Abbott ARCHITECT i series<br>Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL  | 80         | 6.5%  | <5       | 105  |
| Roche cobas 4000/6000/8000<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL  | 113        | 4.1%  | 14       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 82.6       | 8.4%  | 14       | 256  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 67.4       | 9%    | 13       | 113  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 105        | 7.2%  | 12       | 178  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 112        | 5%    | 11       | <100 |
| DiaSorin Liaison /XL<br>Diasorin \ Chemiluminescence Immunoassay (CLIA )                | ng/mL  | 95.9       | 12%   | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 69.2       | 9.8%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Roche \ Immunoturbidimetry (ITA)                          | ng/mL  | 119        | 10%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 68.1       | 9.4%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 80.5       | 12%   | <5       | <100 |
| (F-PSA)游离前列腺特异性抗原                                                                       |        |            |       |          |      |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 2.02       | 2.9%  | 8        | 111  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 2.07       | 6.4%  | 6        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 1.87       | 7.7%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 2.09       | 5%    | 14       | 254  |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL  | 1.81       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 2.06       | 12%   | <5       | <100 |
| (FSH)促卵泡刺激素                                                                             |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/L   | 5.38       | 9%    | 9        | 256  |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | IU/L   | 5.4        | 5.6%  | 7        | 202  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/L   | 5.93       | 2.8%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/L   | 5.55       | 3.4%  | 28       | 806  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/L   | 5.65       | 7.2%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/L   | 5.35       | 8.9%  | 18       | 384  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | IU/L   | 5.56       | 5.3%  | 15       | 305  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | IU/L   | 5.3        | 4.9%  | 13       | 111  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | IU/L   | 4.87       | 7.4%  | 12       | 175  |
| Simens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                        | IU/L   | 5.75       | 6.9%  | 11       | 269  |
| Roche cobas 4000/6000/8000<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | IU/L   | 5.72       | 2.7%  | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | IU/L   | 4.51       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPt<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | IU/L   | 6.47       | 12.1% | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/L   | 5.92       | 12%   | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | IU/L   | 5.47       | 1.9%  | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | IU/L   | 5.81       | 12%   | <5       | <100 |
| (FT3)游离T3                                                                               |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | pmol/L | 4.87       | 7.2%  | 7        | <100 |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | pmol/L | 5.12       | 6.2%  | 7        | 207  |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | pmol/L | 5.5        | 10.1% | 6        | 109  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pmol/L | 5.99       | 8.5%  | 5        | 151  |

质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 6 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        |            |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
|                                                                                         |        | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | pmol/L | 6.67       | 4%    | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | pmol/L | 5.45       | 6.9%  | 35       | 741  |
| Abbott Alinity ci-series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | pmol/L | 5.54       | 4.6%  | <5       | 158  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | pmol/L | 5.94       | 6.4%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | pmol/L | 6.28       | 4.2%  | 22       | 449  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | pmol/L | 6.42       | 6.9%  | 14       | 195  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | pmol/L | 6.69       | 2.9%  | 11       | 274  |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | pmol/L | 4.99       | 12%   | <5       | <100 |
| Beckman UniCel DxI 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | pmol/L | 4.5        | 12%   | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | pmol/L | 6.39       | 2.2%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | pmol/L | 5.18       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | pmol/L | 5.82       | 12%   | <5       | <100 |
| OTHER<br>TELLGEN \ Chemiluminescence Immunoassay (CLIA )                                | pmol/L | 8.36       | 4.8%  | <5       | <100 |
| (FT4 )游离T4                                                                              |        |            |       |          |      |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | pmol/L | 10.9       | 6.1%  | 7        | <100 |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | pmol/L | 12.3       | 7.2%  | 7        | 204  |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | pmol/L | 16.9       | 11.2% | 6        | 108  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pmol/L | 15.9       | 7.2%  | 5        | 143  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | pmol/L | 12.7       | 5.6%  | 35       | 768  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | pmol/L | 15.7       | 4.2%  | <5       | <100 |
| Abbott Alinity ci-series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | pmol/L | 12.8       | 4.8%  | <5       | 156  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | pmol/L | 17.3       | 7.4%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | pmol/L | 16         | 4.4%  | 22       | 452  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | pmol/L | 16.7       | 6.1%  | 14       | 196  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | pmol/L | 14.4       | 4.9%  | 11       | 270  |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | pmol/L | 15.1       | 11.9% | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | pmol/L | 15.1       | 2.8%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | pmol/L | 14.1       | 12%   | <5       | <100 |
| Beckman UniCel DxI 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | pmol/L | 12.3       | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | pmol/L | 12.3       | 12%   | <5       | <100 |
| (hCG)人绒毛膜促性腺激素                                                                          |        |            |       |          |      |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | IU/L   | 3.73       | 5.8%  | 7        | <100 |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | IU/L   | 3.6        | 4.1%  | 6        | 146  |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/L   | 2.64       | 7.3%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ OTHER                                             | IU/L   | 3.93       | 16.6% | 14       | 352  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | IU/L   | 5.94       | 8.8%  | 14       | 301  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/L   | 4.22       | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | IU/L   | 4.32       | 2.7%  | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | IU/L   | 4.41       | 12%   | <5       | <100 |
| (HCY)同型半胱氨酸                                                                             |        |            |       |          |      |
| TOSHIBA TBA<br>Bio-assay \ Enzymatic Cycling Assays Method                              | μmol/L | 11.5       | 7.7%  | 5        | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | μmol/L | 10.8       | 12.5% | <5       | <100 |



|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 7 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        |            |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
|                                                                                         |        | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| Roche cobas 4000/6000/8000<br>Roche \ Enzymatic Cycling Assays Method                   | μmol/L | 13.9       | 10%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | μmol/L | 11.5       | 12%   | <5       | <100 |
| (HGH)人生长激素                                                                              |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | μg/L   | 1.15       | 5.6%  | 6        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | μg/L   | 1.24       | 40.8% | 15       | 291  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | μg/L   | 1.51       | 12%   | <5       | <100 |
| (HS-CTnI) (高) 超敏肌钙蛋白I                                                                   |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 0.066      | 5.6%  | 7        | <100 |
| Beckman Access<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )               | ng/mL  | 0.061      | 12%   | <5       | <100 |
| (HS-CTnT)超敏肌钙蛋白T                                                                        |        |            |       |          |      |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL  | 0.019      | 5.7%  | <5       | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 0.017      | 7.4%  | <5       | <100 |
| (IgE)免疫球蛋白 IgE                                                                          |        |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL  | 182        | 3.1%  | 7        | 206  |
| Beckman Image 800<br>Beckman Coulter \ Immunoturbidimetry (ITA)                         | IU/mL  | 131        | 10%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | IU/mL  | 160        | 12%   | <5       | <100 |
| (IGF1)胰岛素生长因子1                                                                          |        |            |       |          |      |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 217        | 9.5%  | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | nmol/L | 17.7       | 12%   | <5       | <100 |
| (INS)胰岛素                                                                                |        |            |       |          |      |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | mIU/L  | 8.6        | 2.9%  | 7        | 101  |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | mIU/L  | 8.84       | 4%    | 29       | 707  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | mIU/L  | 7.24       | 8.3%  | 20       | 381  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)     | mIU/L  | 12.5       | 5.7%  | 14       | 283  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | mIU/L  | 9.26       | 5%    | 12       | 177  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | mIU/L  | 7.26       | 4.5%  | 11       | 264  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | mIU/L  | 9.01       | 10%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | mIU/L  | 8.83       | 10%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | mIU/L  | 6.95       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | mIU/L  | 7.1        | 12%   | <5       | <100 |
| (IntactPTH)全段甲状旁腺激素                                                                     |        |            |       |          |      |
| Abbott ARCHITECT i series<br>Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pmol/L | 3.17       | 10.3% | 9        | <100 |
| Abbott ARCHITECT i series<br>Abbott \ OTHER                                             | pmol/L | 3.17       | 9.8%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ OTHER                                             | pg/mL  | 3.09       | 13.7% | 11       | 222  |
| (LH)促黄体生成素                                                                              |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/L   | 3.69       | 7.2%  | 9        | 260  |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | IU/L   | 3.83       | 6.4%  | 7        | 203  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/L   | 5.96       | 3.3%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/L   | 5.55       | 3.9%  | 28       | 795  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/L   | 5.71       | 7.5%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/L   | 3.37       | 4.7%  | 17       | 348  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | IU/L   | 5.38       | 6.1%  | 15       | 312  |

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 8 / 12 |       |       |     |            |     |         |       |         |

| 批号：33828CF                                                                             |      |         |       |          |      |
|----------------------------------------------------------------------------------------|------|---------|-------|----------|------|
| 项目\仪器\试剂方法                                                                             | 单位   | 均值 \ CV |       | Lab \ 点数 |      |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | IU/L | 3.78    | 4.9%  | 13       | 111  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | IU/L | 5.26    | 5.6%  | 12       | 177  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                      | IU/L | 4       | 4.7%  | 11       | 262  |
| Roche cobas 4000/6000/8000<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)    | IU/L | 6.2     | 0.7%  | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                        | IU/L | 4.99    | 1.9%  | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)             | IU/L | 5.36    | 2.7%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                    | IU/L | 3.79    | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                             | IU/L | 5.93    | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )        | IU/L | 4.56    | 12.1% | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )           | IU/L | 5.18    | 12%   | <5       | <100 |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )         | IU/L | 4.66    | 11.6% | <5       | <100 |
| Tosoh AIA 900/1800<br>OTHER \ Chemiluminescence Immunoassay (CLIA )                    | IU/L | 3.98    | 4.7%  | <5       | <100 |

(Myo / Mb)肌红蛋白

|                                                                                        |      |      |       |    |      |
|----------------------------------------------------------------------------------------|------|------|-------|----|------|
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)          | μg/L | 33   | 4.2%  | 8  | 320  |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | μg/L | 29.5 | 7.7%  | 7  | <100 |
| TOSHIBA TBA<br>OTHER \ OTHER                                                           | μg/L | 31.5 | 13.6% | 6  | <100 |
| Roche Hitachi/Hitachi Systems<br>Bio-assay \ immunoturbidimetric                       | μg/L | 29.5 | 28.6% | <5 | <100 |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)     | μg/L | 37.3 | 2.7%  | <5 | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)          | μg/L | 30.3 | 6.5%  | 25 | 668  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)             | μg/L | 31   | 7.2%  | 23 | 602  |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )            | μg/L | 39.7 | 4.9%  | <5 | <100 |
| Beckman UniCel DxI 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )  | μg/L | 26.7 | 12%   | <5 | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | μg/L | 29.1 | 10%   | <5 | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                    | μg/L | 34.2 | 12%   | <5 | <100 |

(NT-proBNP)N-末端脑钠肽前体

|                                                                                   |        |      |      |    |      |
|-----------------------------------------------------------------------------------|--------|------|------|----|------|
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pg/mL  | 176  | 5.2% | 25 | 668  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pmol/L | 20.7 | 5.2% | 25 | 668  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | pg/mL  | 184  | 4.3% | 23 | 593  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)        | pmol/L | 21.7 | 4.3% | 23 | 593  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | pg/mL  | 182  | 5.4% | 13 | 194  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | pmol/L | 21.5 | 5.4% | 13 | 194  |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pg/mL  | 195  | 4.2% | 11 | 359  |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pmol/L | 22.9 | 4.2% | 11 | 359  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)         | pg/mL  | 202  | 10%  | <5 | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)         | pmol/L | 23.8 | 10%  | <5 | <100 |

(PRL)催/ 泌乳素

|                                                                                         |       |      |       |    |      |
|-----------------------------------------------------------------------------------------|-------|------|-------|----|------|
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL | 7.3  | 4%    | 7  | 206  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL | 8.66 | 2.3%  | <5 | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL | 8.03 | 3%    | 28 | 794  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 8.57 | 6.6%  | <5 | 144  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 8.93 | 12.3% | 17 | 349  |
| Siemens Immulite1000/2000<br>Siemens \ Electro-Chemiluminescence Immunoassay (ECLIA)    | ng/mL | 6.7  | 5%    | 15 | 312  |



质控品名： 复合免疫分析非定值质控品

|        |       |       |     |            |     |         |       |         |
|--------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF    | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 9 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        |            |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
|                                                                                         |        | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 6.72       | 4%    | 13       | 110  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 7.94       | 5.8%  | 12       | 178  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | ng/mL  | 5.89       | 4%    | 11       | 264  |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL  | 8.3        | 2%    | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL  | 7.89       | 2.9%  | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 7.17       | 12%   | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 9.01       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL  | 6.48       | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL  | 9.12       | 12%   | <5       | <100 |
| (PROG)孕酮                                                                                |        |            |       |          |      |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL  | 1.2        | 11.3% | 7        | 207  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 0.918      | 7.9%  | 6        | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 1.29       | 15.9% | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 0.851      | 7.3%  | 28       | 1218 |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 0.939      | 4.8%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 1.02       | 9.4%  | <5       | 105  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 0.864      | 21.2% | 18       | 361  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | ng/mL  | 1.11       | 10.1% | 15       | 281  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 1.18       | 12.3% | 13       | 111  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 0.911      | 7.4%  | 12       | 178  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | ng/mL  | 0.301      | 9%    | 11       | 258  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL  | 0.671      | 5.8%  | <5       | <100 |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence Immunoassay (CLIA )            | ng/mL  | 1.09       | 11.9% | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL  | 2.54       | 5%    | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL  | 5.05       | 5%    | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL  | 0.92       | 3.2%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 1.01       | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL  | 0.942      | 12%   | <5       | <100 |
| (PSA)前列腺特异性抗原                                                                           |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 4.79       | 5.5%  | 6        | <100 |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 6.51       | 5.6%  | 6        | 107  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 4.78       | 7.2%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 4          | 9.8%  | 16       | 273  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 4.08       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 4.82       | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL  | 4.33       | 12%   | <5       | <100 |
| (PTH)甲状旁腺素                                                                              |        |            |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pmol/L | 2.03       | 8.4%  | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | pmol/L | 2.07       | 6.7%  | 13       | 113  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence enzyme immunoassay (CLEIA)     | pmol/L | 2.41       | 14%   | 12       | 259  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | pmol/L | 1.81       | 6.8%  | 12       | 178  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | pmol/L | 1.94       | 12%   | <5       | <100 |

|         |       |       |     |            |     |         |       |         |
|---------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF     | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 10 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |        |            |       |          |      |
|-----------------------------------------------------------------------------------------|--------|------------|-------|----------|------|
|                                                                                         |        | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位     | 均值 \ CV    |       | Lab \ 点数 |      |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )          | pmol/L | 1.96       | 9.7%  | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | pmol/L | 2.26       | 9.8%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | pmol/L | 2.9        | 9.3%  | <5       | <100 |
| (SCC)鳞状细胞癌抗原                                                                            |        |            |       |          |      |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 0.998      | 15.9% | 14       | 227  |
| (SHBG)性激素结合球蛋白                                                                          |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | nmol/L | 33.2       | 12%   | <5       | <100 |
| (TES)睾酮                                                                                 |        |            |       |          |      |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 5.37       | 5.1%  | 7        | <100 |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL  | 5.25       | 4.6%  | 7        | 205  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )         | ng/mL  | 5.59       | 5.8%  | 7        | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 6.13       | 2.8%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 6.15       | 4.9%  | 28       | 806  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 6.1        | 7.8%  | <5       | <100 |
| Tosoh AIA 900/1800<br>OTHER \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 6.67       | 3.6%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 6.5        | 9.7%  | 18       | 373  |
| Siemens Immulite1000/2000<br>Siemens \ Chemiluminescence enzyme immunoassay (CLEIA)     | ng/mL  | 5.96       | 7.4%  | 15       | 310  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL  | 6.29       | 3.7%  | 12       | 177  |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL  | 5.66       | 12%   | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 5.28       | 12%   | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL  | 5.67       | 3.5%  | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL  | 6.61       | 12%   | <5       | <100 |
| (Tg)甲状腺球蛋白                                                                              |        |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL  | 10.1       | 4.9%  | 8        | <100 |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL  | 7.75       | 3.1%  | 6        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 10.8       | 6.8%  | <5       | 132  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 7.42       | 6.7%  | 11       | 248  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 9.88       | 9.9%  | <5       | <100 |
| (THEO)茶碱                                                                                |        |            |       |          |      |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | µg/mL  | 2.25       | 4%    | <5       | <100 |
| Agilent<br>agilent \ MS                                                                 | µg/mL  | 8.7        | 12%   | <5       | <100 |
| (Tnl)肌钙蛋白I                                                                              |        |            |       |          |      |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)      | ng/mL  | 0.061      | 5.6%  | <5       | <100 |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )             | ng/mL  | 0.075      | 11.9% | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL  | 0.05       | 12%   | <5       | <100 |
| (TnT)肌钙蛋白T                                                                              |        |            |       |          |      |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 0.02       | 4.1%  | 25       | 672  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL  | 0.021      | 5.3%  | 19       | 505  |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL  | 0.017      | 7.1%  | 11       | 359  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL  | 0.02       | 10.2% | <5       | <100 |

质控品名： 复合免疫分析非定值质控品

|         |       |       |     |            |     |         |       |         |
|---------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF     | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 11 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |       | 批号：33828CF |       |          |      |
|-----------------------------------------------------------------------------------------|-------|------------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                             | 单位    | 均值 \ CV    |       | Lab \ 点数 |      |
| (TSH)促甲状腺素                                                                              |       |            |       |          |      |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | mIU/L | 2.81       | 3.3%  | 7        | <100 |
| Beckman AU series<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )            | mIU/L | 2.89       | 4.2%  | 7        | 202  |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | mIU/L | 3.23       | 5.3%  | 6        | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | mIU/L | 3.62       | 11%   | 5        | 129  |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | mIU/L | 3.49       | 2%    | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | mIU/L | 3.41       | 3.4%  | 36       | 869  |
| Abbott Alinity ci-series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | mIU/L | 2.39       | 3.8%  | <5       | 158  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | mIU/L | 2.65       | 6%    | 21       | 458  |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence enzyme immunoassay (CLEIA)      | mIU/L | 2.64       | 6.1%  | 14       | 280  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | mIU/L | 3.41       | 3.5%  | 14       | 194  |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | mIU/L | 3.16       | 4.6%  | 11       | <100 |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | mIU/L | 3          | 4.1%  | 11       | 264  |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | mIU/L | 3.23       | 4%    | <5       | <100 |
| Beckman UniCel DxI 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | mIU/L | 2.9        | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | mIU/L | 3          | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | mIU/L | 2.96       | 12.2% | <5       | <100 |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | mIU/L | 2.87       | 12%   | <5       | <100 |
| OTHER<br>TELLGEN \ Chemiluminescence Immunoassay (CLIA )                                | mIU/L | 3.07       | 2.6%  | <5       | <100 |
| (TT3)总三碘甲状腺原氨酸                                                                          |       |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL | 1.42       | 3.2%  | 8        | <100 |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL | 1.15       | 8.3%  | 7        | 205  |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL | 1.03       | 10.5% | 6        | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL | 1.43       | 8.1%  | 5        | <100 |
| Abbott Alinity ci-series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | ng/mL | 1.21       | 4.7%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 1.16       | 11%   | 21       | 455  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL | 1.4        | 4.2%  | 13       | 194  |
| Siemens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                       | ng/mL | 1.37       | 6.8%  | 11       | 257  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 1.25       | 10%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL | 1.45       | 11.7% | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL | 1.34       | 4.6%  | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL | 1.23       | 3.8%  | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL | 0.956      | 12%   | <5       | <100 |
| Beckman UniCel DxI 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 1.14       | 12%   | <5       | <100 |
| (TT4)总T4                                                                                |       |            |       |          |      |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | ng/mL | 65.9       | 3.2%  | 8        | <100 |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL | 77.6       | 10.2% | 7        | 208  |
| Beckman UniCel DxI 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | ng/mL | 66.9       | 11.9% | 7        | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | ng/mL | 68         | 5%    | <5       | <100 |
| Abbott Alinity ci-series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )              | ng/mL | 58.1       | 3.4%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | ng/mL | 59.5       | 8%    | 21       | 459  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL | 62.7       | 7.5%  | 14       | 195  |

|         |       |       |     |            |     |         |       |         |
|---------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF     | 540-1 | 6×3mL | EXP | 2025-11-30 | LOT | 33828CF | 更新日期： | 2024-05 |
| 12 / 12 |       |       |     |            |     |         |       |         |

|                                                                                         |       |            |       |          |      |
|-----------------------------------------------------------------------------------------|-------|------------|-------|----------|------|
|                                                                                         |       | 批号：33828CF |       |          |      |
| 项目\仪器\试剂方法                                                                              | 单位    | 均值 \ CV    |       | Lab \ 点数 |      |
| Simens IM1600<br>Siemens \ Chemiluminescence Immunoassay (CLIA )                        | ng/mL | 64         | 5.5%  | 11       | 285  |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | ng/mL | 67.5       | 12%   | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | ng/mL | 63.4       | 3.3%  | <5       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | ng/mL | 63.4       | 10%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | ng/mL | 55.4       | 12%   | <5       | <100 |
| Roche cobas e 801<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)              | ng/mL | 50.5       | 3.6%  | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | ng/mL | 69.2       | 12%   | <5       | <100 |
| (VB12)维生素B12                                                                            |       |            |       |          |      |
| maglumi system<br>Snibe \ Electro-Chemiluminescence Immunoassay (ECLIA)                 | pg/mL | 884        | 7.1%  | 9        | 271  |
| YHLO iFlash Analyzer<br>YHL \ Chemiluminescence Immunoassay (CLIA )                     | pg/mL | 795        | 23.8% | 7        | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Electro-Chemiluminescence Immunoassay (ECLIA)     | pg/mL | 373        | 5.5%  | <5       | 104  |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | pg/mL | 477        | 5.4%  | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | pg/mL | 397        | 22.7% | 13       | 203  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | pg/mL | 319        | 7.5%  | 13       | 112  |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | pg/mL | 667        | 4.5%  | 12       | 177  |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | pg/mL | 380        | 8.9%  | 11       | 272  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | pg/mL | 509        | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | pg/mL | 561        | 11.9% | <5       | <100 |
| (β2M)β-2-微球蛋白                                                                           |       |            |       |          |      |
| Beckman AU series<br>Beckman Coulter \ Immunoturbidimetry (ITA)                         | ng/mL | 1171       | 4.7%  | 16       | 213  |
| Tosoh AIA 900/1800<br>TOSOH \ Chemiluminescence Immunoassay (CLIA )                     | ng/mL | 1170       | 12%   | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Roche \ Immunoturbidimetry (ITA)                          | ng/mL | 1020       | 10%   | <5       | <100 |
| (β-hCG)人绒毛膜促性腺激素β亚基                                                                     |       |            |       |          |      |
| Beckman AU series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | IU/mL | 0.004      | 5.6%  | 7        | 180  |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)  | IU/mL | 0.004      | 4.5%  | 7        | 102  |
| Roche cobas e411<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 0.003      | 5.9%  | <5       | <100 |
| Roche cobas e411<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)               | IU/mL | 0.003      | 5.5%  | <5       | <100 |
| Ortho Vitros 3600/5600<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 0.005      | 9.3%  | <5       | <100 |
| Roche cobas 4000/6000/8000<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/mL | 0.003      | 4.5%  | 14       | <100 |
| Autolumo A2000/A2000 PLUS<br>Autobio \ Chemiluminescence Immunoassay (CLIA )            | IU/mL | 0.003      | 7.3%  | 11       | <100 |
| Roche cobas E601/602<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)           | IU/mL | 0.003      | 10%   | <5       | <100 |
| maglumi system<br>Snibe \ Chemiluminescence Immunoassay (CLIA )                         | IU/mL | 0.004      | 2.5%  | <5       | <100 |
| Beckman UniCel Dxl 600/800<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )   | IU/mL | 0.004      | 12%   | <5       | <100 |
| Caris200<br>WanTai \ Chemiluminescence Immunoassay (CLIA )                              | IU/mL | 0.005      | 12%   | <5       | <100 |
| Siemens ADVIA Centaur XP/XPT<br>Siemens \ Chemiluminescence Immunoassay (CLIA )         | IU/mL | 0.004      | 11.9% | <5       | <100 |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )             | IU/mL | 0.003      | 12%   | <5       | <100 |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。