

REF	660-3	6×2mL	EXP	2024-02-20	LOT	3EF92B5	更新日期：	2024-05
1 / 2								

		批号：3EF92B5			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(HCT)血细胞比容/压积					
Maccura F580	%	51.2	2.8%	<5	<100
MAKER \ Electrical impedance					
Sysmex XN series	%	53.4	3.1%	<5	<100
sysmex \ Electrical impedance					
Sysmex XN series	%	60.2	1.7%	<5	<100
sysmex \ Flow cytometry					
OTHER	%	57.2	1%	<5	<100
DYMIND \ Calculate					
Sysmex XN series	%	34.5	1%	<5	<100
sysmex \ Calculate					
Sysmex XE series	%	52.6	9.9%	<5	<100
Dedicated \ Light scattering					
Siemens ADVIA 2120	%	51.9	10.2%	<5	<100
Dedicated \ Light scattering					
Beckman LH series	%	52.4	9.9%	<5	<100
Dedicated \ Electrical impedance					
Mindray BC-5600/BC-5800	%	51.8	10.2%	<5	<100
Dedicated \ Light scattering					
Abbott CELL-DYN 1800/3700	%	52.5	9.9%	<5	<100
Dedicated \ Light scattering					
(HGB/ Hb)血红蛋白					
Maccura F580	g/L	174	2.5%	<5	111
MAKER \ Colorimetry					
OTHER	g/L	181	0.9%	<5	<100
DYMIND \ Colorimetry					
Sysmex XN series	g/L	173	14.1%	10	221
sysmex \ Colorimetry					
Abbott CELL-DYN 1800/3700	g/L	167	10.2%	<5	<100
Dedicated \ Light scattering					
Beckman LH series	g/L	169	10.1%	<5	<100
Dedicated \ Light scattering					
Mindray BC-5600/BC-5800	g/L	167	10.2%	<5	<100
Dedicated \ Light scattering					
Sysmex XE series	g/L	165	10.3%	<5	<100
Dedicated \ Light scattering					
Siemens ADVIA 2120	g/L	168	10.1%	<5	<100
Dedicated \ Light scattering					
(MCH)平均红细胞血红蛋白含量					
Sysmex XN series	pg	30.5	3.8%	6	111
sysmex \ Calculate					
Sysmex XN series	pg	31.4	0.7%	<5	111
sysmex \ Electrical impedance					
Maccura F580	pg	31.7	1%	<5	111
MAKER \ Electrical impedance					
OTHER	pg	30.7	0.9%	<5	<100
DYMIND \ Calculate					
Beckman LH series	pg	31.8	10.1%	<5	<100
Dedicated \ Electrical impedance					
Abbott CELL-DYN 1800/3700	pg	32.9	9.4%	<5	<100
Dedicated \ Electrical impedance					
Sysmex XE series	pg	32.8	9.5%	<5	<100
Dedicated \ Light scattering					
Siemens ADVIA 2120	pg	32.4	9.9%	<5	<100
Dedicated \ Electrical impedance					
Mindray BC-5600/BC-5800	pg	32.4	9.6%	<5	<100
Dedicated \ Light scattering					
(MCHC)平均红细胞血红蛋白浓度					
Sysmex XN series	g/L	332	6.2%	6	110
sysmex \ Calculate					
Maccura F580	g/L	335	2.4%	<5	108
MAKER \ Electrical impedance					
Sysmex XN series	g/L	330	2.7%	<5	107
sysmex \ Electrical impedance					
OTHER	g/L	317	1%	<5	<100
DYMIND \ Calculate					
Mindray BC-5600/BC-5800	g/L	315	10.5%	<5	<100
Dedicated \ Light scattering					
Abbott CELL-DYN 1800/3700	g/L	310	10.3%	<5	<100
Dedicated \ Light scattering					
Siemens ADVIA 2120	g/L	315	10.2%	<5	<100
Dedicated \ Light scattering					
Sysmex XE series	g/L	310	10.6%	<5	<100
Dedicated \ Light scattering					
Beckman LH series	g/L	305	10.5%	<5	<100
Dedicated \ Electrical impedance					
(MCV)平均红细胞体积					
Sysmex XN series	fL	92.6	9.2%	6	109
sysmex \ Calculate					
Sysmex XN series	fL	94.8	2.6%	<5	<100
sysmex \ Electrical impedance					
Maccura F580	fL	94.8	2%	<5	111
MAKER \ Electrical impedance					
OTHER	fL	96.6	1.1%	<5	<100
DYMIND \ Electrical impedance					
Beckman LH series	fL	87	10.3%	<5	<100
Dedicated \ Light scattering					
Abbott CELL-DYN 1800/3700	fL	86	10.5%	<5	<100
Dedicated \ Light scattering					

REF	660-3	6×2mL	EXP	2024-02-20	LOT	3EF92B5	更新日期：	2024-05
2 / 2								

批号：3EF92B5					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Siemens ADVIA 2120 Dedicated \ Light scattering	fL	87	10.3%	<5	<100
Sysmex XE series Dedicated \ Light scattering	fL	88	10.2%	<5	<100
Mindray BC-5600/BC-5800 Dedicated \ Light scattering	fL	88	10.2%	<5	<100
(PLT)血小板					
Maccura F580 MAKER \ Electrical impedance	Giga/L (1E+9/L)	497	5.3%	<5	111
OTHER DYMIND \ Electrical impedance	Giga/L (1E+9/L)	486	4%	<5	<100
Sysmex XN series sysmex \ Electrical impedance	Giga/L (1E+9/L)	438	21.1%	10	222
Sysmex XE series Dedicated \ Light scattering	Giga/L (1E+9/L)	510	10.8%	<5	<100
Beckman LH series Dedicated \ Electrical impedance	Giga/L (1E+9/L)	505	10.9%	<5	<100
Mindray BC-5600/BC-5800 Dedicated \ Light scattering	Giga/L (1E+9/L)	500	11%	<5	<100
Siemens ADVIA 2120 Dedicated \ Electrical impedance	Giga/L (1E+9/L)	490	11.2%	<5	<100
Abbott CELL-DYN 1800/3700 Dedicated \ Electrical impedance	Giga/L (1E+9/L)	510	10.8%	<5	<100
(RBC)红细胞					
Sysmex XN series sysmex \ Electrical impedance	T/L (1E+12/L)	5.33	12.6%	6	151
Sysmex XN series sysmex \ Flow cytometry	T/L (1E+12/L)	6.09	1.7%	<5	<100
Maccura F580 MAKER \ Electrical impedance	T/L (1E+12/L)	5.47	2.4%	<5	110
OTHER DYMIND \ Electrical impedance	T/L (1E+12/L)	5.92	1.1%	<5	<100
Sysmex XE series Dedicated \ Light scattering	T/L (1E+12/L)	5.28	10.6%	<5	<100
Mindray BC-5600/BC-5800 Dedicated \ Light scattering	T/L (1E+12/L)	5.32	10.2%	<5	<100
Abbott CELL-DYN 1800/3700 Dedicated \ Electrical impedance	T/L (1E+12/L)	5.33	10.5%	<5	<100
Siemens ADVIA 2120 Dedicated \ Light scattering	T/L (1E+12/L)	5.34	10.1%	<5	<100
(WBC)白细胞					
Sysmex XN series sysmex \ Flow cytometry	Giga/L (1E+9/L)	17.4	1.3%	<5	<100
Sysmex XN series sysmex \ Electrical impedance	Giga/L (1E+9/L)	21.6	2.1%	<5	111
Maccura F580 MAKER \ Electrical impedance	Giga/L (1E+9/L)	21.4	2.9%	<5	111
OTHER DYMIND \ Electrical impedance	Giga/L (1E+9/L)	23.5	2.1%	<5	<100
Mindray BC-5600/BC-5800 Dedicated \ Light scattering	Giga/L (1E+9/L)	18.1	9.9%	<5	<100
Abbott CELL-DYN 1800/3700 Dedicated \ Electrical impedance	Giga/L (1E+9/L)	18.2	10.4%	<5	<100
Siemens ADVIA 2120 Dedicated \ Light scattering	Giga/L (1E+9/L)	17.9	10.6%	<5	<100
Beckman LH series Dedicated \ Electrical impedance	Giga/L (1E+9/L)	18.3	9.8%	<5	<100
Sysmex XE series Dedicated \ Light scattering	Giga/L (1E+9/L)	17.8	10.1%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。