

REF	671-1L	6×1mL	EXP	2024-05-26	LOT	400C53L	更新日期：	2024-04
1 / 2								

		批号：400C53L			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Werfen ACL series	Seconds	32.8	16.9%	<5	136
Werfen \ Magnetic Beads					
Sekisui Coapresta 3000	Seconds	23.2	32.8%	<5	124
SekiSui \ Light scattering					
Roche cobas T711	Seconds	37.4	9.9%	<5	<100
Roche \ clauss					
Sysmex CS 5100	Seconds	32.2	7.7%	<5	<100
Siemens \ Light scattering					
Sysmex CS 5100	Seconds	32.8	2.8%	<5	<100
sysmex \ Light scattering					
Werfen ACL series	Seconds	34.5	9.3%	<5	<100
Werfen \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	35.2	9.7%	<5	<100
Siemens \ Light scattering					
Mindray EXC810	Seconds	34.5	1.4%	<5	<100
Mindray \ Magnetic Beads					
Stago STA series	Seconds	36.2	8.6%	<5	<100
Dedicated \ Magnetic Beads					
Stago STA-R Evolution System	Seconds	36.8	2.5%	<5	<100
STA-Neoplastine CI \ Magnetic Beads					
Sysmex CA 1500/6000/7000	Seconds	29.3	10%	<5	<100
sysmex \ Light scattering					
(AT-III)抗凝血酶III					
Werfen ACL series	%	87.4	2.1%	<5	103
Werfen \ Nephelometry					
Roche cobas T711	%	89.8	9.9%	<5	<100
Roche \ Nephelometry					
Sysmex CS 5100	%	95	7.5%	<5	<100
Siemens \ Nephelometry					
Sekisui Coapresta 3000	%	171	56.5%	<5	<100
SekiSui \ Nephelometry					
Sekisui Coapresta 3000	g/L	0.273	10%	<5	<100
SekiSui \ Nephelometry					
Werfen ACL series	%	93.6	9.8%	<5	<100
Dedicated \ Nephelometry					
Sysmex CA 1500/6000/7000	%	89	10%	<5	<100
sysmex \ Nephelometry					
Sysmex CA 1500/6000/7000	%	93.8	9.8%	<5	<100
Siemens \ Nephelometry					
Stago STA series	%	104	5.3%	<5	<100
Stago \ Nephelometry					
(D-D(DDU)) D-二聚体(DDU)					
Sekisui Coapresta 3000	mg/L	0.72	50.1%	<5	119
SekiSui \ Nephelometry					
Werfen ACL series	mg/L	0.167	31.1%	<5	111
Werfen \ Nephelometry					
Sysmex CS 5100	mg/L	0.3	10%	<5	<100
Siemens \ Nephelometry					
Stago STA series	mg/L	0.425	3.5%	<5	<100
Stago \ Nephelometry					
Sysmex CS 5100	mg/L	0.31	7.6%	<5	<100
sysmex \ Nephelometry					
(D-D(FEU)) D-二聚体(FEU)					
Roche cobas T711	mg/L	0.3	9.9%	<5	<100
Roche \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	0.39	15.4%	<5	<100
SekiSui \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	0.24	10%	<5	<100
sysmex \ Nephelometry					
Sysmex CS 5100	mg/L	0.32	7.1%	<5	<100
Siemens \ Nephelometry					
(FDP)纤维蛋白(原)降解产物					
Sekisui Coapresta 3000	mg/L	1.85	9.4%	<5	<100
SekiSui \ Nephelometry					
Sysmex CS 5100	mg/L	2.65	11%	<5	<100
Siemens \ Nephelometry					
Roche cobas T711	mg/L	0.16	9.9%	<5	<100
Roche \ Nephelometry					
Werfen ACL series	mg/L	2.37	1.5%	<5	<100
Werfen \ Nephelometry					
Stago STA series	mg/L	2.34	1.1%	<5	<100
Stago \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	2	10%	<5	<100
sysmex \ Nephelometry					
Sysmex CA 1500/6000/7000	mg/L	2.62	12.2%	<5	<100
SekiSui \ Nephelometry					
(FIB)纤维蛋白原					
Werfen ACL series	g/L	3.68	9.8%	<5	126
Werfen \ Magnetic Beads					
Sekisui Coapresta 3000	g/L	3.01	18.9%	<5	101
SekiSui \ Light scattering					
Sysmex CS 5100	g/L	3.58	7%	<5	<100
Siemens \ Light scattering					
Roche cobas T711	g/L	3.61	9.9%	<5	<100
Roche \ clauss					
Sysmex CS 5100	g/L	3.34	5.4%	<5	<100
sysmex \ Light scattering					

REF	671-1L	6×1mL	EXP	2024-05-26	LOT	400C53L	更新日期：	2024-04
2 / 2								

		批号：400C53L			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Werfen ACL series	g/L	3.35	10.1%	<5	<100
Werfen \ Light scattering					
Stago STA series	g/L	3.68	10.3%	<5	<100
Dedicated \ Magnetic Beads					
Sysmex CA 1500/6000/7000	g/L	3.35	9.9%	<5	<100
Siemens \ Light scattering					
Stago STA-R Evolution System	g/L	3.13	3.6%	<5	<100
Stago \ Magnetic Beads					
Mindray EXC810	g/L	3.43	0.3%	<5	<100
Mindray \ Magnetic Beads					
Stago STA series	g/L	3.6	1.7%	<5	<100
Stago \ Magnetic Beads					
(PT INR)PT国际标准化比值					
Werfen ACL series	INR	1.08	6.8%	<5	<100
Werfen \ Calculate					
Roche cobas T711	INR	1.14	9.9%	<5	<100
Roche \ Calculate					
Sekisui Coapresta 3000	INR	1.14	3.1%	<5	<100
SekiSui \ Calculate					
Mindray EXC810	INR	0.99	1%	<5	<100
Mindray \ Calculate					
Stago STA-R Evolution System	INR	1.05	2.1%	<5	<100
STA-Neoplastine CI Plus \ Calculate					
Sysmex CS 5100	INR	1.08	1.6%	<5	<100
Siemens \ Calculate					
(PT)血浆凝血酶原时间					
Sekisui Coapresta 3000	Seconds	12	10.7%	<5	<100
SekiSui \ Light scattering					
Werfen ACL series	Seconds	11.7	15%	<5	167
Werfen \ Magnetic Beads					
Roche cobas T711	Seconds	10.2	9.9%	<5	<100
Roche \ clauss					
Sysmex CS 5100	Seconds	12.3	6.9%	<5	<100
Siemens \ Light scattering					
Stago STA series	Seconds	12.3	4.9%	<5	<100
Stago \ Magnetic Beads					
Stago STA series	Seconds	13.2	9.8%	<5	<100
Dedicated \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	12.3	10%	<5	<100
sysmex \ Light scattering					
Stago STA-R Evolution System	Seconds	13.5	2.3%	<5	<100
STA-Neoplastine CI Plus \ Magnetic Beads					
Sysmex CS 5100	Seconds	12.2	1.5%	<5	<100
sysmex \ Light scattering					
Werfen ACL series	Seconds	11.8	9.3%	<5	<100
Werfen \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	12.2	9.8%	<5	<100
Siemens \ Light scattering					
Mindray EXC810	Seconds	12	1.3%	<5	<100
Mindray \ Magnetic Beads					
(TT)凝血酶时间					
Sekisui Coapresta 3000	Seconds	16.9	7.5%	<5	<100
SekiSui \ Light scattering					
Werfen ACL series	Seconds	13.3	27.4%	<5	<100
Werfen \ Magnetic Beads					
Roche cobas T711	Seconds	19.3	9.9%	<5	<100
Roche \ clauss					
Sysmex CS 5100	Seconds	20.4	8.3%	<5	<100
Siemens \ Light scattering					
Mindray EXC810	Seconds	17.2	0.3%	<5	<100
Mindray \ Magnetic Beads					
Stago STA series	Seconds	16.4	9.8%	<5	<100
Werfen \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	19.2	9.9%	<5	<100
Dedicated \ Light scattering					
Stago STA series	Seconds	16.4	2.1%	<5	<100
Stago \ Magnetic Beads					
Sysmex CS 5100	Seconds	21.4	4.6%	<5	<100
sysmex \ Light scattering					
Sysmex CA 1500/6000/7000	Seconds	21	10%	<5	<100
sysmex \ Light scattering					

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。