

		批号：4000DOE			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex XN series sysmex \ Light scattering	Seconds	30	1.4%	6	229
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	31.2	2.5%	12	464
Stago STA series Stago \ Magnetic Beads	Seconds	33	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Dade Actin FS(APTT) \ Light scattering	Seconds	33	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Pathrombin SL(APTT) \ Light scattering	Seconds	37	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	31.1	12%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	%	101	12%	<5	<100
Werfen ACL series Werfen \ Nephelometry	%	98	12%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Werfen ACL series Werfen \ Nephelometry	mg/L	0.503	12%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex XN series sysmex \ Nephelometry	mg/L	1.63	7.5%	6	224
Sysmex CS 5100 Dedicated \ Nephelometry	mg/L	1.69	8.1%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	1.72	12%	<5	<100
Stago STA series Stago \ Nephelometry	mg/L	0.97	12%	<5	<100
(FIB)纤维蛋白原					
Sysmex XN series sysmex \ Light scattering	g/L	3.58	3.9%	6	226
Sysmex CS 5100 Dedicated \ Light scattering	g/L	3.38	6.1%	12	467
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	3.62	12%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	3.97	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	3.68	12%	<5	<100
(PT INR)PT国际标准化比值					
Sysmex XN series sysmex \ Calculate	INR	1.01	2.6%	6	229
Sysmex CS 5100 Dedicated \ Calculate	INR	1.04	3.8%	12	474
Werfen ACL series Werfen \ Calculate	INR	1.01	12%	<5	<100
Stago STA series Stago \ Calculate	INR	1.02	12%	<5	<100
Sysmex CS 5100 Siemens \ Calculate	INR	1.05	12%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex XN series sysmex \ Light scattering	Seconds	11.8	2.4%	6	230
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	12.2	3.8%	13	449
Stago STA series Stago \ Magnetic Beads	Seconds	14.4	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	13.4	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	11.6	12%	<5	<100
(TT)凝血酶时间					
Sysmex XN series sysmex \ Light scattering	Seconds	19.1	2.6%	6	229
Sysmex CS 5100 Dedicated \ Light scattering	Seconds	20.3	4%	12	470
Sysmex CS 5100 sysmex \ Light scattering	Seconds	19.4	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。