

REF	671-1	6×1mL	EXP	2025-12-27	LOT	4000D72	更新日期：	2024-04
1 / 3								

		批号：4000D72			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CS 5100 Siemens \ OTHER	Seconds	33	2.7%	6	119
Werfen ACL series Werfen \ Magnetic Beads	Seconds	35.4	4.9%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	31.6	30.3%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	69.8	10%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	32.8	2.4%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Dade Actin(APTT) \ Light scattering	Seconds	34	10%	<5	<100
Sysmex CN Series sysmex \ Light scattering	Seconds	29.5	0.9%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	33	10%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Light scattering	Seconds	28	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	Seconds	36.4	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Pathrombin SL(APTT) \ Light scattering	Seconds	37	10%	<5	<100
Sysmex CS-2400 sysmex \ clauss	Seconds	31.5	2.9%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ Light scattering	Seconds	32	10%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ clauss	Seconds	30	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	32.6	0%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Dade Actin FS(APTT) \ Light scattering	Seconds	33	10%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CS 5100 Siemens \ OTHER	%	66.8	15.6%	6	146
Sekisui Coapresta 3000 SekiSui \ Nephelometry	%	210	55.5%	<5	<100
Werfen ACL series Werfen \ Nephelometry	%	94.8	10.7%	<5	<100
Werfen ACL series Werfen IL Antithrombin III Chromogenic (AT□) \ Nephelometry	%	97.7	10%	<5	<100
Sysmex CN Series sysmex \ Nephelometry	%	92.8	1.4%	<5	<100
Werfen ACL series Werfen Hemosil Liquid Antithrombin (AT□) \ Nephelometry	%	101	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Berichrom Antithrombin III (AT□) \ Nephelometry	%	95	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Innovance Antithrombin (AT□) \ Nephelometry	%	98	10%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	%	26	0.6%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Sysmex CS 5100 Siemens \ OTHER	mg/L	1.76	3.1%	6	123
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	1.72	18.5%	<5	109
Werfen ACL series Werfen Hemosil D-Dimer HS \ Nephelometry	mg/L	1.15	57.6%	<5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	1.07	42.6%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	mg/L	1.41	1.4%	<5	<100
Dirui BCA Series DIRUI \ Nephelometry	mg/L	0.88	10%	<5	<100
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	1.99	0.5%	<5	<100
Werfen ACL series Wenfen Hemosil D-Dimer \ Nephelometry	mg/L	0.538	10%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	2.56	10%	<5	<100
Sysmex CN Series sysmex \ Nephelometry	mg/L	0.602	1.2%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CA 1500/6000/7000 Siemens Innovance Antithrombin (AT□) \ Nephelometry	mg/L	1.55	10%	<5	<100
Mindray C Series CHANGDAO \ Nephelometry	mg/L	1.57	10.2%	<5	<100
Sysmex CS-2400 sysmex \ Nephelometry	mg/L	1.62	2.8%	<5	<100
Stago STA series Stago \ Nephelometry	mg/L	0.88	10%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Sysmex CS 5100 Siemens \ OTHER	mg/L	4.78	17.3%	6	142
Werfen ACL series Werfen \ Nephelometry	mg/L	3.68	12.8%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	3.03	14.7%	<5	<100

REF	671-1	6×1mL	EXP	2025-12-27	LOT	4000D72	更新日期：	2024-04
2 / 3								

		批号：4000D72			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	mg/L	5.4	1.4%	<5	<100
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	4000	2.5%	<5	<100
Sysmex CN Series sysmex \ Nephelometry	mg/L	3.44	6.3%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Nephelometry	mg/L	7.87	10%	<5	<100
Sysmex CS 5100 sysmex \ OTHER	mg/L	5.04	8.8%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	15.6	10%	<5	<100
(FIB)纤维蛋白原					
Sysmex CS 5100 Siemens \ OTHER	g/L	3.31	5.9%	6	140
Sekisui Coapresta 3000 SekiSui \ Light scattering	g/L	1.53	93.7%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	g/L	1.67	94.8%	<5	<100
Werfen ACL series Werfen Hemosil QFA Thrombin (FIB) \ clauss	g/L	3.53	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	g/L	3.5	10%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	g/L	3.51	3.1%	<5	<100
Werfen ACL series Werfen Hemosil Fib-C \ Light scattering	g/L	3.37	10%	<5	<100
Sysmex CN Series sysmex \ Light scattering	g/L	3.46	2.8%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Thrombin (FIB) \ Light scattering	g/L	3.3	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Multifibren U (FIB) \ Light scattering	g/L	3.58	10%	<5	<100
Werfen ACL series Werfen Hemosil QFA Thrombin (FIB) \ Light scattering	g/L	3.69	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	3.63	10%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	g/L	2.02	10%	<5	<100
(PT INR)PT国际标准化比值					
Sysmex CS 5100 Siemens \ Calculate	INR	1.04	2.9%	6	144
Dirui BCA Series DIRUI \ Calculate	INR	1	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Calculate	INR	1.24	9.9%	<5	<100
Sysmex CS 5100 sysmex \ Calculate	INR	0.955	0.5%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex CS 5100 Siemens \ OTHER	Seconds	12	2.7%	6	145
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	13	12.4%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	10.9	7.5%	<5	<100
Werfen ACL series Werfen Hemosil PT-Fib HS Plus \ Light scattering	INR	0.93	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ Light scattering	INR	0.91	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	11.2	0.4%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Light scattering	Seconds	10.2	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	INR	0.91	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Thromborel S (PT) \ Light scattering	INR	0.93	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	Seconds	12	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Thromborel S (PT) \ Light scattering	Seconds	11.7	10%	<5	<100
Werfen ACL series Werfen Hemosil PT-Fib HS Plus \ Light scattering	Seconds	14	10%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	11.4	1.1%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Light scattering	INR	0.91	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ clauss	Seconds	11.6	10%	<5	<100
Sysmex CS-2400 sysmex \ clauss	Seconds	11.8	1.8%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	12.6	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ Light scattering	Seconds	10.3	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens PT Innovin \ Light scattering	INR	1.01	10%	<5	<100
(TT)凝血酶时间					
Sysmex CS 5100 Siemens \ OTHER	Seconds	20.6	4.2%	6	121

REF

671-1

6×1mL

EXP

2025-12-27

LOT

4000D72

更新日期：2024-04

3 / 3

		批号：4000D72			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Werfen ACL series	Seconds	13.9	4.4%	<5	<100
Werfen \ Magnetic Beads	Seconds	19.7	21.1%	<5	<100
Sekisui Coapresta 3000	Seconds	18	0.8%	<5	<100
SekiSui \ Light scattering	Seconds	17.4	10%	<5	<100
Sysmex CS 5100	Seconds	19.5	3.2%	<5	<100
sysmex \ Light scattering	Seconds	17.6	10%	<5	<100
Dirui BCA Series	Seconds	11.5	2%	<5	<100
DIRUI \ clauss	Seconds	13.4	10%	<5	<100
Sysmex CN Series	Seconds	16.6	10%	<5	<100
sysmex \ Light scattering	Seconds				
Werfen ACL series	Seconds				
Werfen Hemosil SynthAsil APTT \ clauss	Seconds				
Sysmex CA 1500/6000/7000	Seconds				
sysmex \ Light scattering	Seconds				
Werfen ACL series	Seconds				
Werfen \ Light scattering	Seconds				
Sysmex CS 5100	Seconds				
Siemens \ Light scattering	Seconds				

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。