

质控品名： 尿液定量生化液体非定值质控品

REF	586-1	6×8mL	EXP	2025-10-20	LOT	37E59F2	更新日期：	2024-05
1 / 2								

批号：37E59F2					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(AMY)淀粉酶					
Roche Hitachi/Hitachi Systems DiaSys \ PNP-G7	U/L	31.9	1.9%	<5	114
Beckman AU series Szybio \ Maltotetraose-rate	U/L	35.8	4.9%	10	249
Roche Hitachi/Hitachi Systems Roche \ PNP-G7	U/L	39	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ PNP-G7	U/L	41.3	9.9%	<5	<100
(Ca)钙					
Roche Hitachi/Hitachi Systems DiaSys \ Arsenazo III	mmol/L	1.76	2.2%	<5	112
Roche Hitachi/Hitachi Systems Roche \ Arsenazo III	mmol/L	1.77	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ BAPTA	mmol/L	1.91	12.1%	<5	<100
(Cl)氯					
Roche Hitachi/Hitachi Systems Hitachi \ ISE (indirect)	mmol/L	69	2.1%	<5	113
Roche cobas 4000/6000/8000 Roche \ ISE (indirect)	mmol/L	77.6	9.9%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ ISE (indirect)	mmol/L	75	10%	<5	<100
(Cortisol)皮质醇					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	23.6	11.8%	<5	<100
(Cr)肌酐					
Beckman AU series OTHER \ OTHER	μmol/L	7485	8%	9	210
Roche Hitachi/Hitachi Systems DiaSys \ Enzymatic	μmol/L	6523	1.6%	<5	111
Roche Hitachi/Hitachi Systems Roche \ Enzymatic	μmol/L	6371	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Enzymatic	μmol/L	5914	12.3%	<5	<100
(GLU)葡萄糖					
Roche Hitachi/Hitachi Systems DiaSys \ Glucose oxidase (GOD)	mmol/L	1.19	2.7%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Hexokinase	mmol/L	1.31	12.2%	<5	<100
(hCG)人绒毛膜促性腺激素					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	0.012	10%	<5	<100
(K)钾					
Roche Hitachi/Hitachi Systems Hitachi \ ISE (indirect)	mmol/L	28.5	1.3%	<5	114
Roche cobas 4000/6000/8000 Roche \ ISE (indirect)	mmol/L	29.2	9.9%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ ISE (indirect)	mmol/L	27.9	10%	<5	<100
(MA)微量白蛋白					
Beckman AU series OTHER \ OTHER	mg/L	14.1	23.7%	9	250
Roche Hitachi/Hitachi Systems Roche \ Immunoturbidimetry (ITA)	mg/L	27.2	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Bromocresol green	mg/L	24.3	12.3%	<5	<100
(Mg)镁					
Roche Hitachi/Hitachi Systems DiaSys \ Xylidyl Blue	mmol/L	2.65	6.9%	<5	108
Roche cobas 4000/6000/8000 Roche \ Xylidyl Blue	mmol/L	2.45	12.1%	<5	<100
(Na)钠					
Roche Hitachi/Hitachi Systems Hitachi \ ISE (indirect)	mmol/L	74.6	6%	<5	114
Roche cobas c 501/502/701/702 Roche \ ISE (direct)	mmol/L	73.2	10%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ ISE (indirect)	mmol/L	77.1	10%	<5	<100
(Osm)摩尔渗透压浓度					
WESCOR Osm WESCOR \ Freezing Point Osmometer	mOSM/kg	439	12%	<5	<100
(P)磷					
Roche Hitachi/Hitachi Systems DiaSys \ Phosphomolybdate-UV-end point	mmol/L	8.06	2.2%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Phosphomolybdate UV-rate	mmol/L	8.49	12.8%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ Phosphomolybdate UV-rate	mmol/L	8.43	10%	<5	<100

REF	586-1	6×8mL	EXP	2025-10-20	LOT	37E59F2	更新日期：	2024-05
2 / 2								

		批号：37E59F2			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(pH)酸碱度					
HACH pH Meter OTHER \ Electrode	No Units	5.93	12.7%	<5	<100
(PRO/ TP)(总) 蛋白					
Abbott ARCHITECT i series LEADMAN \ Biuret-endpoint	g/L	0.12	6%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ Biuret-endpoint	g/L	0.137	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Biuret-endpoint	g/L	0.136	12.7%	<5	<100
(SG)比重					
OTHER OTHER \ Dipstick	No Units	1.01	12%	<5	<100
(UA)尿酸					
Roche Hitachi/Hitachi Systems DiaSys \ Uricase colorimetry	μmol/L	449	8.4%	<5	<100
Roche Hitachi/Hitachi Systems Roche \ Uricase colorimetry	μmol/L	383	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Uricase colorimetry	μmol/L	463	13.1%	<5	<100
(Urea)尿素					
Roche Hitachi/Hitachi Systems Roche \ Urease-UV-rate	mmol/L	128	10%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Urease colorimetry	mmol/L	67.4	12.3%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。