

质控品名： 特殊生化与蛋白非定值质控品

REF	583-2	6×3mL	EXP	2024-04-28	LOT	379ECBA	更新日期：	2024-05
1 / 2								

批号：379ECBA					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(5'NT)5`核苷酸酶					
Roche cobas c 501/502/701/702	U/L	34.5	1.7%	<5	<100
Roche \ Enzymatic					
Beckman AU series	U/L	41	12%	<5	<100
Diazyme \ Enzymatic					
(ADA)腺苷脱氨酶					
Roche Hitachi/Hitachi Systems	U/L	150	1.9%	34	719
LEADMAN \ Peroxidase					
Roche cobas c 501/502/701/702	U/L	140	1.2%	<5	<100
Roche \ Trinder reaction					
Roche cobas 4000/6000/8000	U/L	140	1.2%	17	405
BSBE \ Peroxidase					
Beckman AU series	U/L	157	1.8%	15	804
Reebio \ Peroxidase					
Abbott ARCHITECT c series	U/L	149	2.4%	14	285
BSBE \ Peroxidase					
Beckman AU series	U/L	145	2.9%	13	359
KINGSBIO \ OTHER					
Siemens ADVIA 1800/2400	U/L	145	5.4%	12	364
MedicalSystem \ Trinder reaction					
Roche Hitachi/Hitachi Systems	U/L	143	1.9%	10	299
MedicalSystem \ Trinder reaction					
Beckman AU series	U/L	156	12%	<5	<100
Diazyme \ Peroxidase					
(AFU)α-L-岩藻糖苷酶					
Beckman AU series	U/L	33.4	6.2%	9	247
BSBE \ CNPF					
Roche Hitachi/Hitachi Systems	U/L	30.8	6.7%	34	729
BSBE \ CNPF					
Roche cobas c 501/502/701/702	U/L	27.8	4.4%	<5	<100
Roche \ CNPF					
Beckman AU series	U/L	23.4	3.7%	16	804
Reebio \ OTHER					
Beckman AU series	U/L	19.6	2.4%	13	361
chongqiong zybio \ OTHER					
Siemens ADVIA 1800/2400	U/L	35.5	6%	12	364
XJKC \ CNPF					
Beckman AU series	U/L	29	12%	<5	<100
Diazyme \ CNPF					
(anti CCP)抗环瓜氨酸抗体					
Beckman AU series	U/mL	26.6	26.8%	9	192
OTHER \ Immunoturbidimetry (ITA)					
YHLO iFlash Analyzer	U/mL	9.61	4.7%	5	122
YHL \ Chemiluminescence Immunoassay (CLIA)					
Abbott ARCHITECT i series	U/mL	13.3	5.2%	17	380
Abbott \ Chemiluminescence Immunoassay (CLIA)					
OTHER	U/mL	31.8	4.1%	<5	<100
MAKER \ immunoturbidimetric					
Roche cobas c 501/502/701/702	U/mL	34	12%	<5	<100
Roche \ immunoturbidimetric					
Siemens Immulite1000/2000	U/mL	46	4.7%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
(Apo-A I)载脂蛋白A I					
Abbott ARCHITECT c series	g/L	1.46	2.3%	8	162
BSBE \ immunoturbidimetric					
(C1q)补体C1q					
Roche cobas c 501/502/701/702	mg/L	231	3.2%	<5	<100
Roche \ Immunoturbidimetry (ITA)					
Roche Hitachi/Hitachi Systems	mg/L	237	6.1%	10	298
BeiJia(Shanghai) \ Immunoturbidimetry (ITA)					
(CO2 (HCO3-))二氧化碳（碳酸氢根）					
TOSHIBA TBA	mmol/L	28	6.6%	8	218
OTHER \ OTHER					
Abbott ARCHITECT c series	mmol/L	25.7	4.7%	14	311
BSBE \ OTHER					
Beckman AU series	mmol/L	23.6	5.4%	13	340
Ausbio \ OTHER					
Beckman AU series	mmol/L	20	5.1%	11	614
Beckman Coulter \ OTHER					
Roche cobas c 501/502/701/702	mmol/L	27.6	10%	<5	<100
Roche \ Colorimetry					
Beckman AU series	mmol/L	32	12%	<5	<100
OTHER \ PEPC					
(Cys C)半胱氨酸蛋白酶抑制剂C/ 胱抑素C					
Beckman AU series	mg/L	1.92	2.6%	9	253
BSBE \ Immunoturbidimetry (ITA)					
TOSHIBA TBA	mg/L	1.95	4.4%	6	150
OTHER \ Immunoturbidimetry (ITA)					
Roche Hitachi/Hitachi Systems	mg/L	1.94	4.7%	34	723
BSBE \ Immunoturbidimetry (ITA)					
Siemens ADVIA 1800/2400	mg/L	1.91	5.8%	24	704
MAKER \ Immunoturbidimetry (ITA)					
Beckman AU series	mg/L	2.13	3.9%	15	809
Beckman Coulter \ OTHER					
Abbott ARCHITECT c series	mg/L	2	1.6%	14	277
BSBE \ OTHER					

REF	583-2	6×3mL	EXP	2024-04-28	LOT	379ECBA	更新日期：	2024-05
2 / 2								

		批号：379ECBA			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Beckman AU series chongqiong zybio \ OTHER	mg/L	2.29	2%	13	360
Roche cobas c 501/502/701/702 Merit Choice \ Immunoturbidimetry (ITA)	mg/L	1.9	10%	<5	<100
Beckman AU series KaniyaBiomedical \ Immunoturbidimetry (ITA)	mg/L	1.87	12%	<5	<100
(FRUC (GSP))果糖胺 (糖化血清蛋白)					
Beckman AU series BSBE \ NBT	mmol/L	0.373	3.4%	9	253
Beckman AU series Reebio \ NBT	mmol/L	3.23	12.1%	<5	<100
Roche cobas c 501/502/701/702 Roche \ NBT	mmol/L	0.404	12%	<5	<100
(LAP)亮氨酸氨基肽酶					
Roche cobas c 501/502/701/702 OTHER \ OTHER	U/L	27.3	8.4%	<5	<100
Roche Hitachi/Hitachi Systems OTHER \ OTHER	U/L	75.9	3.2%	10	299
(Lp(a))脂蛋白a					
TOSHIBA TBA OTHER \ Immunoturbidimetry (ITA)	mg/L	317	3.8%	6	170
Beckman AU series Beckman Coulter \ OTHER	mg/L	336	3.2%	16	829
Beckman AU series KINGSBIO \ Immunoturbidimetry (ITA)	mg/L	356	2.1%	13	361
Siemens ADVIA 1800/2400 SekiSui \ Immunonephelometry (INA)	mg/L	161	14.1%	12	365
Roche cobas c 501/502/701/702 SekiSui \ Immunoturbidimetry (ITA)	mg/L	227	10%	<5	<100
Beckman AU series KaniyaBiomedical \ Immunonephelometry (INA)	mg/L	326	12%	<5	<100
(NAG)乙酰氨基葡萄糖苷酶					
Roche cobas c 501/502/701/702 Roche \ Enzymatic	No Units	275	2.4%	<5	<100
Roche Hitachi/Hitachi Systems BSBE \ Enzymatic	No Units	223	5.8%	10	297
(TRF)转铁蛋白					
Abbott ARCHITECT c series BSBE \ OTHER	mg/dL	260	1.7%	8	156
(UIBC)不饱和铁结合力					
Roche cobas c 501/502/701/702 Roche \ Ferene	μmol/L	718	5%	<5	<100
(β-Hb)Beta-羟丁酸					
Beckman AU series BSBE \ Dehydrogenase	mmol/L	0.96	3.2%	9	254
Roche cobas c 501/502/701/702 Roche \ Dehydrogenase	mmol/L	1.08	2.4%	<5	<100
Siemens ADVIA 1800/2400 DiaSys \ Dehydrogenase	mmol/L	0.979	13.1%	12	365
Roche Hitachi/Hitachi Systems DiaSys \ Dehydrogenase	mmol/L	0.802	4.5%	10	297
Beckman AU series Beckman Coulter \ Dehydrogenase	mmol/L	1.08	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。