

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
1 / 6								

批号：336CA00					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(25-OH-Vit D)25羟基维生素D					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	43.7	10%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	109	10%	<5	<100
(AFP)甲胎蛋白					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	115	7.1%	6	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	126	6.6%	6	115
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	145	4.3%	5	153
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	136	6.3%	5	141
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	135	4.2%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	134	2.8%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	3.2%	12	306
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	153	4.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	10%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	96.8	2.7%	<5	<100
(ALD)醛固酮					
Agilent aglient \ Chemiluminescence Immunoassay (CLIA)	pg/mL	179	10%	<5	<100
(Anti-Tg)抗甲状腺球蛋白抗体					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	6.7	19.7%	7	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.8	5.4%	5	152
Mindray CL-2000i/1000i/6000i series Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	7.44	16.3%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50.5	4.4%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	24.4	5.8%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	75.7	1.6%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	61.3	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/mL	123	3.3%	<5	<100
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	IU/mL	72.8	2.3%	<5	<100
(Anti-TPO)抗甲状腺过氧化物酶抗体					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	33	13.4%	7	194
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.4	6.4%	6	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/mL	44.6	4.7%	5	152
Mindray CL-2000i/1000i/6000i series Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	33.4	14.8%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	45.2	5%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	39.5	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/mL	43.9	4.3%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	32.1	2.3%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.6	7.5%	<5	<100
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	IU/mL	33.8	3.7%	<5	<100
(ASD)雄烯二酮					
Agilent aglient \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.29	10%	<5	<100
(CA 125)糖类抗原 125					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.6	20.3%	7	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	42.1	6.7%	6	106
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/mL	64.9	2.9%	5	153
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	43.3	2.8%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
2 / 6								

		批号：336CA00			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas e 801	IU/mL	38.2	2.7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	75.7	3%	12	306
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Maccura i3000	IU/mL	73.6	2%	<5	<100
MAKER \ Chemiluminescence Immunoassay (CLIA)					
Lifotronic ecl Series	IU/mL	52	1.7%	<5	<100
Lifotronic \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	79.3	10%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(CA 15-3)糖类抗原 15-3					
Beckman UniCel Dxl 600/800	IU/mL	21.4	21%	7	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	54.9	8.8%	6	124
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
OTHER	IU/mL	71.4	2.7%	5	152
DIRUI \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	IU/mL	54.1	7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e 801	IU/mL	56.4	2.6%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	26.8	5.7%	12	300
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Maccura i3000	IU/mL	54.6	2%	<5	<100
MAKER \ Chemiluminescence Immunoassay (CLIA)					
Lifotronic ecl Series	IU/mL	45.5	1.9%	<5	<100
Lifotronic \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	61.4	10%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(CA 19-9)糖类抗原 19-9					
Beckman UniCel Dxl 600/800	IU/mL	79	20.7%	7	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
OTHER	IU/mL	21.8	3.8%	5	152
DIRUI \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	IU/mL	49.8	3.6%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Mindray CL-2000i/1000i/6000i series	IU/mL	40.4	6.7%	12	291
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Lifotronic ecl Series	IU/mL	43.6	3.9%	<5	<100
Lifotronic \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	45.1	2.7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Maccura i3000	IU/mL	34.2	1.5%	<5	<100
MAKER \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	IU/mL	38.5	10%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(CEA)癌胚抗原					
Beckman UniCel Dxl 600/800	ng/mL	10.2	19.6%	7	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	10.3	5%	6	106
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
OTHER	ng/mL	10.5	5.5%	5	152
DIRUI \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	ng/mL	9.95	7.3%	5	142
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	ng/mL	10.1	4.2%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e 801	ng/mL	10.1	2.4%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Mindray CL-2000i/1000i/6000i series	ng/mL	7.43	9.8%	12	307
Mindray \ Chemiluminescence Immunoassay (CLIA)					
Maccura i3000	ng/mL	18.8	3.1%	<5	<100
MAKER \ Chemiluminescence Immunoassay (CLIA)					
Tosoh AIA 900/1800	ng/mL	20.2	10%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
Lifotronic ecl Series	ng/mL	8.14	4.4%	<5	<100
Lifotronic \ Chemiluminescence Immunoassay (CLIA)					
(CK-MB)肌酸激酶-MB (活性)					
Tosoh AIA 900/1800	ng/mL	15.6	10%	<5	<100
TOSOH \ Chemiluminescence Immunoassay (CLIA)					
(Cortisol)皮质醇					
Tosoh AIA 900/1800	ng/mL	211	10.3%	<5	<100
TOSOH \ Enzyme immunoassay (EIA)					
(C-P)C-肽					
Roche cobas 4000/6000/8000	ng/mL	2.25	1.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e411	ng/mL	2.38	11.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Lifotronic ecl Series	ng/mL	2.59	8.3%	<5	<100
Lifotronic \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	2.12	9.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Maccura i3000	ng/mL	2.2	4.1%	<5	<100
MAKER \ Chemiluminescence Immunoassay (CLIA)					

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
3 / 6								

批号：336CA00					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(DGXN)地高辛					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.89	9.9%	<5	<100
(DHEA-S)硫酸脱氢表雄酮					
Tosoh AIA 900/1800 TOSOH \ Enzyme immunoassay (EIA)	ng/mL	5711	10%	<5	<100
(E2)雌二醇					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	608	55%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1233	2.2%	5	152
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	811	3%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	896	5%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	941	10%	<5	<100
(EPO)促红细胞生成素					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	50	10%	<5	<100
(FA)叶酸					
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.7	3.8%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.19	2.8%	<5	<100
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.89	10%	<5	<100
(Ferritin)铁蛋白					
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	255	5.7%	5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	182	10%	<5	<100
(F-PSA)游离前列腺特异性抗原					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.69	13.9%	6	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.11	5%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.16	3.1%	12	303
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.81	2.3%	<5	<100
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.21	4.9%	<5	<100
(FSH)促卵泡刺激素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.5	4.9%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.9	4.8%	5	152
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	41.4	3.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	30	10%	<5	<100
(FT3)游离T3					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	8.6	18.2%	7	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	12.2	6.7%	7	192
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.2	1.4%	5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	13.2	2.7%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.5	3.1%	5	152
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	10.4	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.38	7%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.9	3.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	22.4	1.1%	<5	<100
(FT4)游离T4					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33.2	8.2%	7	190
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.4	5.7%	6	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	42.5	2.5%	5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
4 / 6								

		批号：336CA00			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	34.9	3.7%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	pmol/L	43.5	4.5%	5	152
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	31.9	2.6%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	pmol/L	45.7	2.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.4	6.5%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	32.8	10%	<5	<100
(hCG)人绒毛膜促性腺激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	199	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	124	9.7%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	138	1.2%	<5	<100
(HCY)同型半胱氨酸					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μmol/L	21.8	0.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Enzyme immunoassay (EIA)	μmol/L	21.5	10%	<5	<100
(HGH)人生长激素					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	3.77	10%	<5	<100
(HS-CTnI) (高) 超敏肌钙蛋白I					
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.309	10%	<5	<100
(IgE)免疫球蛋白 IgE					
Tosoh AIA 900/1800 TOSOH \ Enzyme linked immunosorbent assays (ELISA)	IU/mL	83.8	10%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	77.3	1.4%	<5	<100
(IGF1)胰岛素生长因子1					
Siemens Immulite1000/2000 Siemens \ Chemiluminescence Immunoassay (CLIA)	nmol/L	18	10%	<5	<100
(INS)胰岛素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.7	3.1%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	31	10%	<5	<100
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	mIU/L	33	3.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.8	1.5%	<5	<100
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	mIU/L	28.3	2.5%	<5	<100
(LH)促黄体生成素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.2	2.1%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	IU/L	22.3	4.3%	5	152
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	17.4	10%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	IU/L	58.2	1.4%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	IU/L	15.4	7.8%	<5	<100
(Myo / Mb)肌红蛋白					
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	147	10%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1020	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	120	10%	<5	<100
(PRL)催/ 泌乳素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	14.7	6%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.5	3.6%	5	152
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	26	4.4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.9	10%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.3	4.4%	<5	<100

质控品名： 多专业免疫分析非定值质控品

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
5 / 6								

		批号：336CA00			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(PROG)孕酮					
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.4	5.1%	5	152
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.9	6.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	3.7%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	21.8	2.3%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	11.4	10%	<5	<100
(PSA)前列腺特异性抗原					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.1	4.9%	6	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.6	4.3%	5	152
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.2	5.6%	<5	<100
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.7	4.5%	12	301
Maccura i3000 MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.8	2.4%	<5	<100
Lifotronic ecl Series Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	13.7	4%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	12.6	10%	<5	<100
(PTH)甲状旁腺素					
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.96	1.2%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	8.85	10.1%	<5	<100
(TES)睾酮					
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.429	3.8%	5	152
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.458	6.2%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.54	3.8%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	0.479	10.2%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.452	5.9%	<5	<100
(Tg)甲状腺球蛋白					
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	60.8	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.108	13.6%	<5	<100
(THEO)茶碱					
Roche cobas c 501/502/701/702 Roche \ Chemiluminescence Immunoassay (CLIA)	µg/mL	23.1	10%	<5	<100
(TSH)促甲状腺素					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.47	19.6%	7	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.1	5.8%	7	194
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.7	5.1%	6	117
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	11.5	2.8%	5	<100
OTHER DIRUI \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.6	4.4%	5	152
Mindray CL-2000i/1000i/6000i series Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.2	4.7%	<5	<100
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.54	8.4%	<5	<100
Siemens ADVIA Centaur XP/XPT Siemens \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.1	1.6%	<5	<100
Tosoh AIA 900/1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.2	10%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	10.5	1.4%	<5	<100
(TT3)总三碘甲状腺原氨酸					
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.68	24%	7	<100
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.45	7.1%	7	190
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.17	1.9%	5	<100
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.49	7.6%	5	<100

REF	539-2	6×3mL	EXP	2026-04-10	LOT	336CA00	更新日期：	2024-05
6 / 6								

批号：336CA00					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
OTHER					
DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.06	4.6%	5	153
Mindray CL-2000i/1000i/6000i series					
Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.002	5.1%	<5	<100
Abbott ARCHITECT i series					
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.02	9.2%	<5	<100
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.36	1%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.25	2.7%	<5	<100
(TT4)总T4					
Beckman UniCel DxI 600/800					
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	112	26.6%	7	<100
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	133	3.3%	6	182
Roche cobas 4000/6000/8000					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	129	6%	6	110
Mindray CL-2000i/1000i/6000i series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	136	3.4%	5	<100
OTHER					
DIRUI \ Chemiluminescence Immunoassay (CLIA)	ng/mL	133	3.7%	5	152
Mindray CL-2000i/1000i/6000i series					
Mindray \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.4	5%	<5	<100
Abbott ARCHITECT i series					
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	142	8.4%	<5	<100
Siemens ADVIA Centaur XP/XPT					
Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	181	3.3%	<5	<100
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	134	1.3%	<5	<100
(VB12)维生素B12					
Lifotronic ecl Series					
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	pg/mL	769	1.7%	<5	<100
Tosoh AIA 900/1800					
TOSOH \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	698	10%	<5	<100
Maccura i3000					
MAKER \ Chemiluminescence Immunoassay (CLIA)	pg/mL	630	1.2%	<5	<100
(β2M)β-2-微球蛋白					
Maccura i3000					
MAKER \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1524	2.4%	<5	<100
Tosoh AIA 900/1800					
TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1520	9.9%	<5	<100
Beckman UniCel DxI 600/800					
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1534	2.1%	<5	<100
Lifotronic ecl Series					
Lifotronic \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1258	4.1%	<5	<100
(β-hCG)人绒毛膜促性腺激素β亚基					
Mindray CL-2000i/1000i/6000i series					
Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.114	4.6%	5	<100
Maccura i3000					
MAKER \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.131	8.7%	<5	<100
Roche cobas 4000/6000/8000					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.152	1.7%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。