

REF	588-1	6×3.5mL	EXP	2023-06-16	LOT	381666C	更新日期：	2024-05
1 / 1								

批号：381666C					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Siemens ADVIA Centaur XP/XPT JingYuan(Shanghai) \ Bromocresol green	g/L	0.248	0.5%	<5	<100
(Cl)氯					
Beckman AU series Dedicated \ ISE (direct)	mmol/L	103	1.6%	8	156
Beckman AU series Beckman Coulter \ ISE (indirect)	mmol/L	166	3.1%	12	380
Mindray BS series Mindray \ ISE (indirect)	mmol/L	151	2.6%	12	286
Siemens ADVIA Centaur XP/XPT Siemens \ ISE (indirect)	mmol/L	148	1.9%	<5	<100
Roche cobas 4000/6000/8000 Roche \ ISE (indirect)	mmol/L	132	12%	<5	<100
Roche cobas c 501/502/701/702 Roche \ ISE (indirect)	mmol/L	120	3.1%	<5	<100
(GLU)葡萄糖					
Beckman AU series Dedicated \ Hexokinase	mmol/L	1.83	5.5%	8	139
Mindray BS series Mindray \ Glucose oxidase (GOD)	mmol/L	2.75	3.1%	12	288
Beckman AU series BSBE \ Glucose oxidase (GOD)	mmol/L	2.43	4%	10	289
Siemens ADVIA Centaur XP/XPT Sansure Biotech \ Glucose oxidase (GOD)	mmol/L	2.53	1.1%	<5	<100
Roche cobas 4000/6000/8000 Roche \ Hexokinase	mmol/L	2.57	12%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Hexokinase	mmol/L	2.47	1%	<5	<100
(LA)乳酸					
Beckman AU series BSBE \ Lactate to pyruvate (L-P)	mmol/L	2.61	3.6%	10	291
Mindray BS series BSBE \ Lactate to pyruvate (L-P)	mmol/L	2.44	6.1%	10	222
Roche cobas 4000/6000/8000 Roche \ Lactate to pyruvate (L-P)	mmol/L	2.4	12%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Lactate to pyruvate (L-P)	mmol/L	2.41	2.1%	<5	<100
(LDH)乳酸脱氢酶					
Mindray BS series Mindray \ Lactate to pyruvate (L-P)	U/L	25	3.4%	12	289
Beckman AU series BSBE \ Lactate to pyruvate (L-P)	U/L	20.7	15.8%	10	298
Roche cobas c 501/502/701/702 Roche \ Lactate to pyruvate (L-P)	U/L	21	10%	<5	<100
Siemens ADVIA Centaur XP/XPT JingYuan(Shanghai) \ Lactate to pyruvate (L-P)	U/L	24.7	5.1%	<5	<100
(PRO/ TP)(总) 蛋白					
Siemens ADVIA Chemistry XPT Dedicated \ Biuret-rate	g/L	0.372	8%	5	109
Roche cobas 4000/6000/8000 Roche \ Biuret-endpoint	g/L	0.285	12%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Biuret-endpoint	g/L	0.276	5.6%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。