

		批号：4005ACA			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Succeeder SA series Succeeder \ Magnetic Beads	Seconds	84.5	1.5%	8	175
Werfen ACL series Dedicated \ Light scattering	Seconds	93.8	4.3%	24	598
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	102	5%	23	585
Werfen ACL series Werfen \ Light scattering	Seconds	101	10.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	112	14.3%	<5	<100
Stago STA series Dedicated \ Magnetic Beads	Seconds	105	10.5%	<5	<100
(AT-III)抗凝血酶III					
Werfen ACL series Dedicated \ Nephelometry	%	26.4	9.5%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	%	26.9	9.7%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Succeeder SA series Succeeder \ Immunoturbidimetry (ITA)	mg/L	1.72	1.3%	8	176
Werfen ACL series Dedicated \ Latex agglutination (LA)	mg/L	1.14	3.2%	24	598
Succeeder SF series Succeeder \ Latex agglutination (LA)	mg/L	1.91	2.5%	12	292
Werfen ACL series Dedicated \ Nephelometry	mg/L	2.88	10.1%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	5.26	10.3%	<5	<100
Stago STA series Dedicated \ Latex agglutination (LA)	mg/L	4.38	9.8%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Werfen ACL series Dedicated \ Nephelometry	mg/L	6.67	11.8%	24	235
Sysmex CA 1500/6000/7000 SekiSui \ Nephelometry	mg/L	8.23	10.4%	<5	<100
Stago STA series Dedicated \ Nephelometry	mg/L	7.89	10.3%	<5	<100
(FIB)纤维蛋白原					
Succeeder SA series Succeeder \ Magnetic Beads	g/L	1.33	2.4%	8	175
Werfen ACL series Dedicated \ Light scattering	g/L	1.26	4.5%	24	597
Succeeder SF series Succeeder \ Magnetic Beads	g/L	1.32	10.1%	21	512
Stago STA series Dedicated \ Magnetic Beads	g/L	1.46	12.3%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	1.32	12.1%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	1.35	9.6%	<5	<100
(PT INR)PT国际标准化比值					
Succeeder SF series Succeeder \ Calculate	INR	3.83	2.3%	12	295
(PT)血浆凝血酶原时间					
Succeeder SA series Succeeder \ Magnetic Beads	Seconds	50	4.1%	8	176
Werfen ACL series Dedicated \ Light scattering	Seconds	49.2	4.5%	24	596
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	49.3	3.2%	23	569
Stago STA series Dedicated \ Light scattering	Seconds	43.6	9.9%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	42.5	10.4%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	42.3	10.2%	<5	<100
(TT)凝血酶时间					
Succeeder SA series Succeeder \ Magnetic Beads	Seconds	16.3	2.4%	8	176
Werfen ACL series Dedicated \ Light scattering	Seconds	22.8	4%	24	594
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	15.4	1.6%	23	564
Werfen ACL series Werfen \ Light scattering	Seconds	23.5	9.4%	<5	<100
Sysmex CA 1500/6000/7000 Dedicated \ Light scattering	Seconds	16.8	10.1%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。