

质控品名： 复合免疫分析非定值质控品

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
1 / 6								

		批号：3382991			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(25-OH-Vit D)25羟基维生素D					
Roche cobas E601/602	ng/mL	24.1	10.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	nmol/L	60.2	10.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
YHLO iFlash Analyzer	ng/mL	48.2	10%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
YHLO iFlash Analyzer	nmol/L	120	10%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	ng/mL	18.6	10.6%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	nmol/L	46.5	10.6%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
(AFP)甲胎蛋白					
maglumi system	ng/mL	8.1	10.6%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Beckman AU series	ng/mL	5.79	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	8.8	2.8%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	6.46	1.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	ng/mL	7.23	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
chemclin LiCA Series	ng/mL	9.7	1.5%	<5	<100
Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)					
Hotgen C2000	ng/mL	5.72	1%	<5	<100
Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)					
(ALD)醛固酮					
Autolumo A2000/A2000 PLUS	pg/mL	174	10.3%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(Anti-Tg)抗甲状腺球蛋白抗体					
Roche cobas E601/602	IU/mL	15.2	21.7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e 801	IU/mL	19.6	2.8%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
YHLO iFlash Analyzer	IU/mL	32.5	9.4%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
chemclin LiCA Series	IU/mL	17.4	0.5%	<5	<100
Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)					
(Anti-TPO)抗甲状腺过氧化物酶抗体					
chemclin LiCA Series	IU/mL	17.9	1.1%	<5	<100
Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)					
Roche cobas E601/602	IU/mL	13.8	10.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
YHLO iFlash Analyzer	IU/mL	11.2	2.7%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	IU/mL	16.4	8.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(CA 125)糖类抗原 125					
Roche cobas e 801	IU/mL	26	0.8%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	IU/mL	29.8	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	33.1	3.9%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Hotgen C2000	IU/mL	24.6	4.1%	<5	<100
Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)					
Beckman AU series	IU/mL	29.5	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
(CA 15-3)糖类抗原 15-3					
Roche cobas E601/602	IU/mL	25.2	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e 801	IU/mL	28.1	2.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	24.1	9.9%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Beckman AU series	IU/mL	15	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Hotgen C2000	IU/mL	23.6	3.5%	<5	<100
Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)					
(CA 19-9)糖类抗原 19-9					
Roche cobas e 801	IU/mL	28.8	0.8%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Beckman AU series	IU/mL	38.8	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	IU/mL	23.2	5.5%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	24.9	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
2 / 6								

批号：3382991					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(CEA)癌胚抗原					
Hotgen C2000	ng/mL	1.7	4.6%	<5	<100
Hotgen \ Chemiluminescence enzyme immunoassay (CLEIA)					
maglumi system	ng/mL	3.58	10.6%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Beckman AU series	ng/mL	2.22	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	3.76	10.1%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	ng/mL	2.72	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas e 801	ng/mL	2.79	1.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
chemclin LiCA Series	ng/mL	3.89	2.5%	<5	<100
Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)					
(CK)肌酸激酶					
Roche cobas 4000/6000/8000	U/L	126	10%	<5	<100
Roche \ Urease-UV-rate					
(CK-MB)肌酸激酶-MB (活性)					
Roche cobas 4000/6000/8000	U/L	6.96	9.9%	<5	<100
Roche \ Immunoturbidimetry (ITA)					
Roche cobas e411	ng/mL	4.48	9.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(Cortisol)皮质醇					
Roche cobas E601/602	ng/mL	97.1	2.4%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	ng/mL	129	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	108	1.4%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(C-P)C-肽					
Roche cobas E601/602	ng/mL	1.1	2.9%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Roche cobas E601/602	ng/mL	1.06	10.2%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	ng/mL	0.861	5.6%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
(DGXN)地高辛					
Roche cobas e411	ng/mL	1.04	10.5%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(DHEA-S)硫酸脱氢表雄酮					
Roche cobas E601/602	ng/mL	3651	10%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(E2)雌二醇					
Roche cobas E601/602	pmol/L	283	4.8%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
Autolumo A2000/A2000 PLUS	pmol/L	325	9.6%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	pmol/L	273	0.7%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
maglumi system	pmol/L	74.8	10.4%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	pmol/L	285	10.4%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(EPO)促红细胞生成素					
Siemens Immulite1000/2000	mIU/mL	5.73	10%	<5	<100
Siemens \ Chemiluminescence Immunoassay (CLIA)					
(FA)叶酸					
Roche cobas E601/602	ng/mL	3.81	10.1%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					
(Ferritin)铁蛋白					
Roche cobas E601/602	ng/mL	103	7.3%	<5	<100
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)					
YHLO iFlash Analyzer	ng/mL	69.4	3%	<5	<100
YHL \ Chemiluminescence Immunoassay (CLIA)					
Beckman AU series	ng/mL	49.2	10%	<5	<100
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)					
maglumi system	ng/mL	68.6	10.5%	<5	<100
Snibe \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas 4000/6000/8000	ng/mL	116	10.3%	<5	<100
Roche \ Immunoturbidimetry (ITA)					
(F-PSA)游离前列腺特异性抗原					
Autolumo A2000/A2000 PLUS	ng/mL	1.99	3.6%	<5	<100
Autobio \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas e 801	ng/mL	1.66	0.3%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
3 / 6								

批号：3382991					
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
chemclin LICA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	2.94	0.4%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.11	10%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.79	10%	<5	<100
(FSH)促卵泡刺激素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.26	5.2%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.46	0.3%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.25	9.7%	<5	<100
(FT3)游离T3					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.37	3.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.96	10%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	5.75	2.1%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.12	10%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.75	10%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	5.43	0.3%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.33	1%	<5	<100
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.88	1.7%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.01	10%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	5.09	2.7%	<5	<100
(FT4)游离T4					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	17.4	4.5%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	175	9.6%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	pmol/L	172	10.5%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12	10%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	pmol/L	21.3	2%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	pmol/L	14.7	0.6%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	15.9	1.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	16.6	9.7%	<5	<100
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.7	2.4%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	pmol/L	17.1	6%	<5	<100
(hCG)人绒毛膜促性腺激素					
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.39	10.9%	<5	<100
(HCY)同型半胱氨酸					
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	13.6	9.9%	<5	<100
(HGH)人生长激素					
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	3.71	2.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	1.86	10%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Image 800 Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	148	9.6%	<5	<100
(IGF1)胰岛素生长因子1					
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	228	11.4%	<5	<100
(INS)胰岛素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.93	3.1%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	0.909	10.9%	<5	<100

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
4 / 6								

		批号：3382991			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.61	3.6%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	8.63	9.9%	<5	<100
(LH)促黄体生成素					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.22	4%	<5	<100
Roche cobas e 801					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.22	0.8%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	IU/L	5.51	4.7%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.76	9.3%	<5	<100
(Myo / Mb)肌红蛋白					
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	28.9	9.8%	<5	<100
(NT-proBNP)N-末端脑钠肽前体					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	189	4.2%	<5	<100
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	22.2	4.2%	<5	<100
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	178	9.7%	<5	<100
Roche cobas e411					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	20.9	9.7%	<5	<100
(PRL)催/ 泌乳素					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.67	5.1%	<5	<100
Roche cobas e 801					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.49	1.2%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.43	3.8%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.62	10.1%	<5	<100
(PROG)孕酮					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.18	6.4%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.14	10.1%	<5	<100
Roche cobas e 801					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.19	1%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.19	6.4%	<5	<100
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.33	10.4%	<5	<100
(PSA)前列腺特异性抗原					
Roche cobas e 801					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.41	2.3%	<5	<100
chemclin LiCA Series					
Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	5.34	0.8%	<5	<100
Beckman AU series					
Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5	10%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.92	10%	<5	<100
(PTH)甲状旁腺素					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.74	8%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1.85	10%	<5	<100
YHLO iFlash Analyzer					
YHL \ Chemiluminescence Immunoassay (CLIA)	pmol/L	3.34	3.7%	<5	<100
maglumi system					
Snibe \ Chemiluminescence Immunoassay (CLIA)	pmol/L	0.496	10.2%	<5	<100
(TES)睾酮					
Roche cobas E601/602					
Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.18	3.2%	<5	<100
Autolumo A2000/A2000 PLUS					
Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.43	2.4%	<5	<100
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.53	10.2%	<5	<100
Roche cobas e 801					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	6.64	1.2%	<5	<100
(Tg)甲状腺球蛋白					
Roche cobas E601/602					
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.4	12.6%	<5	<100

质控品名： **复合免疫分析非定值质控品**

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
5 / 6								

批号：3382991					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(TnI)肌钙蛋白I					
Beckman UniCel DxI 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.07	8.6%	<5	<100
(TnT)肌钙蛋白T					
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.025	7.2%	<5	<100
(TSH)促甲状腺素					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.36	4%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.74	10.6%	<5	<100
YHLO iFlash Analyzer YHL \ Chemiluminescence Immunoassay (CLIA)	mIU/L	4.05	1.8%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.58	1.3%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.65	10.2%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.85	8.9%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	3.42	1%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	3.29	1.4%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.89	10%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	mIU/L	3.57	1.4%	<5	<100
(TT3)总三碘甲状原氨酸					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.2	4.4%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.26	10.2%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.02	10%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.28	9.9%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.25	0.2%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.01	0.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.42	9.7%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.28	0.7%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.9	8.1%	<5	<100
(TT4)总T4					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	56.7	3.4%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	ng/mL	70.4	9.9%	<5	<100
chemclin LiCA Series Chemclin \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	57.2	0.3%	<5	<100
Roche cobas e 801 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	60.5	0.9%	<5	<100
Beckman AU series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	76.9	10%	<5	<100
Sysmex HISCL-5000 sysmex \ Chemiluminescence Immunoassay (CLIA)	ng/mL	67.2	9.9%	<5	<100
Mindray BC-6000 series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	70.4	0.9%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	72.1	9.9%	<5	<100
Roche cobas E601/602 Roche \ Chemiluminescence Immunoassay (CLIA)	ng/mL	62.8	3.4%	<5	<100
(VB12)维生素B12					
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	700	3.8%	<5	<100
maglumi system Snibe \ Chemiluminescence Immunoassay (CLIA)	pg/mL	196	10.1%	<5	<100
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	480	9.9%	<5	<100
(β2M)β-2-微球蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	840	10.6%	<5	<100

REF	540-1	6×3mL	EXP	2027-05-18	LOT	3382991	更新日期：	2024-05
6 / 6								

		批号：3382991			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(β-hCG)人绒毛膜促性腺激素β亚基					
Roche cobas E601/602	IU/mL	0.004	5.5%	<5	<100
Roche \ Chemiluminescence Immunoassay (CLIA)					
Roche cobas E601/602	IU/mL	0.004	9.9%	<5	<100
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)					

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。