

REF	671-3	6×1mL	EXP	2024-12-20	LOT	4005B2E	更新日期：	2024-04
1 / 2								

		批号：4005B2E			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Werfen ACL series Dedicated \ Light scattering	Seconds	89.7	6.6%	5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	85	7.7%	<5	169
Sysmex CS 5100 sysmex \ Light scattering	Seconds	111	37.1%	<5	125
Werfen ACL series ACCURDX \ Magnetic Beads	Seconds	61	1.2%	<5	<100
Sysmex CS-2500 sysmex \ Light scattering	Seconds	147	10%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	76.1	0.5%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Light scattering	Seconds	106	3%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	75.8	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	125	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	84.8	12%	<5	<100
(AT-III)抗凝血酶III					
Werfen ACL series Werfen \ Nephelometry	%	27.1	8.6%	<5	110
Sysmex CS 5100 sysmex \ Nephelometry	%	27.8	10.3%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	%	30.8	2.8%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	%	29.3	12%	<5	<100
Stago STA series Stago \ Magnetic Beads	%	31.3	12%	<5	<100
Werfen ACL series ACCURDX \ Nephelometry	%	26.3	0.8%	<5	<100
Sysmex CS-2500 sysmex \ Nephelometry	%	25.3	16.1%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Nephelometry	%	19.7	4.3%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Werfen ACL series Dedicated \ Latex agglutination (LA)	mg/L	1.45	9%	5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	2.44	48.7%	<5	118
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	5.11	8%	<5	<100
Werfen ACL series ACCURDX \ Nephelometry	mg/L	5.12	5%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Nephelometry	mg/L	4.96	1.8%	<5	<100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	5.39	12%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	4.6	7.1%	<5	<100
Succeeder SF series Succeeder \ Nephelometry	mg/L	1.32	0.7%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	4.02	12%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Werfen ACL series Dedicated \ Nephelometry	mg/L	7.1	10.3%	<5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	2.27	169.1%	<5	100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	9.9	12%	<5	<100
Werfen ACL series ACCURDX \ Nephelometry	mg/L	13.8	1.3%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	9.52	5%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Nephelometry	mg/L	7.98	3.6%	<5	<100
Succeeder SF series Succeeder \ Nephelometry	mg/L	3.82	0.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	9.73	12%	<5	<100
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	7.33	35.6%	<5	107
(FIB)纤维蛋白原					
Werfen ACL series Dedicated \ Light scattering	g/L	1.22	5.2%	5	<100
Werfen ACL series Werfen \ Magnetic Beads	g/L	1.25	8.4%	<5	149
Sysmex CS 5100 sysmex \ Light scattering	g/L	1.24	10.9%	<5	145
Werfen ACL series ACCURDX \ Magnetic Beads	g/L	1.25	1.2%	<5	<100
Sysmex CS-2500 sysmex \ Light scattering	g/L	1.23	4.4%	<5	<100

REF	671-3	6×1mL	EXP	2024-12-20	LOT	4005B2E	更新日期：	2024-04
2 / 2								

		批号：4005B2E			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Sekisui Coapresta 3000 SekiSui \ Light scattering	g/L	0.015	2.5%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	1.43	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	1.28	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	1.58	12%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Light scattering	g/L	1.41	2.6%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	g/L	2.5	1.6%	<5	<100
(PT INR)PT国际标准化比值					
Werfen ACL series Werfen \ Calculate	INR	4.31	9.5%	<5	<100
Sysmex CS 5100 sysmex \ Calculate	INR	4.15	4%	<5	112
Werfen ACL series ACCURDX \ Calculate	INR	2.86	22.6%	<5	<100
Stago STA series Stago \ Calculate	INR	5.6	12%	<5	<100
Succeeder SF series Succeeder \ Calculate	INR	2.55	0.4%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Calculate	INR	4.61	7.3%	<5	<100
Sysmex CS 5100 Siemens \ Calculate	INR	4.39	12%	<5	<100
Sysmex CS-2500 sysmex \ Calculate	INR	3.97	5.8%	<5	<100
(PT)血浆凝血酶原时间					
Werfen ACL series Dedicated \ Light scattering	Seconds	49	3.9%	5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	47.7	6.2%	<5	135
Sysmex CS 5100 sysmex \ Light scattering	Seconds	45.8	7%	<5	143
Sysmex CS-2500 sysmex \ Light scattering	Seconds	45.4	5.7%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	49.2	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	44.4	12%	<5	<100
Werfen ACL series ACCURDX \ Magnetic Beads	Seconds	35	10.7%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	47.5	1.1%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	46.6	12%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Light scattering	Seconds	42.6	0.6%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	33.1	0.4%	<5	<100
(TT)凝血酶时间					
Werfen ACL series Dedicated \ Light scattering	Seconds	23.6	4.6%	5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	23.2	7.2%	<5	134
Sysmex CS 5100 sysmex \ Light scattering	Seconds	14.1	9.6%	<5	144
Sysmex CS-2500 sysmex \ Light scattering	Seconds	13.2	2.7%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	15.5	12%	<5	<100
Sekisui Coapresta 3000 SUNBIO \ Light scattering	Seconds	27.7	2.5%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	14.9	3.4%	<5	<100
Werfen ACL series ACCURDX \ Magnetic Beads	Seconds	17.6	2%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	15.7	1.1%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。