

REF	671-3	6×1mL	EXP	2025-09-29	LOT	4005B90	更新日期：	2024-04
1 / 3								

		批号：4005B90			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Succeeder SF series Succeeder \ clauss	Seconds	105	2.3%	5	117
Werfen ACL series Dedicated \ Light scattering	Seconds	38	2.6%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	100	2.7%	<5	129
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	122	9.3%	<5	<100
Roche cobas T711 Roche \ clauss	Seconds	147	9.9%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Light scattering	Seconds	94	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Dade Actin(APTT) \ Light scattering	Seconds	143	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Pathrombin SL(APTT) \ Light scattering	Seconds	148	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	105	3.8%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ Light scattering	Seconds	101	10%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ clauss	Seconds	99	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	85	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	30.4	1.9%	<5	<100
(AT-III)抗凝血酶III					
Werfen ACL series Werfen \ Nephelometry	%	23.5	18.7%	<5	153
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	%	19.8	33.8%	<5	<100
Roche cobas T711 Roche \ Nephelometry	%	26.6	9.9%	<5	<100
Werfen ACL series Werfen IL Antithrombin III Chromogenic (AT□) \ Nephelometry	%	136	10%	<5	<100
Werfen ACL series Werfen Hemosil Liquid Antithrombin (AT□) \ Nephelometry	%	25	10%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	%	28.3	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	%	402	1.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Berichrom Antithrombin III (AT□) \ Nephelometry	%	27	10%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Succeeder SF series Succeeder \ Immunoturbidimetry (ITA)	mg/L	1.7	0.6%	5	116
Werfen ACL series Dedicated \ Nephelometry	mg/L	1.32	2.1%	<5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	1.62	56.5%	<5	169
Werfen ACL series Werfen Hemosil D-Dimer HS \ Nephelometry	mg/L	2.56	54%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	5.17	8.9%	<5	<100
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	4.92	6.3%	<5	<100
Werfen ACL series Wenfen Hemosil D-Dimer \ Nephelometry	mg/L	1.07	10%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	mg/L	3.58	1.5%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Roche cobas T711 Roche \ Nephelometry	mg/L	2	9.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Innovance Antithrombin (AT□) \ Nephelometry	mg/L	3.67	10%	<5	<100
Werfen ACL series MEDCONN \ OTHER	mg/L	3.84	5.4%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	mg/L	4.89	10%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Werfen ACL series Dedicated \ Nephelometry	mg/L	2.29	15.6%	<5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	6.01	11.8%	<5	111
Roche cobas T711 Roche \ Nephelometry	mg/L	15.7	9.9%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	11.8	14.6%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	4.82	2.2%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Nephelometry	mg/L	8.36	0.2%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Nephelometry	mg/L	13.2	10%	<5	<100
(FIB)纤维蛋白原					
Succeeder SF series Succeeder \ clauss	g/L	1.33	1%	5	116

REF	671-3	6×1mL	EXP	2025-09-29	LOT	4005B90	更新日期：	2024-04
2 / 3								

		批号：4005B90			
项目\仪器\试剂方法	单位	均值 \ CV		Lab \ 点数	
Werfen ACL series Dedicated \ Light scattering	g/L	3.96	4.9%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	g/L	0.988	52.9%	<5	154
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	g/L	1.29	25.7%	<5	<100
Roche cobas T711 Roche \ clauss	g/L	1.37	9.9%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	1.36	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Multifibren U (FIB) \ Light scattering	g/L	1.75	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	g/L	1.1	9.2%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	g/L	5.37	0.7%	<5	<100
Werfen ACL series Werfen HemosIL QFA Thrombin (FIB) \ Light scattering	g/L	1.22	10%	<5	<100
Werfen ACL series Werfen HemosIL Fib-C \ Light scattering	g/L	1.11	10%	<5	<100
Werfen ACL series Werfen HemosIL QFA Thrombin (FIB) \ clauss	g/L	1.82	10%	<5	<100
(PT INR)PT国际标准化比值					
Werfen ACL series Dedicated \ Calculate	INR	1.15	2.6%	<5	<100
Roche cobas T711 Roche \ Calculate	INR	3.34	9.9%	<5	<100
Werfen ACL series Werfen \ Calculate	INR	4.31	9.9%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Calculate	INR	4.65	10%	<5	<100
(PT)血浆凝血酶原时间					
Succeeder SF series Succeeder \ clauss	Seconds	55.6	1.4%	5	116
Werfen ACL series Dedicated \ Light scattering	Seconds	12.2	2.6%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	47.9	9%	<5	136
Roche cobas T711 Roche \ clauss	Seconds	32.1	9.9%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	47.4	7.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens PT Innovin \ Light scattering	Seconds	42.8	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	21.3	0.9%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Light scattering	Seconds	52.1	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Light scattering	INR	4.43	9.9%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ Light scattering	Seconds	47.9	10%	<5	<100
Werfen ACL series Werfen Hemosil PT-Fib HS Plus \ Light scattering	INR	5.1	10%	<5	<100
Werfen ACL series Werfen Hemosil PT-Fib HS Plus \ Light scattering	Seconds	62.7	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	42.5	3%	<5	<100
Sysmex CA 1500/6000/7000 Siemens PT Innovin \ Light scattering	INR	4.04	9.9%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Thromborel S (PT) \ Light scattering	Seconds	49.5	10%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Thromborel S (PT) \ Light scattering	INR	4.55	10.1%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	51.9	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ clauss	Seconds	40.5	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	INR	5.8	10%	<5	<100
(TT)凝血酶时间					
Succeeder SF series Succeeder \ clauss	Seconds	16.4	0.6%	5	116
Werfen ACL series Dedicated \ Light scattering	Seconds	15	3%	<5	<100
Roche cobas T711 Roche \ clauss	Seconds	30	9.9%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	25.8	4.5%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	17.9	11.7%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	Seconds	15.6	1.6%	<5	<100
Sysmex CA 1500/6000/7000 sysmex \ Light scattering	Seconds	15.6	0.6%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ clauss	Seconds	31	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	46.8	0.7%	<5	<100

报告类别： 汇总实验室数据报告



质控品名： 凝血七项非定值质控品

REF	671-3	6×1mL	EXP	2025-09-29	LOT	4005B90	更新日期：	2024-04
								3 / 3

		批号：4005B90	
项目\仪器\试剂\方法	单位	均值 \ CV	Lab \ 点数

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。