

质控品名： 特定蛋白液体非定值质控品

REF	521-2	6×3mL	EXP	2024-10-17	LOT	31B51FE	更新日期：	2024-04
1 / 2								

		批号：31B51FE			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(Alb)白蛋白					
Roche cobas 4000/6000/8000 Roche \ Bromocresol green	g/L	39.7	12%	<5	<100
(ASO)抗链球菌溶血素 O 抗体					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	251	1.9%	6	129
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	257	3.8%	5	<100
Siemens ADVIA 1800/2400 chongqiong zybio \ Immunoturbidimetry (ITA)	IU/mL	254	3.6%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	IU/mL	217	2.4%	20	588
Beckman Image 800 OTHER \ Immunonephelometry (INA)	IU/mL	145	3.9%	13	<100
Mindray BS series Mindray \ Immunoturbidimetry (ITA)	IU/mL	252	3.5%	<5	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	142	12%	<5	<100
(AT-III)抗凝血酶III					
Stago STA series Dedicated \ Nephelometry	%	18	12%	<5	<100
(C3)补体 C3					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	164	2.9%	6	129
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	161	3.1%	21	411
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	152	4%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	149	12%	<5	<100
(C4)补体 C4					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	33.5	3.9%	6	129
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	38.6	4.3%	17	<100
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	40.7	9.4%	14	317
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	39	12%	<5	<100
(Cer)铜蓝蛋白					
Beckman Image 800 OTHER \ Immunonephelometry (INA)	g/L	0.274	5.4%	10	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	0.273	12%	<5	<100
(CRP)C反应蛋白					
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	3.55	9.4%	9	132
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/dL	3.46	3.4%	20	462
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	2.91	6.1%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	2.64	12%	<5	<100
(Ferritin)铁蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	ng/mL	184	12%	<5	<100
(IgA)免疫球蛋白 IgA					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	321	3%	6	129
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	267	2.4%	22	440
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	242	3.4%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	242	12%	<5	<100
(IgE)免疫球蛋白 IgE					
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	51.4	12%	<5	<100
(IgG)免疫球蛋白 IgG					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	1536	1%	6	129
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1443	5.5%	21	404
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	1356	3.5%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	1390	12%	<5	<100
(IgM)免疫球蛋白 IgM					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	mg/dL	173	1.4%	6	129
Beckman AU series Beckman Coulter \ Immunonephelometry (INA)	mg/dL	177	3.7%	21	449

REF	521-2	6×3mL	EXP	2024-10-17	LOT	31B51FE	更新日期：	2024-04
2 / 2								

		批号：31B51FE			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Beckman Image 800 OTHER \ Immunonephelometry (INA)	mg/dL	166	3.2%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	mg/dL	166	12%	<5	<100
(Kappa Light Chain)κ轻链抗体					
Beckman Image 800 OTHER \ Immunonephelometry (INA)	g/L	11.4	4.5%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	11.9	12%	<5	<100
(Lambda Light Chain)λ 轻链抗体					
Beckman Image 800 OTHER \ Immunonephelometry (INA)	g/L	5.8	5.1%	17	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	g/L	5.74	12%	<5	<100
(PAB)前白蛋白					
Roche cobas 4000/6000/8000 DiaSys \ Immunonephelometry (INA)	g/L	0.309	12%	<5	<100
(RBP)视黄醇结合蛋白					
Roche cobas 4000/6000/8000 WanTai \ Immunonephelometry (INA)	mg/dL	5.75	12%	<5	<100
(RF)类风湿因子					
Siemens BN ProSpec Siemens \ Immunonephelometry (INA)	IU/mL	59.7	3.3%	6	128
Beckman AU series Beckman Coulter \ Immunoturbidimetry (ITA)	IU/mL	77.2	4%	5	<100
Siemens ADVIA 1800/2400 DiaSys \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	77.1	8.7%	<5	<100
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	IU/mL	80.1	2.9%	20	482
Beckman Image 800 OTHER \ Immunonephelometry (INA)	IU/mL	74.1	3.7%	13	<100
Beckman Image 800 Beckman Coulter \ Immunonephelometry (INA)	IU/mL	77.6	12%	<5	<100
(TRF)转铁蛋白					
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/dL	252	12%	<5	<100
(β2M)β-2-微球蛋白					
Roche Hitachi/Hitachi Systems MedicalSystem \ Immunoturbidimetry (ITA)	ng/mL	6429	7.5%	11	231
Roche cobas 4000/6000/8000 JS (Shanghai) \ Immunonephelometry (INA)	ng/mL	4500	12%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。