

REF	671-2	6×1mL	EXP	2025-09-03	LOT	400347F	更新日期：	2024-04
1 / 3								

		批号：400347F			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sekisui Coapresta 3000 SekiSui \ OTHER	Seconds	58	4.2%	9	241
Mindray EXC810 Mindray \ OTHER	Seconds	89.6	3.8%	8	198
Sysmex CS-2400 sysmex \ OTHER	Seconds	75.6	4.3%	5	175
OTHER OTHER \ Magnetic Beads	Seconds	66.2	3%	5	128
Sekisui Coapresta 3000 SUNBIO \ OTHER	Seconds	79.9	4.7%	<5	111
Sekisui Coapresta 3000 OTHER \ OTHER	Seconds	74.7	4.2%	<5	<100
Werfen ACL series Werfen \ clauss	Seconds	56.4	22.2%	<5	<100
Werfen ACL series Werfen \ OTHER	Seconds	58.1	3.1%	17	362
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	60.8	3.5%	<5	<100
Mindray C Series CHANGDAO \ Magnetic Beads	Seconds	95.2	0.6%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	64.5	2.3%	<5	<100
Sysmex CN Series sysmex \ clauss	Seconds	69.9	1.1%	<5	<100
Dirui BCA Series DIRUI \ clauss	Seconds	91.4	10%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	115	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	25.2	2%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	60.9	1.8%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ clauss	Seconds	73.2	10%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Light scattering	Seconds	64	10%	<5	<100
Sekisui Coapresta 3000 OTHER \ clauss	Seconds	82.8	4.6%	<5	<100
(AT-III)抗凝血酶III					
Sekisui Coapresta 3000 SekiSui \ OTHER	%	38.3	8.5%	8	235
OTHER OTHER \ OTHER	%	51.3	8.5%	5	127
Sysmex CS-2400 sysmex \ OTHER	%	44.3	5.1%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	%	191	85.6%	<5	<100
Werfen ACL series Werfen \ OTHER	%	41.5	8.4%	17	367
Stago STA series Stago \ Nephelometry	%	136	2.6%	<5	<100
Werfen ACL series Werfen IL Antithrombin III Chromogenic (AT□) \ Nephelometry	%	50.2	10%	<5	<100
Werfen ACL series Werfen Hemosil Liquid Antithrombin (AT□) \ Nephelometry	%	40	10%	<5	<100
Werfen ACL series Werfen \ Nephelometry	%	44.5	1.1%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Sekisui Coapresta 3000 SekiSui \ OTHER	mg/L	3.47	7.8%	8	233
Mindray EXC810 Mindray \ OTHER	mg/L	1.61	11.3%	8	197
OTHER OTHER \ Immunoturbidimetry (ITA)	mg/L	2.44	9%	5	139
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	3.45	6.5%	<5	<100
Sysmex CS-2400 Siemens \ Nephelometry	mg/L	3.04	8.1%	<5	<100
Werfen ACL series Werfen \ OTHER	mg/L	0.81	9.3%	19	371
Werfen ACL series Werfen \ Nephelometry	mg/L	0.854	2.8%	<5	<100
Stago STA series Stago \ Nephelometry	mg/L	0.86	3.5%	<5	<100
Sysmex CN Series sysmex \ Nephelometry	mg/L	4.41	4.5%	<5	<100
Werfen ACL series Werfen Hemosil D-Dimer HS \ Nephelometry	mg/L	3.53	10%	<5	<100
Dirui BCA Series DIRUI \ Nephelometry	mg/L	1.3	10%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	3.7	10%	<5	<100
Sekisui Coapresta 3000 OTHER \ Dot-immuogold filtration assay (DIGFA)	mg/L	3.7	7%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Sysmex CS-2400 Siemens \ Latex agglutination (LA)	mg/L	2.97	8.1%	<5	149
Mindray C Series CHANGDAO \ Nephelometry	mg/L	2.93	5.4%	<5	<100

REF	671-2	6×1mL	EXP	2025-09-03	LOT	400347F	更新日期：	2024-04
2 / 3								

		批号：400347F			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	3	10%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Sekisui Coapresta 3000 SekiSui \ OTHER	mg/L	5.54	4.3%	8	226
OTHER OTHER \ Nephelometry	mg/L	12.9	15.9%	5	138
Werfen ACL series Werfen \ OTHER	mg/L	5	6.4%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Nephelometry	mg/L	4.77	17.1%	<5	<100
Werfen ACL series Werfen \ Nephelometry	mg/L	4.01	2.2%	<5	<100
Werfen ACL series Werfen Hemosil APTT-SP \ Nephelometry	mg/L	10.9	10%	<5	<100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	5.61	10.7%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	20	10%	<5	<100
Stago STA series Stago \ Nephelometry	mg/L	3.49	0.6%	<5	<100
Sekisui Coapresta 3000 OTHER \ Nephelometry	mg/L	6.3	3.8%	<5	<100
(FIB)纤维蛋白原					
Sekisui Coapresta 3000 SekiSui \ OTHER	g/L	2.27	5.2%	9	244
Mindray EXC810 Mindray \ OTHER	g/L	1.98	5.7%	8	195
Sysmex CS-2400 sysmex \ OTHER	g/L	1.8	5.7%	5	138
OTHER OTHER \ Magnetic Beads	g/L	1.8	6.3%	5	126
Sekisui Coapresta 3000 SUNBIO \ OTHER	g/L	1.76	8.4%	<5	108
Sekisui Coapresta 3000 OTHER \ OTHER	g/L	2.02	5.9%	<5	<100
Werfen ACL series Werfen \ clauss	g/L	1.68	10.3%	<5	<100
Werfen ACL series Werfen \ OTHER	g/L	2	6.2%	17	383
Sysmex CN Series sysmex \ clauss	g/L	3.53	1.5%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	g/L	4	0.2%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	g/L	2	0.8%	<5	<100
Werfen ACL series Werfen Hemosil QFA Thrombin (FIB) \ Light scattering	g/L	1.94	10%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	g/L	1.3	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	g/L	2.24	10%	<5	<100
Werfen ACL series Werfen Hemosil QFA Thrombin (FIB) \ clauss	g/L	2.17	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	4.18	2%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	g/L	0.023	5.8%	<5	<100
Mindray C Series CHANGDAO \ Magnetic Beads	g/L	2.03	1.8%	<5	<100
Sekisui Coapresta 3000 OTHER \ clauss	g/L	1.89	6.1%	<5	<100
(PT INR)PT国际标准化比值					
Mindray EXC810 Mindray \ Calculate	INR	1.98	4.2%	8	199
Sysmex CS-2400 sysmex \ Calculate	INR	2.65	4.6%	5	174
Sekisui Coapresta 3000 SUNBIO \ Calculate	INR	2.39	3.7%	<5	108
Sekisui Coapresta 3000 OTHER \ Calculate	INR	2.36	8.4%	<5	<100
Werfen ACL series Werfen \ Calculate	INR	2.83	5.3%	16	340
Sekisui Coapresta 3000 SekiSui \ Calculate	INR	2.85	4.2%	10	304
Dirui BCA Series DIRUI \ Calculate	INR	2	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Calculate	INR	2.62	9.9%	<5	<100
Mindray C Series CHANGDAO \ Calculate	INR	1.97	0.6%	<5	<100
(PT)血浆凝血酶原时间					
Sekisui Coapresta 3000 SekiSui \ OTHER	Seconds	35.1	4%	9	252
Mindray EXC810 Mindray \ OTHER	Seconds	24.4	3.3%	8	196
Sysmex CS-2400 sysmex \ OTHER	Seconds	28.3	4.6%	5	172
OTHER OTHER \ Magnetic Beads	Seconds	32.2	9.8%	5	134

REF	671-2	6×1mL	EXP	2025-09-03	LOT	400347F	更新日期：	2024-04
								3 / 3

		批号：400347F			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Sekisui Coapresta 3000 SUNBIO \ OTHER	Seconds	27.7	3.2%	<5	109
Sekisui Coapresta 3000 OTHER \ OTHER	Seconds	26.3	3.4%	<5	<100
Werfen ACL series Werfen \ OTHER	Seconds	31	3.6%	17	359
Werfen ACL series Werfen Hemosil PT Recombiplastin 2G \ Light scattering	Seconds	27.9	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	29.1	2.4%	<5	<100
Sysmex CN Series sysmex \ clauss	Seconds	29.6	2.7%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	15.1	1.6%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	19.7	1%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	25.2	1%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	45.6	10%	<5	<100
Werfen ACL series Werfen Hemosil PT Recombiplastin \ clauss	Seconds	29.4	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	Seconds	23.3	10%	<5	<100
Mindray C Series CHANGDAO \ Magnetic Beads	Seconds	24	0.7%	<5	<100
Werfen ACL series Werfen \ clauss	Seconds	31.6	3.1%	<5	<100
Sekisui Coapresta 3000 OTHER \ clauss	Seconds	28.4	6.3%	<5	<100
(TT)凝血酶时间					
Sekisui Coapresta 3000 SekiSui \ OTHER	Seconds	20.6	2.6%	9	236
Mindray EXC810 Mindray \ OTHER	Seconds	17.7	2.8%	8	196
Sysmex CS-2400 sysmex \ OTHER	Seconds	16.9	4.8%	5	172
OTHER OTHER \ Magnetic Beads	Seconds	17.6	3.2%	5	115
Sekisui Coapresta 3000 SUNBIO \ OTHER	Seconds	18.5	4.4%	<5	110
Sekisui Coapresta 3000 OTHER \ OTHER	Seconds	18.4	7.5%	<5	<100
Werfen ACL series Werfen \ OTHER	Seconds	19.2	5.1%	17	367
Sysmex CN Series sysmex \ clauss	Seconds	18.8	1.2%	<5	<100
Werfen ACL series Wenfen Hemosil D-Dimer \ Light scattering	Seconds	18.1	10.5%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	22.4	0%	<5	<100
Werfen ACL series Werfen \ Magnetic Beads	Seconds	18.7	1.1%	<5	<100
Werfen ACL series Werfen Hemosil SynthAsil APTT \ clauss	Seconds	14.8	10%	<5	<100
Dirui BCA Series DIRUI \ clauss	Seconds	18.4	10%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Magnetic Beads	Seconds	32.6	3.6%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	14.6	10%	<5	<100
Mindray C Series CHANGDAO \ Magnetic Beads	Seconds	18.2	0.4%	<5	<100
Werfen ACL series Werfen \ clauss	Seconds	21.2	4.4%	<5	<100
Sekisui Coapresta 3000 SekiSui \ Light scattering	Seconds	17.3	3.2%	<5	<100
Sekisui Coapresta 3000 OTHER \ clauss	Seconds	12.2	3.9%	<5	<100

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。