

REF	671-2	6×1mL	EXP	2025-04-26	LOT	4003422	更新日期：	2024-04
1 / 2								

批号：4003422					
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
(APTT)活化部分凝血活酶时间					
Sysmex CS 5100 sysmex \ Light scattering	Seconds	150	5%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Dade Actin FS(APTT) \ Light scattering	Seconds	67.8	12%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	74.6	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens Pathrombin SL(APTT) \ Light scattering	Seconds	90.3	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	65.1	12%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	91.1	10%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	63.9	10%	<5	<100
(AT-III)抗凝血酶III					
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	%	46.5	12%	<5	<100
Stago STA series Stago \ Magnetic Beads	%	37.4	12%	<5	<100
Werfen ACL series Werfen \ Nephelometry	%	46.5	12%	<5	<100
(D-D(DDU)) D-二聚体(DDU)					
Werfen ACL series Werfen \ Nephelometry	mg/L	2.07	12%	<5	<100
Sysmex CS 5100 Siemens \ Nephelometry	mg/L	3.81	10%	<5	<100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	3.3	12%	<5	<100
Succeeder SF series Succeeder \ Nephelometry	mg/L	1.86	10%	<5	<100
Sysmex CS 5100 sysmex \ Nephelometry	mg/L	3.05	1.6%	<5	<100
(D-D(FEU)) D-二聚体(FEU)					
Stago STA series Stago \ Nephelometry	mg/L	1.76	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Nephelometry	mg/L	3.42	12%	<5	<100
(FDP)纤维蛋白(原)降解产物					
Succeeder SF series Succeeder \ Nephelometry	mg/L	4.26	10%	<5	<100
Sysmex CS 5100 SekiSui \ Nephelometry	mg/L	5.56	12%	<5	<100
(FIB)纤维蛋白原					
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	g/L	2	12%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	g/L	1.02	10%	<5	<100
Werfen ACL series Werfen \ Light scattering	g/L	2.18	12%	<5	<100
Stago STA series Stago \ Magnetic Beads	g/L	2.02	12%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	g/L	1.41	10%	<5	<100
Sysmex CS 5100 sysmex \ Light scattering	g/L	1.8	11.7%	<5	<100
(PT INR)PT国际标准化比值					
Sysmex CS 5100 Siemens \ Calculate	INR	2.39	26.6%	<5	<100
Stago STA series Stago \ Calculate	INR	2.88	12%	<5	<100
Sysmex CS 5100 sysmex \ Calculate	INR	2.88	3.3%	<5	<100
Werfen ACL series Werfen \ Calculate	INR	2.91	12%	<5	<100
Succeeder SF series Succeeder \ Calculate	INR	4.29	10%	<5	<100
(PT)血浆凝血酶原时间					
Sysmex CS 5100 sysmex \ Light scattering	Seconds	31.4	3%	<5	<100
Sysmex CS 5100 Siemens \ Light scattering	Seconds	17.2	10%	<5	<100
Stago STA series Stago \ Magnetic Beads	Seconds	29.9	12%	<5	<100
Werfen ACL series Werfen \ Light scattering	Seconds	31.7	12%	<5	<100
Sysmex CA 1500/6000/7000 Siemens \ Light scattering	Seconds	32.4	12%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	57	10%	<5	<100
(TT)凝血酶时间					
Sysmex CS 5100 sysmex \ Light scattering	Seconds	17.4	10.8%	<5	<100
Succeeder SF series Succeeder \ Magnetic Beads	Seconds	16.1	10%	<5	<100

REF	671-2	6×1mL	EXP	2025-04-26	LOT	4003422	更新日期：	2024-04
2 / 2								

		批号：4003422			
项目\仪器\试剂\方法	单位	均值 \ CV		Lab \ 点数	
Sysmex CS 5100	Seconds	20.3	10%	<5	<100
Siemens \ Light scattering					
Stago STA series	Seconds	16.2	12%	<5	<100
Stago \ Magnetic Beads					

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。