

|                |       |       |                |            |                |         |       |         |
|----------------|-------|-------|----------------|------------|----------------|---------|-------|---------|
| <div>REF</div> | 673-3 | 6×1mL | <div>EXP</div> | 2026-03-08 | <div>LOT</div> | 40368D5 | 更新日期： | 2024-05 |
| 1 / 1          |       |       |                |            |                |         |       |         |

| 批号：40368D5                                                                          |         |         |       |          |      |
|-------------------------------------------------------------------------------------|---------|---------|-------|----------|------|
| 项目\仪器\试剂\方法                                                                         | 单位      | 均值 \ CV |       | Lab \ 点数 |      |
| (APTT)活化部分凝血活酶时间                                                                    |         |         |       |          |      |
| Werfen ACL series<br>Werfen Hemosil APTT-SP \ Light scattering                      | Seconds | 78      | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil SynthAsil APTT \ Light scattering               | Seconds | 84      | 10%   | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens Dade Actin(APTT) \ Light scattering             | Seconds | 129     | 10%   | <5       | <100 |
| Stago STA series<br>Stago \ Magnetic Beads                                          | Seconds | 71      | 10%   | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens Pathrombin SL(APTT) \ Light scattering          | Seconds | 124     | 10%   | <5       | <100 |
| (AT-III)抗凝血酶III                                                                     |         |         |       |          |      |
| Sysmex CA 1500/6000/7000<br>Siemens Berichrom Antithrombin III (AT□) \ Nephelometry | %       | 28      | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil Liquid Antithrombin (AT□) \ Nephelometry        | %       | 26      | 10%   | <5       | <100 |
| (D-D(DDU)) D-二聚体(DDU)                                                               |         |         |       |          |      |
| Werfen ACL series<br>Werfen Hemosil D-Dimer HS \ Nephelometry                       | mg/L    | 1.12    | 10%   | <5       | <100 |
| Werfen ACL series<br>Wenfen Hemosil D-Dimer \ Nephelometry                          | mg/L    | 1.05    | 10%   | <5       | <100 |
| (D-D(FEU)) D-二聚体(FEU)                                                               |         |         |       |          |      |
| Sysmex CA 1500/6000/7000<br>Siemens Innovance Antithrombin (AT□) \ Nephelometry     | mg/L    | 3.62    | 10%   | <5       | <100 |
| (FIB)纤维蛋白原                                                                          |         |         |       |          |      |
| Sysmex CA 1500/6000/7000<br>Siemens Multifibren U (FIB) \ Light scattering          | g/L     | 1.75    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil Fib-C \ Light scattering                        | g/L     | 1.11    | 10%   | <5       | <100 |
| Stago STA series<br>Stago \ Magnetic Beads                                          | g/L     | 1.36    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil QFA Thrombin (FIB) \ Light scattering           | g/L     | 1.22    | 10%   | <5       | <100 |
| (Protein C)蛋白C                                                                      |         |         |       |          |      |
| Werfen ACL series<br>Werfen \ Light scattering                                      | IU/mL   | 0.25    | 12%   | <5       | <100 |
| (Protein S)蛋白S                                                                      |         |         |       |          |      |
| Werfen ACL series<br>Werfen \ Light scattering                                      | IU/mL   | 0.23    | 8.7%  | <5       | <100 |
| (PT)血浆凝血酶原时间                                                                        |         |         |       |          |      |
| Stago STA series<br>Stago \ Magnetic Beads                                          | Seconds | 46.2    | 10%   | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens PT Innovin \ Light scattering                   | INR     | 3.61    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT Recombiplastin 2G \ Light scattering         | Seconds | 46.4    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT Recombiplastin \ Light scattering            | Seconds | 42.7    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT Recombiplastin \ Light scattering            | INR     | 3.91    | 10%   | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens Thromborel S (PT) \ Light scattering            | INR     | 4.06    | 10.1% | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT Recombiplastin 2G \ Light scattering         | INR     | 3.95    | 10.1% | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens Thromborel S (PT) \ Light scattering            | Seconds | 44      | 10%   | <5       | <100 |
| Sysmex CA 1500/6000/7000<br>Siemens PT Innovin \ Light scattering                   | Seconds | 38.1    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT-Fib HS Plus \ Light scattering               | INR     | 4.55    | 10.1% | <5       | <100 |
| Stago STA series<br>Stago \ Magnetic Beads                                          | INR     | 5.18    | 10%   | <5       | <100 |
| Werfen ACL series<br>Werfen Hemosil PT-Fib HS Plus \ Light scattering               | Seconds | 55.8    | 10%   | <5       | <100 |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。