

REF	586-1	6×8mL	EXP	2021-11-18	LOT	37E598F	更新日期：	2019-01
1 / 1								

批号：37E598F			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AMY)淀粉酶			
Beckman AU series Dedicated \ PNP-G5	U/L	41.4	34.0-48.8
Roche cobas c 501/502/701 Dedicated \ PNP-G7	U/L	39.2	32.4-46.0
(BUN)尿素氮			
Beckman AU series Dedicated \ Dehydrogenase	mg/dL	410	338-482
(Ca)钙			
Roche cobas c 501/502/701 Dedicated \ ISE (direct)	mg/dL	7.76	6.50-9.02
(Cl)氯			
Roche cobas c 501/502/701 Dedicated \ ISE (indirect)	mEq/L	79.7	67.1-92.3
(Cortisol)皮质醇			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/dL	2.65	2.21-3.09
(CRE)肌酐			
Beckman UniCel Dx C 600/800 Dedicated \ Enzymatic	mg/dL	66.8	54.4-79.2
(GLU)葡萄糖			
Beckman AU series MedicalSystem \ Hexokinase	mg/dL	24.5	20.3-28.7
(K)钾			
Roche cobas c 501/502/701 Dedicated \ ISE (indirect)	mEq/L	30.4	29.8-31.0
(MA)微量白蛋白			
Beckman UniCel Dx C 600/800 3V-Diagnostic \ Bromocresol green	mg/L	26.4	22.0-30.8
Roche cobas c 501/502/701 Dedicated \ Biuret-endpoint	mg/L	25.8	21.4-30.2
(Mg)镁			
Beckman UniCel Dx C 600/800 Dedicated \ Calmagite colorimetry	mg/dL	5.96	5.10-6.82
(Na)钠			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	77.4	65.0-89.8
Roche cobas c 501/502/701 Dedicated \ ISE (direct)	mEq/L	75.4	63.2-87.6
(P)磷			
Roche cobas c 501/502/701 Dedicated \ Molybdenum blue colorimetry	mg/dL	27.2	22.6-31.8
(PRO/ TP)(总) 蛋白			
Roche cobas c 501/502/701 Dedicated \ Biuret-rate	mg/dL	15.2	12.6-17.8
(UA)尿酸			
Roche cobas c 501/502/701 Dedicated \ Uricase-UV	mg/dL	8.25	6.79-9.71

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。