

REF	586-2	6×8mL	EXP	2021-11-19	LOT	37E809F	更新日期：	2019-01
								1 / 1

批号：37E809F			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AMY)淀粉酶			
Beckman AU series Dedicated \ PNP-G5	U/L	199	165-233
Roche cobas c 501/502/701 Dedicated \ PNP-G7	U/L	192	158-226
(BUN)尿素氮			
Beckman AU series Dedicated \ Dehydrogenase	mmol/L	839	687-991
Roche cobas 4000/6000/8000 Dedicated \ Urease electrode	mg/dL	824	668-980
(Ca)钙			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	11.7	9.9-13.5
Roche cobas c 501/502/701 Dedicated \ ISE (direct)	mg/dL	11.9	9.7-14.1
(Cl)氯			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	179	151-207
Roche cobas c 501/502/701 Dedicated \ ISE (indirect)	mmol/L	171	143-199
(Cortisol)皮质醇			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/dL	25.7	21.3-30.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/dL	24.8	20.4-29.2
(CRE)肌酐			
Beckman AU series Dedicated \ Enzymatic	mg/dL	163	133-193
Roche cobas c 501/502/701 MedicalSystem \ Enzymatic	mmol/L	15.1	12.5-17.7
(GLU)葡萄糖			
Beckman AU series MedicalSystem \ Hexokinase	mg/dL	277	225-329
Roche cobas c 501/502/701 Dedicated \ Hexokinase	mmol/L	15.4	12.8-18.0
(K)钾			
Beckman AU series Dedicated \ ISE (direct)	mEq/L	78.2	63.8-92.6
Beckman UniCel Dx C 600/800 3V-Diagnostic \ ISE (direct)	mmol/L	76.4	62.6-90.2
(MA)微量白蛋白			
Beckman AU series Dedicated \ Biuret-rate	mg/L	59.5	48.1-70.9
Roche cobas c 501/502/701 Dedicated \ Biuret-endpoint	mg/L	58.4	47.2-69.6
(Mg)镁			
Beckman AU series Dedicated \ Calmagite colorimetry	mg/dL	12.6	10.4-14.8
Roche cobas c 501/502/701 Dedicated \ Calmagite colorimetry	mmol/L	12.4	10.2-14.6
(Na)钠			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	166	140-192
Roche cobas c 501/502/701 Dedicated \ ISE (direct)	mEq/L	164	138-190
(P)磷			
Beckman AU series Dedicated \ Phosphomolybdate-UV-end point	mmol/L	47.4	41.0-53.8
Roche cobas c 501/502/701 Dedicated \ Molybdenum blue colorimetry	mg/dL	47.3	40.9-53.7
(PRO/ TP)(总) 蛋白			
Beckman AU series Dedicated \ Biuret-endpoint	mg/dL	69.1	56.9-81.3
Roche cobas c 501/502/701 Dedicated \ Biuret-rate	mg/L	709	575-843
(UA)尿酸			
Beckman AU series Dedicated \ Uricase colorimetry	mg/dL	14.7	12.1-17.3
Roche cobas c 501/502/701 Dedicated \ Uricase-UV	mmol/L	0.79	0.65-0.93

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。