

REF	510-2	6×3mL	EXP	2022-12-21	LOT	30A8984	更新日期：	2019-01
1 / 2								

批号：30A8984			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AFP)甲胎蛋白			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	94.3	78.1-110.5
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	94.3	78.1-110.5
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	78.7	65.3-92.1
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	334	208-460
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	284	176-392
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	182	114-250
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	U/mL	217	135-299
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	114.1	70.7-157.5
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	97.8	60.6-135.0
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	75.8	47.0-104.6
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	U/mL	125.1	77.7-172.5
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	482	414-550
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	237	203-271
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	100	86-114
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	U/mL	204	176-232
(CA 72-4)糖类抗原 72-4			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	13.4	8.4-18.4
(CEA)癌胚抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	68.1	58.5-77.7
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	46.7	40.1-53.3
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	51.1	43.9-58.3
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	53.9	46.3-61.5
(Cyfra 21-1)Cyfra 21-1(细胞角蛋白19片段)			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.2	7.6-16.8
(Ferritin)铁蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	488	422-554
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	293	253-333
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	304	262-346
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	307	265-349
(F-PSA)游离PSA			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.64	3.20-6.08
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.1	4.9-9.3

REF	510-2	6×3mL	EXP	2022-12-21	LOT	30A8984	更新日期：	2019-01
2 / 2								

		批号：30A8984	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.88	3.36-6.40
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.12	3.52-6.72
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	19.7	16.5-22.9
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	20.1	16.9-23.3
Roche Modular E170 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.9	15.9-21.9
Siemens ADVIA Centaur XP Siemens \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.1	13.5-18.7
(SCC)鳞状细胞癌抗原			
Abbott AxSYM Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	15.5	9.5-21.5
(TPA)组织多肽抗原			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	8030	4,980-11,080
(TPS)组织多肽特异抗原			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/L	262	162-362
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	153	131-175
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	110	94-126

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。