

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 510-1 | 6×3mL | <b>EXP</b> | 2022-12-31 | <b>LOT</b> | 30A6274 | 更新日期： | 2019-01 |
| 1 / 2      |       |       |            |            |            |         |       |         |

|                                                                                      |       |      |            |
|--------------------------------------------------------------------------------------|-------|------|------------|
| 批号：30A6274                                                                           |       |      |            |
| 项目\仪器\试剂\方法                                                                          | 单位    | 均值   | + / - 2 SD |
| <b>(AFP)甲胎蛋白</b>                                                                     |       |      |            |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )      | ng/mL | 4.43 | 3.67-5.19  |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | ng/mL | 4.1  | 3.4-4.8    |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | ng/mL | 4    | 3-5        |
| <b>(CA 125)糖类抗原 125</b>                                                              |       |      |            |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | U/mL  | 44.5 | 27.5-61.5  |
| DiaSorin Liaison /XL<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )            | U/mL  | 33.3 | 20.5-46.1  |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | U/mL  | 22.3 | 13.7-30.9  |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | U/mL  | 27   | 17-37      |
| <b>(CA 15-3)糖类抗原 15-3</b>                                                            |       |      |            |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | U/mL  | 17.4 | 10.8-24.0  |
| DiaSorin Liaison /XL<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )            | U/mL  | 17.1 | 10.5-23.7  |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | U/mL  | 15.9 | 9.9-21.9   |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | U/mL  | 19   | 12-26      |
| <b>(CA 19-9)糖类抗原 19-9</b>                                                            |       |      |            |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | U/mL  | 48.4 | 41.6-55.2  |
| DiaSorin Liaison /XL<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )            | U/mL  | 28.9 | 24.9-32.9  |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | U/mL  | 14.6 | 12.6-16.6  |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | U/mL  | 24.9 | 21.3-28.5  |
| <b>(CA 72-4)糖类抗原 72-4</b>                                                            |       |      |            |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | U/mL  | 2.86 | 1.76-3.96  |
| <b>(CEA)癌胚抗原</b>                                                                     |       |      |            |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | ng/mL | 2.92 | -1.28-7.12 |
| Beckman UniCel Dxl 600/800<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )      | ng/mL | 2.07 | -0.93-5.07 |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | ng/mL | 2.6  | -0.4-5.6   |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA) | ng/mL | 2.2  | 1.6-2.8    |
| <b>(Cyfra 21-1)Cyfra 21-1(细胞角蛋白19片段)</b>                                             |       |      |            |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | ng/mL | 1.5  | 0.9-2.1    |
| <b>(Ferritin)铁蛋白</b>                                                                 |       |      |            |
| Abbott ARCHITECT i series<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | ng/mL | 76.1 | 65.7-86.5  |
| Beckman Access<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )                  | ng/mL | 49.5 | 42.7-56.3  |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)      | ng/mL | 70.3 | 60.7-79.9  |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )        | ng/mL | 57   | 49-65      |
| <b>(F-PSA)游离PSA</b>                                                                  |       |      |            |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )          | ng/mL | 0.82 | 0.56-1.08  |
| Beckman Access<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )            | ng/mL | 1.43 | 0.99-1.87  |

|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 510-1 | 6×3mL | EXP | 2022-12-31 | LOT | 30A6274 | 更新日期： | 2019-01 |
| 2 / 2 |       |       |     |            |     |         |       |         |

|                                                                                 |       |            |            |
|---------------------------------------------------------------------------------|-------|------------|------------|
|                                                                                 |       | 批号：30A6274 |            |
| 项目\仪器\试剂\方法                                                                     | 单位    | 均值         | + / - 2 SD |
| Roche Modular E170<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL | 1.03       | 0.71-1.35  |
| Siemens ADVIA Centaur XP<br>Siemens \ Chemiluminescence Immunoassay (CLIA )     | ng/mL | 1.03       | 0.71-1.35  |
| (PSA)前列腺特异性抗原                                                                   |       |            |            |
| Abbott ARCHITECT i series<br>Abbott \ Chemiluminescence Immunoassay (CLIA )     | ng/mL | 2.66       | 2.24-3.08  |
| Beckman Access<br>Beckman Coulter \ Chemiluminescence Immunoassay (CLIA )       | ng/mL | 2.59       | 2.19-2.99  |
| Roche Modular E170<br>Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)     | ng/mL | 2.76       | 2.32-3.20  |
| Siemens ADVIA Centaur XP<br>Siemens \ Chemiluminescence Immunoassay (CLIA )     | ng/mL | 2.36       | 1.98-2.74  |
| (SCC)鳞状细胞癌抗原                                                                    |       |            |            |
| Abbott AxSYM<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)       | ng/mL | 1.7        | 1.0-2.4    |
| (TPA)组织多肽抗原                                                                     |       |            |            |
| DiaSorin Liaison /XL<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )       | U/L   | 377        | 233-521    |
| (TPS)组织多肽特异抗原                                                                   |       |            |            |
| Siemens Immulite1000/2000<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )  | U/L   | 28         | 17-39      |
| (β-hCG)人绒毛膜促性腺激素Beta亚基                                                          |       |            |            |
| Roche Modular E170<br>Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA) | IU/L  | 15.3       | 12.7-17.9  |
| Siemens ADVIA Centaur XP<br>Dedicated \ Chemiluminescence Immunoassay (CLIA )   | IU/L  | 12.5       | 10.3-14.7  |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。