

REF

539-2

6×3mL

EXP

2021-06-30

LOT

336C93A

更新日期： 2019-01

1 / 3

批号：336C93A			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AFP)甲胎蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	86.2	71.0-101.4
Beckman UniCel Dxi 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	117	95-139
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	109	87-131
(Anti-Tg)抗甲状腺球蛋白抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	37.6	30.2-45.0
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	16.7	13.5-19.9
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	kU/L	116	98-134
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	86.7	69.5-103.9
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	55.2	45.0-65.4
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	36.7	30.3-43.1
Mindray CL-2000i/1000i/6000i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	40.2	31.8-48.6
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	kU/L	64.3	52.1-76.5
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	78.3	63.1-93.5
Beckman UniCel Dxi 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	65.7	53.3-78.1
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	61.3	49.7-72.9
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	58.8	48.2-69.4
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	38.6	31.4-45.8
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	44.8	36.0-53.6
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	48.2	39.6-56.8
Beckman UniCel Dxi 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	82.7	66.9-98.5
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	48.2	39.0-57.4
(CEA)癌胚抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	57.9	47.3-68.5
Beckman UniCel Dxi 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.7	8.9-12.5
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	50.6	40.8-60.4
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	11.8	9.6-14.0
(C-P)C-肽			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.24	1.82-2.66
Roche cobas E601/602 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.19	1.77-2.61
(E2)雌二醇			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	754	610-898

REF

539-2

6×3mL

EXP

2021-06-30

LOT

336C93A

更新日期： 2019-01

2 / 3

批号：336C93A

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	898	720-1,076
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	678	554-802
(FA)叶酸			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.96	7.28-10.64
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	8.42	6.80-10.04
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.78	8.12-11.44
(Ferritin)铁蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	232	190-274
(F-PSA)游离PSA			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	6.22	4.98-7.46
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.03	4.01-6.05
(FT3)游离T3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.2	9.0-13.4
Mindray CL-2000i/1000i/6000i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.3	10.1-14.5
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	9.46	7.62-11.30
(FT4)游离T4			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.2	23.4-35.0
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.7	26.9-38.5
Mindray CL-2000i/1000i/6000i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.1	22.5-33.7
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	36.3	29.7-42.9
(HCY)同型半胱氨酸			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	umol/L	19.6	16.0-23.2
(LH)促黄体生成素			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	15.9	12.7-19.1
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	19.8	16.0-23.6
(PRL)催/泌乳素			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	17.9	14.5-21.3
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	18.8	15.2-22.4
(PROG)孕酮			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.4	8.2-12.6
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	32.6	26.4-38.8
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.32	7.46-11.18
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	28.4	23.2-33.6
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.4	9.2-13.6
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	14.5	12.1-16.9

REF	539-2	6×3mL	EXP	2021-06-30	LOT	336C93A	更新日期：	2019-01
								3 / 3

		批号：336C93A	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche Modular E170 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.8	9.6-14.0
(TES)睾酮			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.31	0.25-0.37
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.67	0.55-0.79
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	1.39	1.13-1.65
(Tg)甲状腺球蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.2	8.4-12.0
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	8.09	6.47-9.71
(TSH)促甲状腺素			
Abbott ARCHITECT i series Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.1	7.9-12.3
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.86	7.24-10.48
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.2	9.0-13.4
(TT3)总T3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.47	2.03-2.91
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	4.52	3.68-5.36
(TT4)总T4			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	125	103-147
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	141	113-169
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.32	0.26-0.38

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。