

REF	540-2	6×3mL	EXP	2023-01-25	LOT	338509D	更新日期：	2019-09
1 / 6								

批号：338509D			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	5.45	4.33-6.57
(AFP)甲胎蛋白			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	109	85-133
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	121	97-145
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	122	98-146
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	U/mL	92.6	74.0-111.2
(ALD)醛固酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	498	404-592
(Anti-Tg)抗甲状腺球蛋白抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	IU/mL	29.5	23.9-35.1
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	43.2	35.0-51.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	62.1	50.3-73.9
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	28.6	23.8-33.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	22.8	18.0-27.6
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	29.3	23.5-35.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	35.6	29.0-42.2
(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	10.5	8.1-12.9
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	58.5	45.9-71.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	60.4	47.8-73.0
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	60.4	48.8-72.0
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	U/mL	79.9	63.7-96.1
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	40.6	31.4-49.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	44.9	35.7-54.1
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	49.6	40.0-59.2
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	U/mL	59.5	47.3-71.7
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	34.9	26.7-43.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	49.1	38.9-59.3
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	51.2	41.0-61.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	U/mL	38.1	30.3-45.9

REF

540-2

6×3mL

EXP

2023-01-25

LOT

338509D

更新日期： 2019-09

2 / 6

批号：338509D

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(CEA)癌胚抗原			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	9.87	7.91-11.83
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	10.7	8.3-13.1
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/L	17.5	13.7-21.3
(CK-MB)肌酸激酶-MB			
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	9.0-13.4
Roche cobas c 501/502/701 Roche \ Chemiluminescence Immunoassay (CLIA)	µg/L	5.86	4.70-7.02
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.6	13.0-20.2
(Cortisol)皮质醇			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/dL	24.2	19.2-29.2
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/dL	23.9	18.9-28.9
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µmol/L	0.62	0.50-0.74
(C-P)C-肽			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	2.59	2.07-3.11
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µg/L	2.54	2.02-3.06
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.67	2.14-3.20
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.51	2.03-2.99
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.98	1.58-2.38
(DGXN)地高辛			
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	2.63	2.07-3.19
(DHEA-S)硫酸脱氢表雄酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	µmol/L	16.4	12.8-20.0
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	µmol/L	14.9	12.1-17.7
(E2)雌二醇			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	976	780-1,172
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	1098	882-1,314
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	pg/mL	168	134-202
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	1015	813-1,217
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	935	749-1,121
(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	48.6	38.8-58.4
(FA)叶酸			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.22	6.58-9.86
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.67	3.75-5.59
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	nmol/L	16.9	13.3-20.5

REF

540-2

6×3mL

EXP

2023-01-25

LOT

338509D

更新日期： 2019-09

3 / 6

批号：338509D

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(Ferritin)铁蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	206	164-248
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	146	118-174
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	224	180-268
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	139	111-167
(F-PSA)游离PSA			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.87	3.81-5.93
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	4.02-6.14
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.98	4.00-5.96
(FSH)促卵泡刺激素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	33.1	26.1-40.1
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	34.2	27.2-41.2
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	28.9	23.1-34.7
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	24.9	20.3-29.5
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	26.7	20.9-32.5
(FT3)游离T3			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.6	8.2-13.0
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.5	9.1-13.9
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.4	10.0-14.8
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	pg/mL	8.73	6.98-10.48
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	13.6	11.0-16.2
(FT4)游离T4			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	30.5	23.7-37.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.8	23.8-35.8
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	32.3	26.5-38.1
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/dL	2.55	2.04-3.06
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	36.9	29.7-44.1
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	34.6	27.4-41.8
(HCY)同型半胱氨酸			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	24.2	19.2-29.2
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.6	17.8-27.4
(HGH)人生长激素			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mU/L	11.4	9.0-13.8
(HS-CTnT)超敏肌钙蛋白T			
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.48	2.78-4.18

REF

540-2

6×3mL

EXP

2023-01-25

LOT

338509D

更新日期： 2019-09

4 / 6

批号：338509D

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(IgE)免疫球蛋白 IgE			
Beckman Immage 800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	153	123-183
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	kU/L	78.4	63.8-93.0
(IGF1)胰岛素生长因子1			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	22.8	18.2-27.4
(INS)胰岛素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	34.2	27.0-41.4
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.1	27.9-42.3
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	32.9	26.3-39.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	35.2	28.4-42.0
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.5	25.7-39.3
(LH)促黄体生成素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.7	16.5-24.9
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	19.7	15.5-23.9
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	14.3	11.4-17.2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	19.8	16.0-23.6
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.1	14.3-21.9
(Myo / Mb)肌红蛋白			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	121	99-143
HYBIOME EVERESYS Series Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	109	87-131
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	96.6	77.3-115.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	142	116-168
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	142	114-170
(NT-proBNP)N-末端脑钠肽前体			
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	pg/mL	1588	1,272-1,904
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	121	97-145
(PRL)催/ 泌乳素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	13.6-20.0
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.8	13.4-20.2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.2	14.6-21.8
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	392	314-470
(PROG)孕酮			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	32.5	25.3-39.7
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	34.2	27.0-41.4
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.9	9.5-14.3

REF

540-2

6×3mL

EXP

2023-01-25

LOT

338509D

更新日期： 2019-09

5 / 6

批号：338509D

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.09	7.27-10.91
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	nmol/L	37.2	29.4-45.0
(PSA)前列腺特异性抗原			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.5	9.7-15.3
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.7	9.9-15.5
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	11.2	8.6-13.8
(PTH)甲状旁激素			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	78.4	61.8-95.0
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	81.3	64.7-97.9
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	pg/mL	112	90-134
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.42	7.60-11.24
(TES)睾酮			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	1.88	1.50-2.26
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	1.84	1.46-2.22
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/dL	238	190-286
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.44	0.36-0.52
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	nmol/L	1.72	1.36-2.08
(Tg)甲状腺球蛋白			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.4	8.2-12.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.23	5.67-8.79
(THEO)茶碱			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	137	111-163
(TnI)肌钙蛋白I			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	0.43	0.33-0.53
(TnT)肌钙蛋白T			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	1024	820-1,228
(TSH)促甲状腺素			
Abbott ARCHITECT i series Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	9.67	7.71-11.63
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.1	7.7-12.5
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	11.6	9.2-14.0
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	10.7	8.6-12.8
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	8.7-13.1
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.81	7.81-11.81
(TT3)总三碘甲状腺原氨酸			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.37	2.71-4.03
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.65	2.89-4.41

REF	540-2	6×3mL	EXP	2023-01-25	LOT	338509D	更新日期：	2019-09
6 / 6								

批号：338509D			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.29	1.83-2.75
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	2.52	2.08-2.96
(TT4)总T4			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	147	121-173
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	137	111-163
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	µg/dL	13.7	11.0-16.4
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	121	99-143
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.1	8.7-13.5
(VB12)维生素B12			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	583	471-695
Siemens ADVIA Centaur CP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	pg/mL	402	320-484
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	467	371-563
(Vit D)维生素D			
HYBIOME EVERESYS Series HYBIOME \ Chemiluminescence Immunoassay (CLIA)	ng/mL	28.3	22.6-34.0
(β2M)β-2-微球蛋白			
Roche cobas 4000/6000/8000 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/L	1.16	0.90-1.42
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mg/L	1.48	1.16-1.80
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Abbott ARCHITECT i series Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	U/mL	126	98-154
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	133	105-161

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。