

REF	671-2	6×1mL	EXP	2023-01-13	LOT	4003357	更新日期：	2020-12
1 / 1								

批号：4003357			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(APTT)活化部分凝血活酶时间			
Stago STA series	Seconds	75	57-93
Stago \ Magnetic Beads			
Sysmex CS 5100	Seconds	98.6	80.2-117
Siemens \ Light scattering			
Sysmex CA 1500/6000/7000	Seconds	96.3	77.9-115
Siemens \ Light scattering			
IL ACL series	Seconds	68.4	55.2-81.6
IL \ Light scattering			
(AT-III)抗凝血酶III			
IL ACL series	%	46.9	38.3-55.5
Dedicated \ Nephelometry			
(D-D(DDU)) D-二聚体(DDU)			
Stago STA series	mg/L	2.68	2.04-3.32
SekiSui \ Nephelometry			
IL ACL series	mg/L	0.689	0.545-0.833
Dedicated \ Nephelometry			
(D-D(FEU)) D-二聚体(FEU)			
Sysmex CS 5100	mg/L	2.88	2.32-3.44
Dedicated \ Nephelometry			
Sysmex CA 1500/6000/7000	mg/L	2.96	2.4-3.52
SekiSui \ Nephelometry			
(FDP)纤维蛋白(原)降解产物			
Sysmex CA 1500/6000/7000	mg/L	5.89	4.71-7.07
SekiSui \ Nephelometry			
IL ACL series	mg/L	4970	3990-5950
SekiSui \ Nephelometry			
Stago STA series	mg/L	7.24	5.5-8.98
SekiSui \ Nephelometry			
(FIB)纤维蛋白原			
Stago STA series	g/L	2.54	1.93-3.15
Stago \ Magnetic Beads			
Sysmex CS 5100	g/L	1.87	1.53-2.21
Siemens \ Light scattering			
Sysmex CA 1500/6000/7000	g/L	1.96	1.6-2.32
Siemens \ Light scattering			
IL ACL series	g/L	2.21	1.83-2.59
IL \ Light scattering			
(PT INR)PT国际标准化比值			
Stago STA series	INR	2.9	2.2-3.6
Stago \ Calculate			
(PT)血浆凝血酶原时间			
Stago STA series	Seconds	29.6	22.5-36.7
Stago \ Magnetic Beads			
Sysmex CS 5100	Seconds	32.3	26.5-38.1
Siemens \ Light scattering			
Sysmex CA 1500/6000/7000	Seconds	32.1	25.9-38.3
Siemens \ Light scattering			
IL ACL series	Seconds	30.6	24.8-36.4
IL \ Light scattering			
(TT)凝血酶时间			
Stago STA series	Seconds	14.6	11.1-18.1
Stago \ Magnetic Beads			
Sysmex CA 1500/6000/7000	Seconds	17.6	14.2-21
Dedicated \ Light scattering			
IL ACL series	Seconds	17.2	13.8-20.6
SUNBIO \ Light scattering			

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。