

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2021-11
1 / 5								

批号：3384F15			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ng/mL	1.8	1.43-2.17
(AFP)甲胎蛋白			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	132	100-164
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	119	96.6-141
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	105	83-127
(ALD)醛固酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	pg/mL	179	146-212
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	pg/mL	2900	2320-3480
(Anti-Tg)抗甲状腺球蛋白抗体			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	57.1	45.9-68.3
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	20.9	16.7-25.1
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	23.7	18.9-28.5
(ASD)雄烯二酮			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.92	2.23-3.61
(CA 125)糖类抗原 125			
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	62.1	49.5-74.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	53.3	42.7-63.9
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	86.6	70.4-103
(CA 15-3)糖类抗原 15-3			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	58.7	46.5-70.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	51.4	40.6-62.2
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	45.8	36.6-55
(CA 19-9)糖类抗原 19-9			
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/mL	46.5	36.3-56.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	50.9	40.7-61.1
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	39.8	32-47.6
(CEA)癌胚抗原			
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.7	8.3-13.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.8	8.2-13.4
(CK)肌酸激酶			
Roche cobas 4000/6000/8000 Roche \ Urease-UV-rate	U/L	57.1	43.5-70.7
(CK-MB)肌酸激酶-MB (活性)			
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	U/L	9.78	7.58-12
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.4	8-12.8
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	16.2	12.8-19.6

REF

540-2

6×3mL

EXP

2023-08-26

LOT

3384F15

更新日期：

2021-11

2 / 5

批号：3384F15

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(Cortisol)皮质醇			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	207	164-250
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	233	177-289
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	231	181-281
(C-P)C-肽			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.33	1.87-2.79
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.23	1.71-2.75
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.76	1.36-2.16
(DGXN)地高辛			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.42	1.86-2.98
Beckman Access Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	2.09	1.65-2.53
(DHEA-S)硫酸脱氢表雄酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	μmol/L	17.4	13.8-21
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μmol/L	16.5	13.3-19.7
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	15.7	12.5-18.9
(E2)雌二醇			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	812	648-976
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	998	800-1190
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	978	782-1170
(EPO)促红细胞生成素			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/mL	56.1	44.5-67.7
(FA)叶酸			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.46	5.87-9.05
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.92	4.51-7.33
(Ferritin)铁蛋白			
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	μg/L	180	144-216
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	139	111-167
(F-PSA)游离前列腺特异性抗原			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.08	4.02-6.14
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	5.07	3.99-6.15
(FSH)促卵泡刺激素			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	33.7	26.7-40.7
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	31.8	24.2-39.4
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	27.1	21.3-32.9
(FT3)游离T3			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	14.1	10.7-17.5

REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2021-11
							3 / 5	

		批号：3384F15	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	11.5	9.1-13.9
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	10.6	8.47-12.7
(FT4)游离T4			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pmol/L	34.6	27.4-41.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	33.2	25.4-41
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pmol/L	29.8	23.8-35.8
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	pmol/L	26.9	21.8-32
(HCY)同型半胱氨酸			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.9	18.3-27.5
Roche cobas 4000/6000/8000 Roche \ Enzymatic Cycling Assays Method	μmol/L	25.1	19.1-31.1
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μmol/L	22.4	17.6-27.2
(HGH)人生长激素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	4.43	3.53-5.33
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	3.8	3-4.6
(IgE)免疫球蛋白 IgE			
Beckman Image 800 Dedicated \ Immunoturbidimetry (ITA)	IU/mL	85.8	70.6-101
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/mL	80.9	64.5-97.3
(IGF1)胰岛素生长因子1			
Autolumo A2000/A2000 PLUS Autobio \ Chemiluminescence Immunoassay (CLIA)	ng/mL	207	159-255
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	22.8	18.2-27.4
(INS)胰岛素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	31.4	24-38.8
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	mIU/L	35.1	27.9-42.3
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	32.5	25.7-39.3
(LH)促黄体生成素			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	IU/L	18.1	14.3-21.9
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	22.6	18.6-26.6
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	20.4	16.2-24.6
(Myo / Mb)肌红蛋白			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	μg/L	139	111-167
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	145	110-180
(NT-proBNP)N-末端脑钠肽前体			
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1240	944-1530
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	146	111-181
(PRL)催/ 泌乳素			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	18.4	14.7-22.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.2	14.6-21.8

REF

540-2

6×3mL

EXP

2023-08-26

LOT

3384F15

更新日期：

2021-11

4 / 5

批号：3384F15

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.3	13.7-20.9
(PROG)孕酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.8	8.54-13.1
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.8	8.22-13.4
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.7	9.25-14.2
(PSA)前列腺特异性抗原			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	13.6	10.3-16.9
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	11.2	8.6-13.8
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	12.7	9.9-15.5
(PTH)甲状旁腺素			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	84.2	67-101
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	pg/mL	81.3	64.7-97.9
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	91.1	72.1-110
(TES)睾酮			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.531	0.421-0.641
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.476	0.366-0.586
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.487	0.383-0.591
(Tg)甲状腺球蛋白			
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.6	8.4-12.8
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.4	10-14.8
(THEO)茶碱			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	µg/mL	24.7	20-29.4
(TnI)肌钙蛋白I			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.46	0.34-0.58
(TnT)肌钙蛋白T			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.05	0.846-1.25
Roche cobas e411 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1200	912-1480
(TSH)促甲状腺素			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	mIU/L	8.81	7.01-10.6
Beckman UniCel Dxl 600/800 Beckman Coulter \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	10.9	8.5-13.3
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.9	9.1-14.7
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.84	8.12-11.6
(TT3)总三碘甲状腺原氨酸			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.35	2.67-4.03
Beckman UniCel Dxl 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	4.62	3.86-5.38
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	3.76	2.86-4.66



REF	540-2	6×3mL	EXP	2023-08-26	LOT	3384F15	更新日期：	2021-11
5 / 5								

批号：3384F15			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(TT4)总T4			
Abbott ARCHITECT i series Abbott \ Chemiluminescence Immunoassay (CLIA)	nmol/L	170	137-203
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	nmol/L	138	112-164
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	140	106-174
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	nmol/L	0.011	0.009-0.013
(VB12)维生素B12			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	pg/mL	665	530-800
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	688	522-854
(Vit D)维生素D			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	85.1	67.3-103
(β2M)β-2-微球蛋白			
Tosoh AIA 1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	mg/L	1.46	1.14-1.78
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	mg/L	1.37	1.09-1.65
(β-hCG)人绒毛膜促性腺激素β亚基			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	154	121-187
Beckman UniCel Dxi 600/800 Beckman Coulter \ Chemiluminescence Immunoassay (CLIA)	IU/L	138	110-166

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。