



REF	662-3	6×4.5mL	EXP	2023-11-18	LOT	3F29FF4	更新日期：	2023-08
								1 / 3

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-3

6×4.5mL

EXP

2023-11-18

LOT

3F29FF4

更新日期：2023-08

2 / 3

批号：3F29FF4			
项目\仪器\试剂方法	单位	均值	+ / - 2 SD
(BASO %)嗜碱性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	50	30-70
Sysmex XN series sysmex \ Calculate	%	10	0-20
Sysmex XS series sysmex \ Calculate	%	10	0-20
Sysmex XT series sysmex \ Calculate	%	15	0-30
(EOS %)嗜酸性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	15	5-25
Sysmex XE series sysmex \ Calculate	%	15	0-30
Sysmex XN series sysmex \ Calculate	%	20	0-40
Sysmex XS series sysmex \ Calculate	%	15	5-25
Sysmex XT series sysmex \ Calculate	%	15	0-30
(HCT)血细胞比容/压积			
Mindray BC-6000 series Mindray \ Light scattering	%	53.8	43.4-64.2
Sysmex XE series sysmex \ Light scattering	%	53.5	43.1-63.9
Sysmex XN series sysmex \ Light scattering	%	53.7	43.5-63.9
Sysmex XS series sysmex \ Light scattering	%	53.2	42.8-63.6
Sysmex XT series sysmex \ Light scattering	%	53.4	42.8-64
(HGB/ Hb)血红蛋白			
Mindray BC-6000 series Mindray \ Colorimetry	g/L	187	151-223
Sysmex XE series sysmex \ Light scattering	g/L	185	149-221
Sysmex XN series sysmex \ Light scattering	g/L	186	150-222
Sysmex XS series sysmex \ Light scattering	g/L	184	150-218
Sysmex XT series sysmex \ Light scattering	g/L	187	151-223
(LYM %)淋巴细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	20	10-30
Sysmex XE series sysmex \ Calculate	%	25	10-40
Sysmex XN series sysmex \ Calculate	%	25	10-40
Sysmex XS series sysmex \ Calculate	%	20	10-30
Sysmex XT series sysmex \ Calculate	%	25	10-40
(MCH)平均红细胞血红蛋白含量			
Mindray BC-6000 series Mindray \ Light scattering	pg	33.5	27.3-39.7
Sysmex XE series sysmex \ Light scattering	pg	32.4	26.2-38.6
Sysmex XN series sysmex \ Light scattering	pg	32.6	26.2-39
Sysmex XS series sysmex \ Light scattering	pg	32.5	26.3-38.7
Sysmex XT series sysmex \ Light scattering	pg	32.6	26.2-39
(MCHC)平均红细胞血红蛋白浓度			
Mindray BC-6000 series Mindray \ Light scattering	g/L	326	262-390
Sysmex XE series sysmex \ Light scattering	g/L	323	259-387
Sysmex XN series sysmex \ Light scattering	g/L	325	259-391
Sysmex XS series sysmex \ Light scattering	g/L	323	259-387
Sysmex XT series sysmex \ Light scattering	g/L	322	258-386
(MCV)平均红细胞体积			
Mindray BC-6000 series Mindray \ Light scattering	fL	95	77-113
Sysmex XE series sysmex \ Light scattering	fL	96	78-114

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-3

6×4.5mL

EXP

2023-11-18

LOT

3F29FF4

更新日期：

2023-08

3 / 3

Sysmex XN series sysmex \ Light scattering	fL	95	77-113
Sysmex XS series sysmex \ Light scattering	fL	96	78-114
Sysmex XT series sysmex \ Light scattering	fL	94	76-112
(MONO %)单核细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	10	0-20
Sysmex XN series sysmex \ Calculate	%	5	0-10
Sysmex XS series sysmex \ Calculate	%	10	0-20
Sysmex XT series sysmex \ Calculate	%	5	0-10
(NEUT %)嗜中性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	60	40-80
Sysmex XE series sysmex \ Calculate	%	50	30-70
Sysmex XN series sysmex \ Calculate	%	55	35-75
Sysmex XS series sysmex \ Calculate	%	55	35-75
Sysmex XT series sysmex \ Calculate	%	60	40-80
(PLT)血小板			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	450	360-540
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E +9/L)	450	350-550
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	455	355-555
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E +9/L)	455	345-565
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	450	350-550
(RBC)红细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E +12/L)	5.43	4.37-6.49
Sysmex XE series sysmex \ Light scattering	T/L (1E +12/L)	5.46	4.4-6.52
Sysmex XN series sysmex \ Light scattering	T/L (1E +12/L)	5.48	4.36-6.6
Sysmex XS series sysmex \ Light scattering	T/L (1E +12/L)	5.51	4.45-6.57
Sysmex XT series sysmex \ Light scattering	T/L (1E +12/L)	5.47	4.35-6.59
(WBC)白细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	17.5	13.9-21.1
Sysmex XE series sysmex \ Light scattering	Giga/L (1E +9/L)	17.6	14-21.2
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	17.8	14.2-21.4
Sysmex XS series sysmex \ Light scattering	Giga/L (1E +9/L)	17.7	14.1-21.3
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	17.7	14.1-21.3

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。