

REF	500-1	6×3mL	EXP	2022-10-01	LOT	2FB1F06	更新日期：	2019-11
								1 / 2

批号：2FB1F06			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(C-P)C-肽			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.96	0.78-1.14
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	0.62	0.48-0.76
(E2)雌二醇			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	336	270-402
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	211	171-251
(FA)叶酸			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	4.31	3.45-5.17
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	8.96	7.12-10.80
(FSH)促卵泡刺激素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.33	3.45-5.21
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	6.53	5.21-7.85
(FT3)游离T3			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/dL	5.66	4.56-6.76
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.67	3.69-5.65
(FT4)游离T4			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/dL	16.9	13.7-20.1
(hCG)人绒毛膜促性腺激素			
Tosoh AIA1800 TOSOH \ Chemiluminescence Immunoassay (CLIA)	U/L	4.57	3.63-5.51
(HGH)人生长激素			
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	4.85	3.85-5.85
(INS)胰岛素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	64.6	52.0-77.2
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	μIU/mL	7.75	6.15-9.35
(LH)促黄体生成素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.56	3.64-5.48
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.18	2.50-3.86
(PRL)催/ 泌乳素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	8.26	6.62-9.90
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	200	160-240
(PROG)孕酮			
Roche cobas E601/602 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/dL	3.04	2.42-3.66
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.52	2.80-4.24
(PTH)甲状旁激素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	11.8	9.6-14.0
(TES)睾酮			
Genrui-PA120 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/dL	23.6	19.0-28.2

REF	500-1	6×3mL	EXP	2022-10-01	LOT	2FB1F06	更新日期：	2019-11
2 / 2								

批号：2FB1F06			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	27.1	21.5-32.7
(Tg)甲状腺球蛋白			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	10.3	8.1-12.5
(TSH)促甲状腺素			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	1.94	1.56-2.32
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	1.56	1.22-1.90
(TT3)总三碘甲状原氨酸			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	1.76	1.40-2.12
(TT4)总T4			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	76.4	60.2-92.6
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	15.2	12.2-18.2
(VB12)维生素B12			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/L	263	211-315
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	309	247-371
(Vit D)维生素D			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	43.5	34.9-52.1
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	2.74	2.20-3.28
Tosoh AIA1800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	4.44	3.48-5.40

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。