

REF	539-1	6×3mL	EXP	2022-02-28	LOT	336A28A	更新日期：	2018-08
								1 / 3

批号：336A28A			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AFP)甲胎蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	2.38	1.94-2.82
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	2.98	2.34-3.62
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	kU/L	3.08	2.44-3.72
(Anti-Tg)抗甲状腺球蛋白抗体			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	12.7	10.1-15.3
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	3.98	3.14-4.82
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	kU/L	32.6	26.2-39.0
(Anti-TPO)抗甲状腺过氧化物酶抗体			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	kU/L	10.8	8.6-13.0
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	kU/L	13.9	11.3-16.5
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	kU/L	14.3	8.1-20.5
(CA 125)糖类抗原 125			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	10.9	8.7-13.1
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	13.8	11.2-16.4
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	11.4	9.2-13.6
(CA 15-3)糖类抗原 15-3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	10.9	8.7-13.1
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	10.7	8.9-12.5
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	U/mL	10.3	8.1-12.5
(CA 19-9)糖类抗原 19-9			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	3.57	2.93-4.21
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	U/mL	5.86	4.80-6.92
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/mL	6.79	5.53-8.05
(CEA)癌胚抗原			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	1.48	1.20-1.76
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	1.28	1.02-1.54
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	ug/L	0.85	0.67-1.03
(C-P)C-肽			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.96	0.78-1.14
Roche Modular E170 Hyssen \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.14	0.92-1.36
(E2)雌二醇			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	318	256-380
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	252	204-300
(FA)叶酸			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	3.36	2.72-4.00

REF	539-1	6×3mL	EXP	2022-02-28	LOT	336A28A	更新日期：	2018-08
								2 / 3

		批号：336A28A	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	3.96	3.18-4.74
(Ferritin)铁蛋白			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	36.8	29.0-44.6
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	34.2	27.4-41.0
(F-PSA)游离PSA			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.152	0.112-0.192
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.132	0.092-0.172
(FT3)游离T3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.31	3.45-5.17
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.34	3.36-5.32
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	4.63	3.67-5.59
(FT4)游离T4			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.9	10.3-15.5
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	13.3	13.0-13.6
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	14.2	11.4-17.0
(HCY)同型半胱氨酸			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	umol/L	15.9	12.7-19.1
(LH)促黄体生成素			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	3.46	2.74-4.18
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.74	3.90-5.58
(PRL)催/ 泌乳素			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	8.11	6.37-9.85
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ug/L	7.94	6.30-9.58
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	8.66	6.88-10.44
(PROG)孕酮			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	2.18	1.74-2.62
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	2.35	1.89-2.81
Siemens ADVIA Centaur XP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	nmol/L	2.49	2.03-2.95
(PSA)前列腺特异性抗原			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.09	0.87-1.31
Roche Modular E170 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.81	0.63-0.99
(TES)睾酮			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.38	4.26-6.50
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.27	4.15-6.39
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	22.9	18.1-27.7

REF	539-1	6×3mL	EXP	2022-02-28	LOT	336A28A	更新日期：	2018-08
								3 / 3

批号：336A28A			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(Tg)甲状腺球蛋白			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	8.84	7.20-10.48
Roche cobas e411 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	7.92	6.30-9.54
(TSH)促甲状腺素			
Abbott ARCHITECT i series Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	2.29	1.81-2.77
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	2.62	2.08-3.16
(TT3)总T3			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	2.08	1.64-2.52
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	1.34	1.08-1.60
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	1.84	1.48-2.20
(TT4)总T4			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	49.3	38.9-59.7
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	58.7	46.5-70.9
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	52.8	41.0-64.6
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.046	0.034-0.058
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	0.013	0.007-0.019

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。