

REF	530-3	6×3mL	EXP	2021-06-25	LOT	32934AA	更新日期：	2020-09
1 / 1								

批号：32934AA			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(CK-MB)肌酸激酶-MB			
Roche Modular E170			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	U/L	38.8	31-46.6
Beckman AU series			
Beckman Coulter \ Enzyme immunoassay (EIA)	U/L	36.6	27.8-45.4
(D-D(DDU)) D-二聚体(DDU)			
Sysmex CA 1500/6000/7000			
sysmex \ Latex agglutination (LA)	mg/L	4.32	3.56-5.08
(FDP)纤维蛋白(原)降解产物			
Sysmex CA 1500/6000/7000			
sysmex \ Nephelometry	mg/L	12.3	9.9-14.7
(hs-CRP)超敏C反应蛋白			
Beckman AU series			
Beckman Coulter \ Immunoturbidimetry (ITA)	mg/L	9.05	6.89-11.2
(Myo / Mb)肌红蛋白			
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	μg/L	326	270-382
(NT-proBNP)N-末端脑钠肽前体			
Roche Modular E170			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pg/mL	1320	1070-1560
Roche Modular E170			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	155	127-183
(TnI)肌钙蛋白I			
Abbott ARCHITECT i series			
Abbott \ Chemiluminescence Immunoassay (CLIA)	ng/mL	4.88	4.04-5.72
(TnT)肌钙蛋白T			
Roche cobas E601/602			
Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.73	0.59-0.87

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。