

REF	500-2	6×3mL	EXP	2022-03-31	LOT	2FB473B	更新日期：	2018-11
								1 / 2

		批号：2FB473B	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(17-OHP)17-羟孕酮			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	3.87	3.11-4.63
(Cortisol)皮质醇			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/dL	27.9	22.7-33.1
(C-P)C-肽			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.88	1.56-2.20
(E2)雌二醇			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	786	642-930
(FA)叶酸			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	7.67	6.15-9.19
(FSH)促卵泡刺激素			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	21.8	17.6-26.0
(FT3)游离T3			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	28.7	23.5-33.9
(FT4)游离T4			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	9.98	8.02-11.94
(IgA)免疫球蛋白 IgA			
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	g/L	183	145-221
(INS)胰岛素			
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	32.4	26.2-38.6
(LH)促黄体生成素			
DiaSorin Liaison /XL Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/mL	18.2	15.0-21.4
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	17.6	14.4-20.8
(PRL)催/ 泌乳素			
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	346	282-410
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	15.8	12.6-19.0
(PROG)孕酮			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.09	5.73-8.45
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	9.76	7.92-11.60
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	26.3	21.5-31.1
(TES)睾酮			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.66	0.54-0.78
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.64	0.52-0.76
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	2.56	2.12-3.00
(TSH)促甲状腺素			
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	9.57	7.97-11.17
Siemens ADVIA Centaur CP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	9.08	7.56-10.60



REF	500-2	6×3mL	EXP	2022-03-31	LOT	2FB473B	更新日期：	2018-11
2 / 2								

批号：2FB473B			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(TT3)总T3			
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	3.35	2.69-4.01
(TT4)总T4			
Beckman UniCel DxI 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	136	114-158
(β-hCG)人绒毛膜促性腺激素Beta亚基			
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.14	0.10-0.18
Roche cobas 4000/6000/8000 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	229	187-271

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。