

REF

512-2

6×3mL

EXP

2021-11-23

LOT

30D965F

更新日期：

2020-11

1 / 2

批号：30D965F

项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(AFP)甲胎蛋白			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	149	119-179
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	142	114-170
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	ng/mL	124	99.2-149
(CA 125)糖类抗原 125			
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	IU/mL	163	130-196
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	84.7	67.7-102
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	115	93-137
(CA 15-3)糖类抗原 15-3			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	42.9	34.3-51.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	93.1	74.5-112
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	IU/mL	61.5	49.2-73.8
(CA 19-9)糖类抗原 19-9			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	157	127-187
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	81.7	65.3-98.1
(CA 72-4)糖类抗原 72-4			
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	IU/mL	35.6	28.5-42.7
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	IU/mL	27.9	21.3-34.5
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	31.5	25.1-37.9
Roche Hitachi/Hitachi Systems Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/mL	31.5	25.1-37.9
(CEA)癌胚抗原			
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	ng/mL	24.3	19.4-29.2
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	31.2	24.8-37.6
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	21.8	17.4-26.2
(Cyfra 21-1)细胞角蛋白19片段			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	39.9	31.9-47.9
Roche Hitachi/Hitachi Systems Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	39.7	31.7-47.7
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	40.5	30.7-50.3
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	ng/mL	41.4	33.1-49.7
(Ferritin)铁蛋白			
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	μg/L	274	219-329
Roche cobas 4000/6000/8000 Roche \ Immunoturbidimetry (ITA)	μg/L	301	241-361
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	μg/L	264	212-316
(F-PSA)游离前列腺特异性抗原			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	3.99	3.19-4.79
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	2.51	1.99-3.03

REF	512-2	6×3mL	EXP	2021-11-23	LOT	30D965F	更新日期：	2020-11
								2 / 2

		批号：30D965F	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	ng/mL	1.61	1.29-1.93
(hCG)人绒毛膜促性腺激素			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/L	116	92-140
(NSE)神经元特异性烯醇化酶			
Mindray CL-2000i/1000i/6000i series Mindray \ Chemiluminescence Immunoassay (CLIA)	ng/mL	17.4	13.2-21.6
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	11.7	9.3-14.1
Roche Hitachi/Hitachi Systems Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	12.2	9.6-14.8
(PSA)前列腺特异性抗原			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	22.3	17.7-26.9
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	18.7	14.9-22.5
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	ng/mL	10.4	8.32-12.5
(Tg)甲状腺球蛋白			
Beckman Access Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	50.2	40.2-60.2
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	49.5	39.5-59.5
(β2M)β-2-微球蛋白			
Roche cobas c 501/502/701/702 Roche \ Immunoturbidimetry (ITA)	mg/L	16.1	12.7-19.5
Siemens BN II Dedicated \ Chemiluminescence Immunoassay (CLIA)	mg/L	22.3	17.7-26.9
Caris200 WanTai(BeiJin) \ Chemiluminescence Immunoassay (CLIA)	mg/L	16.1	12.9-19.3
(β-hCG)人绒毛膜促性腺激素β亚基			
Roche cobas E601/602 Roche \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	123	99-147

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。