

REF	500-1	Level 1	6×3mL	EXP	2020-12-31	LOT	Level1	500131	更新日期：	2017-04
	500-2	Level 2	6×3mL				Level2	500132		
	500-3	Level 3	6×3mL				Level3	500133		
									1 / 2	

		批号：500131		批号：500132		批号：500133	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD	均值	+ / - 2 SD	均值	+ / - 2 SD
(17-OHP)17-羟孕酮							
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	1.56	1.34-1.78	3.86	3.20-4.52	10.2	8.4-12.0
(Cortisol)皮质醇							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/dL	11	9-13	29	25-33	34	28-40
(C-P)C-肽							
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	0.89	0.73-1.05	1.88	1.56-2.20	3.48	2.90-4.06
(E2)雌二醇							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	pmol/L	280	240-320	800	660-940	1800	1,480-2,120
(FA)叶酸							
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	3.93	3.11-4.75	7.93	6.41-9.45	10.8	9.0-12.6
(FSH)促卵泡刺激素							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.32	4.36-6.28	22.3	18.1-26.5	42.8	35.0-50.6
(FT3)游离T3							
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	12.1	9.9-14.3	29.9	24.7-35.1	45.6	37.8-53.4
(FT4)游离T4							
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	4.95	4.21-5.69	10.5	8.3-12.7	30.2	24.6-35.8
(IgA)免疫球蛋白 IgA							
Roche Modular E170 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	g/L	208	170-246	186	148-224	192	154-230
(INS)胰岛素							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/L	11.3	9.1-13.5	33.2	27.0-39.4	127	105-149
(LH)促黄体生成素							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	4.1	3.5-4.7	18	15-21	39	33-45
DiaSorin Liaison /XL Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mIU/mL	3.76	3.04-4.48	18.2	15.0-21.4	66.8	55.0-78.6
(PRL)催/ 泌乳素							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ng/mL	9.2	7.4-11.0	16	13-19	35	29-41
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	204	168-240	357	293-421	758	620-896
(PROG)孕酮							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	2.6	2.2-3.0	27	23-31	74	62-86
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.69	0.55-0.83	6.98	5.60-8.36	17.2	14.0-20.4
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	0.82	0.66-0.98	9.82	7.98-11.66	23.7	19.3-28.1
(TES)睾酮							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	nmol/L	25.3	20.9-29.7	2.52	2.12-2.92	10.3	8.1-12.5
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	7.62	6.34-8.90	0.64	0.52-0.76	3.46	2.82-4.10
DiaSorin Liaison /XL Dedicated \ Chemiluminescence Immunoassay (CLIA)	ng/mL	5.16	4.18-6.14	0.57	0.45-0.69	2.41	1.95-2.87
(TSH)促甲状腺素							
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	mIU/L	2.8	2.2-3.4	9.5	7.9-11.1	26	22-30
Siemens ADVIA Centaur CP Dedicated \ Chemiluminescence enzyme immunoassay (CLEIA)	mIU/L	2.53	2.07-2.99	8.93	7.41-10.45	25.8	21.4-30.2

REF	500-1	Level 1 6×3mL	EXP	2020-12-31	LOT	Level1	500131	更新日期：2017-04 2 / 2
	500-2	Level 2 6×3mL				Level2	500132	
	500-3	Level 3 6×3mL				Level3	500133	

		批号：500131		批号：500132		批号：500133	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD	均值	+ / - 2 SD	均值	+ / - 2 SD
(TT3)总T3							
Beckman UniCel Dxl 600/800 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	1.48	1.24-1.72	3.74	3.08-4.40	8.98	7.42-10.54
(TT4)总T4							
Beckman UniCel Dxl 600/800 Dedicated \ Electro-Chemiluminescence Immunoassay (CLIA)	nmol/L	64.4	53.6-75.2	138	116-160	167	141-193
(β-hCG)人绒毛膜促性腺激素Beta亚基							
Roche cobas 4000/6000/800 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	IU/L	5.23	4.27-6.19	233	195-271	665	547-783
Abbott ARCHITECT i series Dedicated \ Chemiluminescence Immunoassay (CLIA)	IU/mL	0.012	0.006-0.018	0.136	0.110-0.162	0.369	0.303-0.435

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。