

REF	581-1	6×5mL	EXP	2022-01-31	LOT	376B8C9	更新日期：	2018-05
1 / 3								

批号：376B8C9			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(a-HBDH)α-羟丁酸脱氢酶			
Beckman AU series Dedicated \ immunoturbidimetric	U/L	136	112-160
(ALB)白蛋白			
Beckman AU series Dedicated \ Bromocresol green	g/L	42.2	34.0-50.4
(ALP)碱性磷酸酶			
Roche cobas c 501/502/701 Dedicated \ Rate -AMP buffer	U/L	79.7	65.5-93.9
(ALT/ SGPT)丙氨酸氨基转移酶			
Beckman AU series Dedicated \ UV with P5P	U/L	44.3	36.1-52.5
(AMY)淀粉酶			
Beckman UniCel DxC 600/800 Dedicated \ Maltotetraose-rate	U/L	103	85-121
(Apo-A1)载脂蛋白A1			
Beckman AU series LEADMAN \ immunoturbidimetric	g/L	1.78	1.38-2.18
(Apo-B)载脂蛋白B			
Beckman AU series LEADMAN \ immunoturbidimetric	g/L	0.66	0.44-0.88
(AST/ SGOT)天门冬氨酸氨基转移酶			
Beckman AU series Dedicated \ UV with P5P	U/L	41.2	33.8-48.6
(Ca)钙			
Beckman AU series Dedicated \ Arsenazo III	mmol/L	2.2	1.8-2.6
(CHE)胆碱酯酶			
Siemens ADVIA 1800/2400 LEADMAN \ Butyrylthiocholine Method	U/L	8094	6,736-9,452
(CHOL/ TC)总胆固醇			
Beckman AU series Dedicated \ Cholesterol oxidase	mmol/L	3.98	3.32-4.64
(CK)肌酸激酶			
Beckman AU series Beckman Coulter \ Rate -AMP buffer	U/L	90	72-108
(CK-MB)肌酸激酶-MB			
Roche cobas c 501/502/701 Dedicated \ Immune suppression	U/L	16.1	13.1-19.1
(Cl)氯			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	117.6	99.6-135.6
(CRE)肌酐			
Beckman AU series Dedicated \ Enzymatic	umol/L	77.4	64.2-90.6
(CRP)C反应蛋白			
Beckman AU series Dedicated \ Immunonephelometry (INA)	mg/L	7.4	6.0-8.8
(Cu)铜			
J&J Vitros Fusion 5.1 FS Dedicated \ dry chemistry	umol/L	20.4	16.2-24.6
(DBIL)直接胆红素			
Roche cobas c 501/502/701 Dedicated \ Diazotization -Blank Compensation	umol/L	5.5	4.4-6.6
(Digoxin)地高辛			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	0.79	0.65-0.93
(E2)雌二醇			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	pmol/L	86.6	70.4-102.8
(Fe)铁			
Beckman AU series Dedicated \ Ferrozine colorimetry	umol/L	17.6	14.2-21.0
(GGT)γ-谷氨酰基转移酶			
Beckman AU series Dedicated \ γ-glutamyl -3-carboxy-4-nitroaniline	U/L	36	29-43
(GLDH)谷氨酸脱氢酶			
Beckman UniCel DxC 600/800 LEADMAN \ Direct Measure	U/L	5.32	4.38-6.26

REF	581-1	6×5mL	EXP	2022-01-31	LOT	376B8C9	更新日期：	2018-05
								2 / 3

批号：376B8C9			
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(GLU)葡萄糖			
Beckman AU series Dedicated \ Hexokinase	mmol/L	4.11	3.47-4.75
(HCY)同型半胱氨酸			
ADC Elisa Series MAKER \ Enzymatic Cycling Assays Method	umol/L	13.42	10.62-16.22
(HDL-C)高密度脂蛋白胆固醇			
Beckman AU series Dedicated \ Direct Measure	mmol/L	0.91	0.73-1.09
(IgA)免疫球蛋白 IgA			
Beckman AU series Dedicated \ Immunonephelometry (INA)	g/L	2.2	1.5-2.9
(IgG)免疫球蛋白 IgG			
Beckman AU series Dedicated \ Immunonephelometry (INA)	g/L	11.08	8.68-13.48
(IgM)免疫球蛋白 IgM			
Beckman AU series Dedicated \ Immunonephelometry (INA)	g/L	0.9	0.6-1.2
(K)钾			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	3.62	3.04-4.20
(LA)乳酸			
Beckman UniCel Dx C 600/800 Dedicated \ Lactate to pyruvate (L-P)	mmol/L	2.12	1.74-2.50
(LDH)乳酸脱氢酶			
Beckman AU series Dedicated \ Lactate to pyruvate (L-P)	U/L	130	106-154
(LDL-C)低密度脂蛋白胆固醇			
Beckman AU series Dedicated \ Direct Measure	mmol/L	2.73	2.27-3.19
(Li)锂			
Beckman UniCel Dx C 600/800 Dedicated \ ISE (indirect)	mmol/L	0.79	0.65-0.93
(LIP)脂肪酶			
Siemens ADVIA 1800/2400 LEADMAN \ Enzymatic Cycling Assays Method	U/L	48	38-58
(Mg)镁			
Beckman AU series Dedicated \ Xylidyl Blue	mmol/L	0.89	0.71-1.07
(MYO / Mb)肌红蛋白			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	ug/L	24.13	19.93-28.33
(Na)钠			
Beckman AU series Dedicated \ ISE (indirect)	mmol/L	140	126-154
(NEFA)游离脂肪酸			
Siemens ADVIA 1800/2400 Co-Health \ ACS - ACOD	umol/L	602	487-717
(P)磷			
Beckman AU series Dedicated \ Phosphomolybdate UV-rate	mmol/L	1.02	0.86-1.18
(PRO/ TP)(总) 蛋白			
Beckman AU series Dedicated \ Biuret-endpoint	g/L	65.4	53.2-77.6
(TBA)总胆汁酸			
Beckman AU series Dedicated \ Enzymatic Cycling Assays Method	umol/L	17.2	8.2-26.2
(TBIL)总胆红素			
Beckman AU series Dedicated \ Diazotization	umol/L	17.9	15.1-20.7
(TES)睾酮			
Siemens Immulite1000/2000 Dedicated \ Chemiluminescence Immunoassay (CLIA)	nmol/L	10.3	8.1-12.5
(TG)甘油三酯			
Beckman AU series Dedicated \ GPO-POD	mmol/L	1.11	0.93-1.29



REF	581-1	6×5mL	EXP	2022-01-31	LOT	376B8C9	更新日期：	2018-05
3 / 3								

		批号：376B8C9	
项目\仪器\试剂\方法	单位	均值	+ / - 2 SD
(THEO)茶碱			
Roche cobas E601/602 Dedicated \ Electro-Chemiluminescence Immunoassay (ECLIA)	mg/L	13.81	11.41-16.21
(TRF)转铁蛋白			
Beckman AU series Dedicated \ Immunonephelometry (INA)	g/L	2.95	1.75-4.15
(UA)尿酸			
Beckman AU series Beckman Coulter \ Uricase colorimetry	umol/L	308	252-364
(UIBC)不饱和铁结合力			
Beckman AU series Dedicated \ Ferrozine colorimetry	umol/L	45.3	44.5-46.1
(Urea)尿素			
Beckman AU series Dedicated \ Urease-UV-rate	mmol/L	4.21	3.43-4.99
(Zn)锌			
J&J Vitros Fusion 5,1 FS Dedicated \ dry chemistry	umol/L	27.5	22.1-32.9

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。