

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2022-02-28 | <b>LOT</b> | 376DFDA | 更新日期： | 2018-05 |
|            |       |       |            |            |            |         |       | 1 / 3   |

|                                                                                      |        |       |              |
|--------------------------------------------------------------------------------------|--------|-------|--------------|
| 批号：376DFDA                                                                           |        |       |              |
| 项目\仪器\试剂\方法                                                                          | 单位     | 均值    | + / - 2 SD   |
| <b>(a-HBDH)α-羟丁酸脱氢酶</b>                                                              |        |       |              |
| Beckman AU series<br>Dedicated \ immunoturbidimetric                                 | U/L    | 336   | 278-394      |
| <b>(ALB)白蛋白</b>                                                                      |        |       |              |
| Beckman AU series<br>Dedicated \ Bromocresol green                                   | g/L    | 59.4  | 48.2-70.6    |
| <b>(ALP)碱性磷酸酶</b>                                                                    |        |       |              |
| Roche cobas c 501/502/701<br>Dedicated \ Rate -AMP buffer                            | U/L    | 248   | 208-288      |
| <b>(ALT/ SGPT)丙氨酸氨基转移酶</b>                                                           |        |       |              |
| Beckman AU series<br>Dedicated \ UV with P5P                                         | U/L    | 151   | 125-177      |
| <b>(AMY)淀粉酶</b>                                                                      |        |       |              |
| Beckman UniCel DxC 600/800<br>Dedicated \ Maltotetraose-rate                         | U/L    | 391   | 321-461      |
| <b>(Apo-A1)载脂蛋白A1</b>                                                                |        |       |              |
| Beckman AU series<br>LEADMAN \ immunoturbidimetric                                   | g/L    | 2.7   | 2.2-3.2      |
| <b>(Apo-B)载脂蛋白B</b>                                                                  |        |       |              |
| ADC Elisa Series<br>LEADMAN \ immunoturbidimetric                                    | g/L    | 0.99  | 0.75-1.23    |
| <b>(AST/ SGOT)天门冬氨酸氨基转移酶</b>                                                         |        |       |              |
| Beckman AU series<br>Dedicated \ UV with P5P                                         | U/L    | 259   | 221-297      |
| <b>(Ca)钙</b>                                                                         |        |       |              |
| Beckman AU series<br>Dedicated \ Arsenazo III                                        | mmol/L | 3.15  | 2.59-3.71    |
| <b>(CHE)胆碱酯酶</b>                                                                     |        |       |              |
| Siemens ADVIA 1800/2400<br>LEADMAN \ Butyrylthiocholine Method                       | U/L    | 10513 | 8,561-12,465 |
| <b>(CHOL/ TC)总胆固醇</b>                                                                |        |       |              |
| Beckman AU series<br>Dedicated \ Cholesterol oxidase                                 | mmol/L | 6.29  | 5.21-7.37    |
| <b>(CK)肌酸激酶</b>                                                                      |        |       |              |
| Beckman AU series<br>Beckman Coulter \ Rate -AMP buffer                              | U/L    | 302   | 240-364      |
| <b>(CK-MB)肌酸激酶-MB</b>                                                                |        |       |              |
| Roche cobas c 501/502/701<br>Dedicated \ Immune suppression                          | U/L    | 51.8  | 42.0-61.6    |
| <b>(Cl)氯</b>                                                                         |        |       |              |
| Beckman AU series<br>Dedicated \ ISE (indirect)                                      | mmol/L | 133.7 | 115.7-151.7  |
| <b>(CRE)肌酐</b>                                                                       |        |       |              |
| Beckman AU series<br>Dedicated \ Enzymatic                                           | umol/L | 254   | 208-300      |
| <b>(CRP)C反应蛋白</b>                                                                    |        |       |              |
| Beckman AU series<br>Dedicated \ Immunonephelometry (INA)                            | mg/L   | 2.7   | 2.3-3.1      |
| <b>(Cu)铜</b>                                                                         |        |       |              |
| J&J Vitros Fusion 5.1 FS<br>Dedicated \ dry chemistry                                | umol/L | 33.2  | 27.4-39.0    |
| <b>(DBIL)直接胆红素</b>                                                                   |        |       |              |
| Roche cobas c 501/502/701<br>Dedicated \ Diazotization -Blank Compensation           | umol/L | 15.4  | 12.6-18.2    |
| <b>(Digoxin)地高辛</b>                                                                  |        |       |              |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence<br>Immunoassay (ECLIA) | ug/L   | 1.76  | 1.46-2.06    |
| <b>(E2)雌二醇</b>                                                                       |        |       |              |
| Siemens Immulite1000/2000<br>Dedicated \ Chemiluminescence Immunoassay<br>(CLIA )    | pmol/L | 142   | 116-168      |
| <b>(Fe)铁</b>                                                                         |        |       |              |
| Beckman AU series<br>Dedicated \ Ferrozine colorimetry                               | umol/L | 32    | 26-38        |
| <b>(GGT)γ-谷氨酰基转移酶</b>                                                                |        |       |              |
| Beckman AU series<br>Dedicated \ γ-glutamyl -3-carboxy-4-nitroaniline                | U/L    | 152   | 124-180      |
| <b>(GLDH)谷氨酸脱氢酶</b>                                                                  |        |       |              |
| Beckman UniCel DxC 600/800<br>LEADMAN \ Direct Measure                               | U/L    | 9.66  | 8.02-11.30   |

|            |       |       |            |            |            |         |       |         |
|------------|-------|-------|------------|------------|------------|---------|-------|---------|
| <b>REF</b> | 581-2 | 6×5mL | <b>EXP</b> | 2022-02-28 | <b>LOT</b> | 376DFDA | 更新日期： | 2018-05 |
|            |       |       |            |            |            |         |       | 2 / 3   |

|                                                                                      |        |       |             |
|--------------------------------------------------------------------------------------|--------|-------|-------------|
| 批号：376DFDA                                                                           |        |       |             |
| 项目\仪器\试剂\方法                                                                          | 单位     | 均值    | + / - 2 SD  |
| <b>(GLU)葡萄糖</b>                                                                      |        |       |             |
| Beckman AU series<br>Dedicated \ Hexokinase                                          | mmol/L | 11.28 | 9.48-13.08  |
| <b>(HCY)同型半胱氨酸</b>                                                                   |        |       |             |
| Beckman AU series<br>MAKER \ Enzymatic Cycling Assays Method                         | umol/L | 18.88 | 15.08-22.68 |
| <b>(HDL-C)高密度脂蛋白胆固醇</b>                                                              |        |       |             |
| Beckman AU series<br>Dedicated \ Direct Measure                                      | mmol/L | 1.52  | 1.26-1.78   |
| <b>(IgA )免疫球蛋白 IgA</b>                                                               |        |       |             |
| Beckman AU series<br>Dedicated \ Immunonephelometry (INA)                            | g/L    | 3.1   | 2.3-3.9     |
| <b>(IgG)免疫球蛋白 IgG</b>                                                                |        |       |             |
| Beckman AU series<br>Dedicated \ Immunonephelometry (INA)                            | g/L    | 15.67 | 12.27-19.07 |
| <b>(IgM)免疫球蛋白 IgM</b>                                                                |        |       |             |
| Beckman AU series<br>Dedicated \ Immunonephelometry (INA)                            | g/L    | 1.4   | 0.9-1.9     |
| <b>(K)钾</b>                                                                          |        |       |             |
| Beckman AU series<br>Dedicated \ ISE (indirect)                                      | mmol/L | 5.85  | 4.83-6.87   |
| <b>(LA)乳酸</b>                                                                        |        |       |             |
| Beckman UniCel Dx C 600/800<br>Dedicated \ Lactate to pyruvate (L-P)                 | mmol/L | 6.03  | 4.91-7.15   |
| <b>(LDH)乳酸脱氢酶</b>                                                                    |        |       |             |
| Beckman AU series<br>Dedicated \ Lactate to pyruvate (L-P)                           | U/L    | 319   | 263-375     |
| <b>(LDL-C)低密度脂蛋白胆固醇</b>                                                              |        |       |             |
| Beckman AU series<br>Dedicated \ Direct Measure                                      | mmol/L | 4.3   | 1.9-6.7     |
| <b>(Li)锂</b>                                                                         |        |       |             |
| Beckman UniCel Dx C 600/800<br>Dedicated \ ISE (indirect)                            | mmol/L | 1.90  | 1.58-2.22   |
| <b>(LIP)脂肪酶</b>                                                                      |        |       |             |
| Siemens ADVIA 1800/2400<br>LEADMAN \ Enzymatic Cycling Assays Method                 | U/L    | 118   | 100-136     |
| <b>(Mg)镁</b>                                                                         |        |       |             |
| Beckman AU series<br>Dedicated \ Xylidyl Blue                                        | mmol/L | 1.52  | 1.24-1.80   |
| <b>(MYO / Mb)肌红蛋白</b>                                                                |        |       |             |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence<br>Immunoassay (ECLIA) | ug/L   | 34.57 | 27.57-41.57 |
| <b>(Na)钠</b>                                                                         |        |       |             |
| Beckman AU series<br>Dedicated \ ISE (indirect)                                      | mmol/L | 152   | 138-166     |
| <b>(NEFA)游离脂肪酸</b>                                                                   |        |       |             |
| Siemens ADVIA 1800/2400<br>Co-Health \ ACS - ACOD                                    | umol/L | 994   | 816-1,172   |
| <b>(P)磷</b>                                                                          |        |       |             |
| Beckman AU series<br>Dedicated \ Phosphomolybdate UV-rate                            | mmol/L | 3.06  | 2.54-3.58   |
| <b>(PRO/ TP)(总) 蛋白</b>                                                               |        |       |             |
| Beckman AU series<br>Dedicated \ Biuret-endpoint                                     | g/L    | 91.2  | 73.0-109.4  |
| <b>(TBA)总胆汁酸</b>                                                                     |        |       |             |
| Beckman AU series<br>Dedicated \ Enzymatic Cycling Assays Method                     | umol/L | 103.9 | 49.9-157.9  |
| <b>(TBIL)总胆红素</b>                                                                    |        |       |             |
| Beckman AU series<br>Dedicated \ Diazotization                                       | umol/L | 75.7  | 61.5-89.9   |
| <b>(TES)睾酮</b>                                                                       |        |       |             |
| Siemens Immulite1000/2000<br>Dedicated \ Chemiluminescence Immunoassay<br>(CLIA )    | nmol/L | 6.9   | 5.6-8.2     |
| <b>(TG)甘油三酯</b>                                                                      |        |       |             |
| Beckman AU series<br>Dedicated \ GPO-POD                                             | mmol/L | 4.8   | 3.8-5.8     |



|       |       |       |     |            |     |         |       |         |
|-------|-------|-------|-----|------------|-----|---------|-------|---------|
| REF   | 581-2 | 6×5mL | EXP | 2022-02-28 | LOT | 376DFDA | 更新日期： | 2018-05 |
| 3 / 3 |       |       |     |            |     |         |       |         |

|                                                                                      |        |       |             |
|--------------------------------------------------------------------------------------|--------|-------|-------------|
| 批号：376DFDA                                                                           |        |       |             |
| 项目\仪器\试剂\方法                                                                          | 单位     | 均值    | + / - 2 SD  |
| (THEO)茶碱                                                                             |        |       |             |
| Roche cobas E601/602<br>Dedicated \ Electro-Chemiluminescence<br>Immunoassay (ECLIA) | mg/L   | 24.32 | 20.12-28.52 |
| (TRF)转铁蛋白                                                                            |        |       |             |
| Beckman AU series<br>Dedicated \ Immunonephelometry (INA)                            | g/L    | 4.2   | 2.4-6.0     |
| (UA)尿酸                                                                               |        |       |             |
| Beckman AU series<br>Beckman Coulter \ Uricase colorimetry                           | umol/L | 651   | 529-773     |
| (UIBC)不饱和铁结合力                                                                        |        |       |             |
| Beckman AU series<br>Dedicated \ Ferrozine colorimetry                               | umol/L | 58.1  | 48.5-67.7   |
| (Urea)尿素                                                                             |        |       |             |
| Beckman AU series<br>Dedicated \ Urease-UV-rate                                      | mmol/L | 13.1  | 10.7-15.5   |
| (Zn)锌                                                                                |        |       |             |
| J&J Vitros Fusion 5,1 FS<br>Dedicated \ dry chemistry                                | umol/L | 33.1  | 26.7-39.5   |

注：此报告所提供的数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来，仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。