



REF	662-3	6×4.5mL	EXP	2023-09-17	LOT	3F29FF2	更新日期：	2023-07
								1 / 3

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-3

6×4.5mL

EXP

2023-09-17

LOT

3F29FF2

更新日期：2023-07

2 / 3

批号：3F29FF2			
项目\仪器\试剂方法	单位	均值	+ / - 2 SD
(BASO %)嗜碱性粒细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series			
sysmex \ Calculate	%	50	30-70
Sysmex XN series			
sysmex \ Calculate	%	10	0-20
Sysmex XS series			
sysmex \ Calculate	%	10	0-20
Sysmex XT series			
sysmex \ Calculate	%	50	25-75
(EOS %)嗜酸性粒细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	15	0-30
Sysmex XE series			
sysmex \ Calculate	%	15	0-30
Sysmex XN series			
sysmex \ Calculate	%	15	5-25
Sysmex XS series			
sysmex \ Calculate	%	15	5-25
Sysmex XT series			
sysmex \ Calculate	%	15	0-30
(HCT)血细胞比容/压积			
Mindray BC-6000 series			
Mindray \ Light scattering	%	51.6	41.2-62
Sysmex XE series			
sysmex \ Light scattering	%	52.7	42.3-63.1
Sysmex XN series			
sysmex \ Light scattering	%	51.9	41.5-62.3
Sysmex XS series			
sysmex \ Light scattering	%	52.3	41.9-62.7
Sysmex XT series			
sysmex \ Light scattering	%	51.8	41.4-62.2
(HGB/ Hb)血红蛋白			
Mindray BC-6000 series			
Mindray \ Colorimetry	g/L	171	139-203
Sysmex XE series			
sysmex \ Light scattering	g/L	168	134-202
Sysmex XN series			
sysmex \ Light scattering	g/L	169	137-201
Sysmex XS series			
sysmex \ Light scattering	g/L	167	133-201
Sysmex XT series			
sysmex \ Light scattering	g/L	173	139-207
(LYM %)淋巴细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	20	10-30
Sysmex XE series			
sysmex \ Calculate	%	20	10-30
Sysmex XN series			
sysmex \ Calculate	%	15	0-30
Sysmex XS series			
sysmex \ Calculate	%	20	10-30
Sysmex XT series			
sysmex \ Calculate	%	15	0-30
(MCH)平均红细胞血红蛋白含量			
Mindray BC-6000 series			
Mindray \ Light scattering	pg	32.2	25.8-38.6
Sysmex XE series			
sysmex \ Light scattering	pg	32.4	26.2-38.6
Sysmex XN series			
sysmex \ Light scattering	pg	32.3	25.9-38.7
Sysmex XS series			
sysmex \ Light scattering	pg	32.4	26-38.8
Sysmex XT series			
sysmex \ Light scattering	pg	32.5	26.1-38.9
(MCHC)平均红细胞血红蛋白浓度			
Mindray BC-6000 series			
Mindray \ Light scattering	g/L	328	264-392
Sysmex XE series			
sysmex \ Light scattering	g/L	332	268-396
Sysmex XN series			
sysmex \ Light scattering	g/L	326	260-392
Sysmex XS series			
sysmex \ Light scattering	g/L	337	273-401
Sysmex XT series			
sysmex \ Light scattering	g/L	335	271-399
(MCV)平均红细胞体积			
Mindray BC-6000 series			
Mindray \ Light scattering	fL	92	74-110
Sysmex XE series			
sysmex \ Light scattering	fL	94	76-112

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-3

6×4.5mL

EXP

2023-09-17

LOT

3F29FF2

更新日期：

2023-07

3 / 3

Sysmex XN series sysmex \ Light scattering	fL	93	75-111
Sysmex XS series sysmex \ Light scattering	fL	94	76-112
Sysmex XT series sysmex \ Light scattering	fL	92	74-110
(MONO %)单核细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	5	0-10
Sysmex XN series sysmex \ Calculate	%	10	0-20
Sysmex XS series sysmex \ Calculate	%	5	0-10
Sysmex XT series sysmex \ Calculate	%	10	0-20
(NEUT %)嗜中性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	55	35-75
Sysmex XE series sysmex \ Calculate	%	55	35-75
Sysmex XN series sysmex \ Calculate	%	55	35-75
Sysmex XS series sysmex \ Calculate	%	60	40-80
Sysmex XT series sysmex \ Calculate	%	55	35-75
(PLT)血小板			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	450	360-540
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E +9/L)	460	360-560
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	455	365-545
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E +9/L)	445	355-535
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	450	360-540
(RBC)红细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E +12/L)	5.67	4.53-6.81
Sysmex XE series sysmex \ Light scattering	T/L (1E +12/L)	5.68	4.56-6.8
Sysmex XN series sysmex \ Light scattering	T/L (1E +12/L)	5.68	4.54-6.82
Sysmex XS series sysmex \ Light scattering	T/L (1E +12/L)	5.69	4.57-6.81
Sysmex XT series sysmex \ Light scattering	T/L (1E +12/L)	5.71	4.67-6.75
(WBC)白细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	18.3	14.7-21.9
Sysmex XE series sysmex \ Light scattering	Giga/L (1E +9/L)	18.7	15.1-22.3
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	18.6	15-22.2
Sysmex XS series sysmex \ Light scattering	Giga/L (1E +9/L)	18.3	14.7-21.9
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	18.4	14.8-22

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。