



REF	662-2	6×4.5mL	EXP	2023-09-17	LOT	3F278E2	更新日期：	2023-07
								1 / 3

本报告含以下项目的示范实验室数据

参数 (Parameters)	参数 (Parameters)	参数 (Parameters)
(BASO %)嗜碱性粒细胞 %	(EOS %)嗜酸性粒细胞 %	(HCT)血细胞比容/压积
(HGB/ Hb)血红蛋白	(LYM %)淋巴细胞 %	(MCH)平均红细胞血红蛋白含量
(MCHC)平均红细胞血红蛋白浓度	(MCV)平均红细胞体积	(MONO %)单核细胞 %
(NEUT %)嗜中性粒细胞 %	(PLT)血小板	(RBC)红细胞
(WBC)白细胞		

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-2

6×4.5mL

EXP

2023-09-17

LOT

3F278E2

更新日期： 2023-07

2 / 3

批号：3F278E2			
项目\仪器\试剂方法	单位	均值	+ / - 2 SD
(BASO %)嗜碱性粒细胞 %			
Mindray BS series			
Mindray \ Flow cytometry	%	10	0-20
Sysmex XE series			
sysmex \ Calculate	%	50	30-70
Sysmex XS series			
sysmex \ Calculate	%	10	0-20
Sysmex XT series			
sysmex \ Calculate	%	55	35-75
Sysmex XN series			
sysmex \ Calculate	%	5	0-10
(EOS %)嗜酸性粒细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	15	0-30
Sysmex XE series			
sysmex \ Calculate	%	10	0-20
Sysmex XS series			
sysmex \ Calculate	%	15	0-30
Sysmex XT series			
sysmex \ Calculate	%	10	0-20
Sysmex XN series			
sysmex \ Calculate	%	10	0-20
(HCT)血细胞比容/压积			
Mindray BC-6000 series			
Mindray \ Light scattering	%	36.8	29.6-44
Sysmex XE series			
sysmex \ Light scattering	%	36.9	29.3-44.5
Sysmex XS series			
sysmex \ Light scattering	%	36.8	29.4-44.2
Sysmex XT series			
sysmex \ Light scattering	%	37.4	29.8-45
Sysmex XN series			
sysmex \ Light scattering	%	37.3	29.7-44.9
(HGB/ Hb)血红蛋白			
Mindray BC-6000 series			
Mindray \ Colorimetry	g/L	132	106-158
Sysmex XE series			
sysmex \ Light scattering	g/L	129	103-155
Sysmex XS series			
sysmex \ Light scattering	g/L	125	99-151
Sysmex XT series			
sysmex \ Light scattering	g/L	127	101-153
Sysmex XN series			
sysmex \ Light scattering	g/L	128	102-154
(LYM %)淋巴细胞 %			
Mindray BC-6000 series			
Mindray \ Flow cytometry	%	30	15-45
Sysmex XE series			
sysmex \ Calculate	%	30	15-45
Sysmex XS series			
sysmex \ Calculate	%	35	20-50
Sysmex XT series			
sysmex \ Calculate	%	35	20-50
Sysmex XN series			
sysmex \ Calculate	%	35	15-55
(MCH)平均红细胞血红蛋白含量			
Mindray BC-6000 series			
Mindray \ Light scattering	pg	26.8	21.6-32
Sysmex XE series			
sysmex \ Light scattering	pg	26.7	21.3-32.1
Sysmex XS series			
sysmex \ Light scattering	pg	26.8	21.6-32
Sysmex XT series			
sysmex \ Light scattering	pg	26.4	21.2-31.6
Sysmex XN series			
sysmex \ Light scattering	pg	26.6	21.4-31.8
(MCHC)平均红细胞血红蛋白浓度			
Mindray BC-6000 series			
Mindray \ Calculate	g/L	329	261-397
Sysmex XE series			
sysmex \ Light scattering	g/L	323	261-385
Sysmex XS series			
sysmex \ Light scattering	g/L	324	262-386
Sysmex XT series			
sysmex \ Light scattering	g/L	328	264-392
Sysmex XN series			
sysmex \ Light scattering	g/L	336	274-398
(MCV)平均红细胞体积			
Mindray BC-6000 series			
Mindray \ Light scattering	fL	84	68-100
Sysmex XE series			
sysmex \ Light scattering	fL	85	69-101

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。

REF

662-2

6×4.5mL

EXP

2023-09-17

LOT

3F278E2

更新日期：

2023-07

3 / 3

Sysmex XS series sysmex \ Light scattering	fL	87	71-103
Sysmex XT series sysmex \ Light scattering	fL	84	68-100
Sysmex XN series sysmex \ Light scattering	fL	86	70-102
(MONO %)单核细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	5	0-10
Sysmex XE series sysmex \ Calculate	%	5	0-10
Sysmex XS series sysmex \ Calculate	%	10	0-20
Sysmex XT series sysmex \ Calculate	%	5	0-10
Sysmex XN series sysmex \ Calculate	%	5	0-10
(NEUT %)嗜中性粒细胞 %			
Mindray BC-6000 series Mindray \ Flow cytometry	%	55	35-75
Sysmex XE series sysmex \ Calculate	%	55	35-75
Sysmex XS series sysmex \ Calculate	%	50	30-70
Sysmex XT series sysmex \ Calculate	%	50	30-70
Sysmex XN series sysmex \ Calculate	%	50	30-70
(PLT)血小板			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	250	200-300
Sysmex XE series sysmex \ Electrical impedance	Giga/L (1E +9/L)	255	205-305
Sysmex XS series sysmex \ Electrical impedance	Giga/L (1E +9/L)	250	200-300
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	245	195-295
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	245	195-295
(RBC)红细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	T/L (1E +12/L)	4.45	3.59-5.31
Sysmex XS series sysmex \ Light scattering	T/L (1E +12/L)	4.39	3.53-5.25
Sysmex XE series sysmex \ Light scattering	T/L (1E +12/L)	4.43	3.57-5.29
Sysmex XT series sysmex \ Light scattering	T/L (1E +12/L)	4.46	3.58-5.34
Sysmex XN series sysmex \ Light scattering	T/L (1E +12/L)	4.38	3.52-5.24
(WBC)白细胞			
Mindray BC-6000 series Mindray \ Electrical impedance	Giga/L (1E +9/L)	7.35	5.93-8.77
Sysmex XE series sysmex \ Light scattering	Giga/L (1E +9/L)	7.36	5.92-8.8
Sysmex XT series sysmex \ Light scattering	Giga/L (1E +9/L)	7.29	5.85-8.73
Sysmex XN series sysmex \ Light scattering	Giga/L (1E +9/L)	7.38	5.92-8.84
Sysmex XS series sysmex \ Light scattering	Giga/L (1E +9/L)	7.27	5.83-8.71

注：此报告所提供的项目与数据均基于检测相同批号质控品的若干实验室的数据汇总统计而来。此报告不可代替产品说明书。此报告仅供学习、参考之用。因所用技术、仪器和试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据。根据良好实验室规范的要求，实验室须遵循相关技术规范确立自己的均值和可接受范围。